



Sep 24, 2021

CD34+ isolation from human bone marrow

Mohsen Khosravi-Maharlooei¹, Markus Holzl¹, Austin Chen¹, Megan Sykes¹¹Columbia Center for Translational Immunology, Columbia University, New York

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.bvsw6fe

Human Islet Research Network



Lili Liang

ABSTRACT

This protocol details the steps for isolating CD34+ cells from human bone marrow. The CD34+ cells isolated from this protocol can be used for generating humanized mice through reconstitution of immune cells via IV injection after bone marrow ablation. These cells can also be used for mixed lymphocyte reaction experiments.

Corresponding Authors

Mohsen Khosravi-Maharlooei

Email: mkm2182@cumc.columbia.edu

Austin Chen

Email: ac4274@cumc.columbia.edu

Tel: 425-283-6900

DOI

dx.doi.org/10.17504/protocols.io.bvsw6fe

PROTOCOL CITATION

Mohsen Khosravi-Maharlooei, Markus Holzl, Austin Chen, Megan Sykes 2021. CD34+ isolation from human bone marrow. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bvsw6fe>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 14, 2021

LAST MODIFIED

Sep 24, 2021

PROTOCOL INTEGER ID








50742

MATERIALS TEXT

Required material





- [Bone Marrow Medium](#)
- MACS buffer(degassed)
- Sterile flask for BM rinse
- [15 and 50 mL Falcons](#)
- Histopaque
- Sterile pipetts
- CD34+ MACS kit ([130-046-702](#))
- [Human serum](#)
- MACS supply
- antibodies

Required Buffers

- BM Medium ( **500 mL** Media 199,  **5 mL** Hepes,  **5 mL** DNase,  **40 µl** Gentamycin)
- MACS buffer ( **500 mL** PBS,  **5 g** BSA,  **2 mL** EDTA, sterile filtrated and degassed)
- Cryomedium ( **90 mL** PBS,  **10 mL** FBS,  **10 mL** DMSO)

BEFORE STARTING

Human bone marrow is a rich source for CD34+ hematopoietic stem cells. CD34+ cells can be easily isolated and further processed.

- 1 Transfer the content of the collection bag into a sterile flask
- 2 Add  **250 mL** BM Medium to the bag and rinse it thoroughly
- 3 Transfer the content of the collection bag into the sterile flask
- 4 Layer  **35 mL** of the suspension over  **15 mL** of Histopaque
- 5 Centrifuge the tubes for 30 minutes  **500 g** without brake at RT

- 6 Collect the leukocyte ring in **50 mL** Falcons and fill up with BM medium
- 7 Wash once by centrifuging 6 minutes **500 g**
- 8 Resuspend the cells in MACS buffer and count (*take **50 µl** for FACS confirmation = PRE*)
- 9 Wash down again and resuspend the pellet according to the protocol ([130-046-702](#) MACS Human CD34+ kit)
- 10 Add **300 µl** of MACS buffer per 10^8 cells
- 11 Add **100 µl** of FcR-B reagent per 10^8 cells, mix it and incubate in fridge for 15 minutes
- 12 Add **100 µl** of CD34 beads per 10^8 cells to the suspension and mix and let it sit for 30 min (fridge)
- 13 Fill up with **50 mL** MACS buffer and strain through a blue strainer (40 µm)
- 14 Wash cells (**500 g** 6 min) and resuspend in MACS buffer 500 µL/200,000,000 cells. If you have more cells, increase volume accordingly. E.g 3×10^9 cells = 7,5mL. Aliquot this volume to more than one (with **3 mL** prerinsed) MACS column.
- 15 Wash with buffer **3 mL** 3 times and keep negative fraction (*Take **50 µl** for FCM = POST neg*)
- 16 Put the column out of the magnet and push out positive fraction with **5 mL** Buffer and the plunger
- 17 Collect the positive fractions. (*Take **50 µl** for FCM = Post pos*)
- 18 Process the cells as desired (injection in mouse or cryopreservation)

19 Check the purity with FACS

20

<u>FACS panel</u>		
CD45	FITC	5 ul
CD3	PercpC5.5	5 ul
CD14	PacBlue	5 ul
CD19	APC	5 ul
CD34	PE	5 ul
CD38	PeC7	5 ul
Human Serum		5 ul
<u>FACSbuffer</u>		<u>15 ul</u>
Total		50 ul/sample