



JUN 19, 2023

TSS competent cells and transformation

David Ruano Gallego¹

¹Universidad Autónoma de Madrid



David Ruano Gallego

Universidad Autónoma de Madrid

ABSTRACT

This protocol allows the efficient transformation of plasmid into E. coli and related strains.

OPEN  ACCESS

DOI:

dx.doi.org/10.17504/protocols.io.5jyl8jzbrg2w/v1

Protocol Citation: David Ruano Gallego 2023. TSS competent cells and transformation .

protocols.io

<https://dx.doi.org/10.17504/protocols.io.5jyl8jzbrg2w/v1>

MANUSCRIPT CITATION:

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: Apr 11, 2023

Last Modified: Jun 19, 2023

PROTOCOL integer ID:

80304

MATERIALS

Transformation Storage Solution (TSS) buffer

- PEG3350 10%. Sigma. Ref: 202444
- DMSO 5%.



DMSO (dimethyl sulfoxide) Merck MilliporeSigma (Sigma-Aldrich) Catalog #D8418

- Glycerol 10%



Glycerol - for molecular biology, $\geq 99\%$ Merck MilliporeSigma (Sigma-Aldrich) Catalog #G5516

- MgCl_2 20 mM.



1 M Magnesium Chloride (MgCl_2) Merck MilliporeSigma (Sigma-Aldrich) Catalog #M8266

- LB medium 2x.

Adjust to pH 6.1 with 6 M HCl.

Autoclave.

Add previously sterilized MnCl_2 4 M to a final concentration of 140 mM.

Add water to adjust the modified LB medium to 1x.

Prepare aliquots of 1 ml and store at -20°C .

5x KCM solution

- KCl 0.5 M.



Potassium Chloride Merck MilliporeSigma (Sigma-Aldrich) Catalog #P9541

- CaCl_2 150 mM.



Calcium chloride dihydrate Merck MilliporeSigma (Sigma-Aldrich) Catalog #C7902

- MgCl_2 250 mM.



Magnesium chloride Merck MilliporeSigma (Sigma-Aldrich) Catalog #M8266

Sterilize at 121°C for 20 min. Store at 4°C .

BEFORE START INSTRUCTIONS

Prepare solution TSS and KCM.


When prepared, TSS has a cloudy appearance that disappears once autoclaved. Also, over time, TSS can seem to be contaminated without being so, since the polyethylene glycol tends to precipitate.

Preparing competent cells

2d

- 1 Grow o/n a 5 mL culture of the strain to be transformed.

12h

 5 mL

 37 °C

 Overnight

- 2 Dilute 1:100 of the preinoculum in 50 mL prewarmed LB ($OD_{600} = 0,05$).

 50 mL

 37 °C


- 3 Grow at 37°C with agitation (170 rpm) until it reaches $OD_{600} = 0,5$.

1h 30m


 170 rpm, 37°C


- 4 Centrifuge the culture at 4000 g and 4°C for 10 min.

10m

 4000 x g, 4°C,
00:10:00


- 5 Discard supernatant and wash with 1 mL of chilled TSS.

 1 mL

 On ice

- 6 Centrifuge at 4000 g and 4°C for 3 min.

3m

 4000 x g, 4°C,
00:03:00

- 7 Discard supernatant and resuspend in 1 mL of TSS.

🧪 1 mL

🌡️ On ice

8 Incubate 10 min on ice.

10m

🌡️ On ice

⌚ 00:10:00

9 Store 30 µL aliquots at -80°C.

🧪 30 µL

🌡️ -80 °C

Transformation protocol

1h 40m

10 Mix 5 µL of 5x KCM solution with DNA (around 1-5 µL) and H₂O to a total of a 25 µL mixture.

🧪 25 µL

11 Mix the mixture with 25 µL of competent cells.

🧪 25 µL

12 Incubate 30 min on ice.

30m

🌡️ On ice

⌚ 00:30:00

13 Thermal shock: incubate at 42°C for 90 s.

1m 30s

🌡️ 42 °C

⌚ 00:01:30

14 Incubate 2 min on ice.

2m

🌡️ On ice

⌚ 00:02:00

15 Recovery: add 250 µl of LB and incubate at 37°C for 1 h.

1h

 37 °C  01:00:00

- 16** Plate 100 µl of the solution on LB agar plates supplemented with the indicated antibiotic and incubate o/n at 37°C.

12h

 37 °C  Overnight