



SEP 23, 2023

OPEN ACCESS



DOI:
dx.doi.org/10.17504/protocols.io.n2bvj3y5blk5/v1

Protocol Citation: Elias Adriaenssens 2023. Cell culture. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n2bvj3y5blk5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
 We use this protocol and it's working

Created: Jul 07, 2023

Last Modified: Sep 23, 2023

PROTOCOL integer ID:
 84663

Keywords: Cell culture

Cell culture

Elias Adriaenssens¹

¹Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna

ASAP Collaborative Research Network

Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

ABSTRACT

This protocol describes routine cell culture maintenance.

MATERIALS

Reagents

- Dulbecco Modified Eagle Medium (DMEM, Thermo Fisher)
- Fetal Bovine Serum (FBS, Thermo Fisher)
- HEPES (Thermo Fisher)
- non-essential amino acids (NEAA, Thermo Fisher)
- Penicillin-Streptomycin (Thermo Fisher)
- Iscove's Modified Dulbecco's Medium (IMDM, Thermo Fisher)

Cell culture

5m



- 1 Grow HeLa and HEK293T cells in Dulbecco Modified Eagle Medium (DMEM, Thermo Fisher) supplemented with 10% (v/v) Fetal Bovine Serum (FBS, Thermo Fisher), **[M] 25 millimolar (mM)** HEPES (Thermo Fisher), 1% (v/v) non-essenMal amino acids (NEAA, Thermo Fisher), and 1% (v/v) Penicillin-Streptomycin (Thermo Fisher).
- 2 Culture HAP1 cells in Iscove's Modified Dulbecco's Medium (IMDM, Thermo Fisher), supplemented with 10% (v/v) Fetal Bovine Serum

(FBS, Thermo Fisher) and 1% (v/v) Penicillin-Streptomycin (Thermo Fisher).

3 Regularly test all cell lines for mycoplasma contaminations.

4 Split the cells by trypsinization once confluent.

5 To this end, remove medium with Pasteur pipette.

6 Wash cells with 1x PBS and then add trypsin for  00:05:00 at  37 °C .

5m



7 Afterwards, resuspend cells in medium and retain 10% for further culturing and expansion.

8 Discard the rest or collect for seeding for experiments.