



Jul 17, 2024

Agrobacterium-Mediated Transient Expression in *Nicotiana benthamiana* Leaves



DOI

dx.doi.org/10.17504/protocols.io.rm7vzj7prlx1/v1

Maryam Rahmati Ishka¹

¹Boyce Thompson Institute (BTI)



Maryam Rahmati Ishka

Boyce Thompson Institute (BTI)

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.rm7vzj7prlx1/v1

External link: <http://n/a>

Protocol Citation: Maryam Rahmati Ishka 2024. Agrobacterium-Mediated Transient Expression in *Nicotiana benthamiana* Leaves. protocols.io <https://dx.doi.org/10.17504/protocols.io.rm7vzj7prlx1/v1>

Manuscript citation:
n/a

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 17, 2024

Last Modified: July 17, 2024

Protocol Integer ID: 103602

Keywords: Agrobacterium infiltration, *Nicotiana benthamiana*, plant construct delivery, transient protein expression



Disclaimer

n/a

Abstract

This procedure outlines the *Agrobacterium* infiltration of tobacco leaves for the transient expression of fluorescently labeled proteins.


Image Attribution

n/a

Guidelines

n/a

Materials

Infiltration medium: 10 mM MgCl₂, 10 mM MES (pH 5.6). Autoclave for 20 min. Add acetosyringone immediately before using it  MgCl₂ **Applied Biosystems (ThermoFisher Scientific)**

 20X MES Buffer **Thermo Fisher Scientific Catalog #NP0002**

 Acetosyringone **Merck MilliporeSigma (Sigma-Aldrich) Catalog #2478-38-8**

 1g Rifampicin USP Grade **G-Biosciences Catalog #RC-191** to a final concentration of 200 µM.

Safety warnings

 n/a

Ethics statement

n/a

Before start

n/a

Transform constructs into *Agrobacterium tumefaciens* strains that are compatible with vectors. Strain GV3101 is typically used (GV2260 can also be used). Make glycerol stock of.

- 1 Refresh bacteria containing the construct of interest from glycerol stock on agar YEP (or any other bacterial media such as LB) medium with an appropriate antibiotic. Incubate at 28°C for two days. If *A. tumefaciens* GV3101 is used add **rifampicin**.
- 2 Inoculate **5 mL** of liquid YEP medium with appropriate **antibiotic** in a test tube. Grow at 28°C for **24 hours** with shaking at 250-300 rpm. Be sure to grow bacterial strain with **p19** as a silencing inhibitor as well!
- 3 Measure the OD₆₀₀ using 100 µL of bacterial culture.
- 4 Calculate the culture volume needed for an OD₆₀₀ of 0.2 - 0.5 per construct for a 4 mL total volume of Infiltration medium. 4 mL are used for test spot infiltration. The specific OD₆₀₀ to use depends on how well expressed your construct is in *N. benth*, the general range is 0.1 – 0.5. Use an OD₆₀₀ of 0.3 for P19.
(I use 0.3 for all constructs and it works well).
- 5 Pipette the appropriate volume of each strain together in 2 mL tubes.
- 6 Centrifuge combined bacteria at 8000×*g* for 2 min at room temperature. Carefully remove all supernatant.
- 7 Add 200 mM **acetosyringone** to the **Infiltration medium** at a 1:1000 ratio.
Note: acetosyringone is light-sensitive, keep out of light until ready for use.
- 8 Add 2 mL of infiltration medium with **acetosyringone** to each pellet, and resuspend. Leave in the dark at room temperature for 1–2 hrs (or even longer but not overnight).
- 9 Infiltrate 4-week-old *N. benthamiana* (pre-flowering) leaf using 1 mL syringe.
- 10 Grow plants for 2-3 days after infiltration (up to 5 days, depending on expression level of constructs).



- 11 Punch **1 cm** leaf discs from leaves and observe under the microscope to visualize the fluorescence.

Protocol references

n/a