



MAR 29, 2024

# 🌐 Immunofluorescent Staining of Fixed Mouse Brain Tissue Sections

Katerina Rademacher<sup>1</sup>, Ken Nakamura<sup>1</sup>

<sup>1</sup>Gladstone Institute of Neurological Disease



Haru Yamamoto  
UCSF

## ABSTRACT

This protocol describes steps for immunofluorescent staining of free floating fixed mouse brain tissue sections.

## MATERIALS

- 1X PBS (UCSF Stem Cell Core)
- Triton X-100 (Sigma Aldrich 93443)
- Fetal bovine serum (UCSF Stem Cell Core)
- Primary antibody
- Secondary antibody
- Blocking buffer: 0.2% Triton, 4% FBS in 1X PBS
- Glass microscope slides (Fisher 22-038-103)
- Anti-fade mounting medium with or without DAPI (Vector Laboratories H140010)

OPEN  ACCESS



DOI:

[dx.doi.org/10.17504/protocols.io.kxygx38owg8j/v1](https://dx.doi.org/10.17504/protocols.io.kxygx38owg8j/v1)

**Protocol Citation:** Katerina Rademacher, Ken Nakamura 2024. Immunofluorescent Staining of Fixed Mouse Brain Tissue Sections. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kxygx38owg8j/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Jan 07, 2024

**Last Modified:** Mar 29, 2024

**PROTOCOL integer ID:** 93032

**Keywords:** ASAPCRN

**Funders Acknowledgement:**

ASAP

Grant ID: 020529

**1** Free floating sections in 12-well plates, all steps performed at room temperature.

**1.1** 2 x 10min PBS washes.

**1.2** 2 x 10min 0.2% PBS-Triton (PBS-T) washes.

**1.3** Block 1hr in blocking buffer.

**1.4** Incubate in primary antibody diluted in blocking buffer overnight.

**1.5** 3 x 10min 0.2% PBS-T washes.

**1.6** Incubate in secondary antibody diluted in blocking buffer for 2hrs.

**1.7** 2 x 10min PBS-T washes.

**1.8** 2 x 10min PBS washes.

**1.9** Mount sections on glass slides, add mounting media and coverslip.