

May 10, 2024

## AAV Craniotomy

DOI

[dx.doi.org/10.17504/protocols.io.5jyl829q7l2w/v1](https://dx.doi.org/10.17504/protocols.io.5jyl829q7l2w/v1)

Santiago Unda<sup>1</sup>, Roberta Marongiu<sup>1</sup>, Mihaela Stavarache<sup>1</sup>, Michael G. Kaplitt<sup>1</sup>

<sup>1</sup>Weill Cornell Medical College

ASAP Collaborative Rese...

Kaplitt Protocols



Eileen Ruth Torres

Weill Cornell Medicine

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.5jyl829q7l2w/v1](https://dx.doi.org/10.17504/protocols.io.5jyl829q7l2w/v1)

**Protocol Citation:** Santiago Unda, Roberta Marongiu, Mihaela Stavarache, Michael G. Kaplitt 2024. AAV Craniotomy. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.5jyl829q7l2w/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 08, 2024

**Last Modified:** May 10, 2024

**Protocol Integer ID:** 99452

**Keywords:** ASAPCRN, surgery, mouse, AAV

**Funders Acknowledgement:**

**Aligning Science Across**

**Parkinson's**

**Grant ID:** 020608

## Abstract

This protocol is for injection of AAV into a mouse brain using a stereotactic frame and Hamilton syringe.

## Guidelines

Note that sterile technique must be followed during this procedure to avoid health complications.

## Anesthesia

- 1 Weigh mice.
- 2 Administer intraperitoneal mixture of ketamine (ButlerAnimal Health Supply) and xylazine (Lloyd Laboratories) at concentrations of 110 and 10 mg/kg of body weight, respectively.
- 3 Confirm induction of anesthesia.

## Injection

- 4 Place mouse into stereotactic frame (David Kopf Instruments) and place eye lubricant.
- 5 Remove fur from skull area using either electric shaver or depilatory cream.
- 6 Sterilize skull area and then cut open with scalpel.
- 7 Adjust brain so that it is level. Locate bregma and zero X and Y axes on stereotaxic frame.
- 8 Move to injection coordinates and drill small hole into the skull.
- 9 Zero the Z axis on the stereotaxic frame.
- 10 Prepare AAV solution by diluting AAV in PBS with 500mM Mg<sub>2</sub>/Cl to a final dilution of total of  $2 \times 10^{11}$  genomic particles/mL.
- 11 Using a 10  $\mu$ L stereotactic syringe with a 33 G needle attached to a micro-infusion pump (World Precision Instruments), withdraw the AAV solution.



- 12 With the loaded syringe, navigate the needle to the desired X, Y, Z coordinates.
- 13 Infuse the brain with AAV at a flow rate of 0.1 - 0.4  $\mu\text{L}/\text{min}$ .
- 14 To prevent reflux, after each infusion, leave the injection needle in place for 5 min, withdraw a short distance (0.3 - 0.5 mm), and then leave in the new position for an additional 2 min before complete removal of the needle.
- 15 Close skin with non-absorbable sutures.

## AAV Preparation

- 16 Inject a total of  $2 \times 10^{11}$  genomic particles/mL of vector is injected into the appropriate coordinates.

## Post-procedure

- 17 Give meloxicam 2mg/kg subcutaneously.
- 18 Perform health checks every day for the first 72 hrs post-injection.
- 19 Remove sutures 15 days post-injection.