

Aug 18, 2022

# Growth Curve Stress Test

Brian Teague<sup>1</sup><sup>1</sup>University of Wisconsin - Stout

1 Works for me

Share

This protocol is published without a DOI.

Yeast ORFans CURE

Brian Teague  
University of Wisconsin - Stout

## ABSTRACT

This is one (fairly quantitative and reproducible) way to run a yeast growth study to analyze the effect of a stressor on yeast. It is not the only way! It is also less of a "protocol" and more of a "template" -- it measures the growth of a single yeast strain at one concentration of "chemical X." Makes sure to adapt it to your own needs!

## PROTOCOL CITATION

Brian Teague 2022. Growth Curve Stress Test. **protocols.io**  
<https://protocols.io/view/growth-curve-stress-test-cffjtjkn>



## KEYWORDS

yeast, phenotype, growth curve

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## IMAGE ATTRIBUTION

Chalik1, CC BY 3.0 , via Wikimedia Commons

## CREATED

Aug 17, 2022

## LAST MODIFIED

Aug 18, 2022

## PROTOCOL INTEGER ID

68811

## PARENT PROTOCOLS

In steps of

### Growth Curve Stress Test (Instructor Protocol)

#### GUIDELINES

Remember, if you see a difference between the knockout and the wild-type under stress, it could be that the stress is slowing down the knockout strain. However, there's another possible explanation: it could just be that the knockout strain is growing more slowly under all conditions. Hence the no-stress controls.

#### MATERIALS TEXT

##### Equipment

- Vortexer
- Microplate reader (with incubation and shaking)
- Sterile workspace (such as a biosafety cabinet or laminar flow hood)

##### Materials

- Corning COSTAR 96-well clear flat-bottom assay plate
- Wild-type yeast in double-strength synthetic complete media, OD600 ~ 0.1-0.2
- Knockout yeast in double-strength synthetic complete media, OD600 ~ 0.1-0.2
- A solution of Chemical X in water

#### SAFETY WARNINGS

Yeast and yeast media are not hazardous.

However, you may be using a hazardous "chemical X." Make sure to use appropriate PPE and engineering controls and dispose of it appropriately.

The yeast we are using are genetically modified. Make sure to inactivate them (by bleach or autoclaving, for example) before disposing of them.

#### Assay setup

1

If possible, perform steps 1 and 2 in a sterile environment such as a biosafety cabinet or laminar flow hood.

Vortex the yeast culture briefly to resuspend the yeast cells.


2 1. Set up four wells according to the following table:

A	B	C	D	E
	Well 1	Well 2	Well 3	Well 4
Yeast culture (wild-type)	50 $\mu$ l	50 $\mu$ l	–	–
Yeast culture (knockout)	–	–	50 $\mu$ l	50 $\mu$ l
Chemical X, 100 mM	10 $\mu$ l	–	10 $\mu$ l	–
Water	40 $\mu$ l	50 $\mu$ l	40 $\mu$ l	50 $\mu$ l

- 3 Set up the plate reader as follows:
  - Temperature: 30°C
  - Mode: Kinetic
  - Wavelength: 600 nm
  - Interval: 5 minutes
  - Total run time: 24 hours
  - Shake before read: 30 seconds
- 4 Transfer the assay plate to the reader and read for 24 hours.

#### Data analysis

- 5 Plot each growth curve with the Y axis scale set to "logarithmic."



- 6 Find the part of the curve that is linear.



In this case, I might choose from 250 minutes to 350 minutes.

- 7 Find the OD600 reading at these two times.



At 250 minutes, my OD600 is 0.2086; at 350 minutes, my OD600 is 0.3583.

- 8 Compute the doubling time using a bit of math:

$$T_d = (t_2 - t_1) \cdot \frac{\ln(2)}{\ln\left(\frac{q_2}{q_1}\right)}.$$

Where  $t_1$  and  $t_2$  are the start and end times you found in step 2, and  $q_1$  and  $q_2$  are the OD600 at those times.



Thus, for my BY4735 strain growing in YPD, the doubling time is 128 minutes.