



FEB 27, 2024

## Quantitative PCR, 384 well format

Malu G Tansey<sup>1</sup>

<sup>1</sup>College of Medicine |University of Florida



Senthilkumar Karuppagounder

### ABSTRACT

Quantitative PCR, 384 well format

### MATERIALS

#### Materials:

- PCR 8-tube strips or 96 well plates depending on sample# and gene# (*General stores 4868*)
- PCR strip tube lids (*General stores 1400-0800*)
- ABI Sybr Green (*ABI 4301955*)
- QPCR plate (*ABI 4309849*)
- Adhesive plate film (*ABI 4311971*)
- Brown plastic adhesive helper
- Filtered 200 µl pipette tips (2 boxes for one complete plate)
- Electronic Multichannel Pipette (referred to as Hilda from here on out)
- Primers – final concentration of 1.25 µM each
- Autoclaved MilliQ H<sub>2</sub>O

### OPEN ACCESS



#### DOI:

[dx.doi.org/10.17504/protocols.io.yxmvm341bl3p/v1](https://dx.doi.org/10.17504/protocols.io.yxmvm341bl3p/v1)

**Protocol Citation:** Malu G Tansey 2024. Quantitative PCR, 384 well format. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.yxmvm341bl3p/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Jan 29, 2024

Last Modified: Feb 27, 2024

PROTOCOL integer ID: 94312

Keywords: ASAPCRN

**Funders Acknowledgement:**

Aligning Science Across  
Parkinson's (ASAP)  
Grant ID: ASAP-020527

## Set-up:

1 Sample **cDNA**: thaw your cDNA on ice (*vortex and quick spin before plating*)

1.1 Primer stock (+/-): Add 12.5 µl of (+) and 12.5 µl (-) to 1 mL MilliQ H<sub>2</sub>O

1.2 PCR Master Mix (10ul rxn; Triplicate)

A	B
MilliQ H <sub>2</sub> O	10.20ul
SYBR Green	20ul
primers (2.5 µM each set)	4.80ul
Total master mix	35.00ul

Equipment: ABI 7900 Prism

## Procedure:

2 Place 8-tube PCR strips in PCR tube racks (each single tube runs 1 sample and 1 gene)

- 3** Add 5  $\mu$ l cDNA to the bottom of each tube, use 20  $\mu$ L pipetor (keep on ice)
- 4** Add 25  $\mu$ l (duplicates) or 35  $\mu$ l (triplicates) of Sybr Mastermix to the 8-well tube (keep on ice)
  - 4.1** Use 200 $\mu$ L multi-channel pipette.
- 5** Mix using the multichannel and quick spin.
- 6** Dispense 10  $\mu$ l of cDNA/Mastermix (20 $\mu$ L pipetor) into each well on 384 well plate according to plate layout made in advance (keep plate on ice)
- 7** Gently blot top of plate with kimwipe (to keep samples from transferring to other wells)
- 8** Place clear Adhesive plate cover over the plate.
  - 8.1** use brown 'helper' to smooth out

**8.2** pay attention to edges

**8.3** work from center of the plate out

**9** Spin plate for 5 min at 3500 rpm (4° C)

**9.1** During spin: set up ABI SDS program (keep plate in centrifuge until ready to run)

**10** Seal plate with sticky film. Vortex and spin down plate 3500rpm for 5 min at 4C

**11** Open SDS 2.3 program

**11.1** File -> new

**11.2** One instrument tab: real time -> Connect to machine -> open/close door

**11.3** Insert plate, aligning A1 to A1

**11.4** Close door

**11.5** On layout tab; highlight unused wells, click “omit wells”

**11.6** Highlight used wells and click “add detector” for each specific gene

**11.7** Set to 10uL Rxn VL

**11.8** Check cycle times and temperatures

**11.9** Add dissociation stage (SYBR primers only)

	A	B	C	D	E	F	G	H
Temp C	50	95	95	60	95	60	95	
Time	2:00	10:00	0:15	1:00	0:15	0:15	0:55	

Stage C&D 40 time

- 12 Run plate
- 13 Primer Validation Procedure: Set-up is same as above plus cDNA standard curve for each gene in an extra set of 8-tube PCR strips (see workflow file)