

# Freezing cancer cell lines

Emily Souster<sup>1</sup>, Verity Goodwin<sup>1</sup>, Adam Jackson<sup>1</sup>, Charlotte Beaver<sup>1</sup>, Rizwan Ansari<sup>1</sup>, Fiona Behan<sup>1</sup>, Mathew Garnett<sup>1</sup>

<sup>1</sup>Wellcome Sanger Institute

1 Works for me dx.doi.org/10.17504/protocols.io.bgtyjwpw


Cellular Generation and Phenotyping

Emily Souster

## ABSTRACT

This protocol outlines routine banking of cancer cell lines and Ca9 transduced cancer lines.

Process diagram:



DOI

dx.doi.org/10.17504/protocols.io.bgtyjwpw

## PROTOCOL CITATION

Emily Souster, Verity Goodwin, Adam Jackson, Charlotte Beaver, Rizwan Ansari, Fiona Behan, Mathew Garnett 2020. Freezing cancer cell lines. **protocols.io**  
dx.doi.org/10.17504/protocols.io.bgtyjwpw

## COLLECTIONS

Whole genome CRISPR screening in cancer cell lines

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

May 25, 2020

## LAST MODIFIED

Jul 23, 2020

PROTOCOL INTEGER ID

37464

PARENT PROTOCOLS

Part of collection

[Whole genome CRISPR screening in cancer cell lines](#)

GUIDELINES

- As a guideline, we usually bank 5 cryovials from a 70% confluent T150 flask, each containing 1ml cell suspension.


MATERIALS

NAME	CATALOG #	VENDOR
<a href="#">Falcon™ 15mL Conical Centrifuge Tubes</a>	14-959-53A	Fisher Scientific
<a href="#">TrypLE™ Express Enzyme (1X), no phenol red</a>	12604021	Thermo Fisher
<a href="#">Nunc™ Biobanking and Cell Culture Cryogenic Tubes, 1.8mL, 48mm, external thread, printed</a>	375418	Thermo Fisher
DMSO	D2650	Sigma Aldrich
DPBS	14190	Invitrogen - Thermo Fisher

MATERIALS TEXT


Select an appropriate culture media for your cell line. Common culture medias used for cancer cell lines are serum supplemented Advanced DMEM F-12 or RPMI, in the presence of pen-strep.

### Equipment

Light Microscope  
Microbiological Safety Cabinet (MSC)  
Pipette Boy  
Stripettes  
Pipettes and tips  
 **37 °C** , 5% CO<sub>2</sub> incubator  
Centrifuge  
CoolCell or appropriate freezing container  
-80C freezer  
Liquid Nitrogen storage

BEFORE STARTING

- Pre-warm complete culture media to room-temperature.
- Check the cells under the microscope and record percentage confluency. Cancer cells should be banked when ~70% confluent.

1 Prepare  **1 mL** freezing media per vial as follows: complete culture media + 10% DMSO.

2 Detach and collect cells from a flask, by following Steps 1-6 of the protocol: [Passaging adherent cancer cell lines](#).

- 3 Aspirate the supernatant, taking care to avoid disturbing the cell pellet. Resuspend the pellet in an appropriate volume of **freezing media**- depending on the number of vials being frozen. Mix well to ensure a single cell suspension.

For example, if 5 vials are being frozen from a T150, resuspend the cell pellet in **5 mL** of freezing media.

- 4 Transfer **1 mL** aliquots of the cell suspension to pre-labelled 1.8ml cryovials.
- 5 Place vials in a 'CoolCell' or appropriate freezing container and store at **-80 °C** overnight.



Appropriate freezing containers will ensure that the liquid freezes at a controlled rate of around **-1 °C** per minute at **-80 °C**.

- 6 Transfer the vials to liquid nitrogen for long-term storage.