



Jul 20, 2022

ADP Glo Max measuring ATP13A2 ATPase activity in microsomes

Sue Sim¹, eunyong_park¹¹University of California, Berkeley

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.6qpvr65ozvmk/v1 Sue Sim

ABSTRACT

Measuring ATP13A2 ATPase activity in microsomes using commercial luminescence-based ADP detection assay (ADP Glo Max; Promega)

DOI

dx.doi.org/10.17504/protocols.io.6qpvr65ozvmk/v1

PROTOCOL CITATION

Sue Sim, eunyong_park 2022. ADP Glo Max measuring ATP13A2 ATPase activity in microsomes. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.6qpvr65ozvmk/v1>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 17, 2022

LAST MODIFIED

Jul 20, 2022

PROTOCOL INTEGER ID

66882

MATERIALS TEXT

Luminescence-based ADP detection assay (ADP Glo Max assay; Promega)

Reaction Buffer

50 mM MOPS-KOH pH 7.0

100 mM KCl

11 mM MgCl₂

1 mM DTT

0.02% DDM

0.001-10 mM spermine

- 1 Thaw microsomes on ice
- 2 Use 1 µg microsomes per reaction (1 µL of 5 mg/mL stock)
- 3 Dilute to 4 µL in Reaction Buffer at 4°C
- 4 Equilibrate microsomes in reaction buffer for ⌚01:30:00 🌡️ On ice and for ⌚00:10:00^{1h 40m} at 🌡️ 37 °C
- 5 Add 1 µL of 25 mM ATP (final concentration 5 mM ATP) to start reaction (total reaction volume 5 µL)
- 6 Incubate for ⌚00:15:00 at 🌡️ 37 °C upon addition of ATP 15m
- 7 Terminate reaction by heating samples for ⌚00:05:00 at 🌡️ 80 °C 5m
- 8 Mix with 5 µL of ADP-Glo Reagent for 60-80 min at 🌡️ 23 °C

- 9 Mix with 10 uL of ADP-Glo Max Detection Reagent for 🕒 01:00:00
- 10 Measure luminescence in 384 multi-well plate using a luminometer (NOVOstar; BMG Labtech)