



Aug 19, 2021

PEI/Laminin Coating

In 1 collection

Zoe¹¹Khurana Lab, Harvard Medical School, Brigham and Women's Hospital

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.bu5iny4e

hendersa

ABSTRACT

This protocol is used to prepare plates for iPSC-derived neuron maturation. It is a part of the NGN2 cortical neuron differentiation protocol: coat plates in preparation for Day 7 dissociation and replating.

ATTACHMENTS

[dh4fbia7.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.bu5iny4e

PROTOCOL CITATION

Zoe 2021. PEI/Laminin Coating. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bu5iny4e>

COLLECTIONS ⓘ

HCNT-Cortical Neuron Transdifferentiation Overview

KEYWORDS

PEI, Laminin Coating, iPSC-derived neuron maturation

LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

May 18, 2021

LAST MODIFIED

Aug 19, 2021

OWNERSHIP HISTORY

May 18, 2021 Urmilas

Jun 23, 2021 hendersa



PROTOCOL INTEGER ID

50058











Materials

-  **50 mL** Falcon tube
[SteriFlip-GP Sterile Centrifuge Tube Top Filter](#)
- [Unit Millipore Catalog #SCGP00525](#)
- **96 square-well plates (Brooks, #MGB096-1-2-LG-L)**
- Cell culture hood
- Multichannel pipet
- **Reagent Reservoir (Fisher Scientific, #21385104)**

Reagents

- [Pierce™ 20X Borate Buffer Thermo](#)
- [Fisher Catalog #28341](#)
- Distilled water ()
- 5% PEI (poly-ethylenimine, stored at  **4 °C**)
- [Laminin Sigma Catalog #L2020](#) , ^[M] **1 mg/ml**
[DPBS with Ca and Mg Gibco - Thermo](#)
- [Fischer Catalog #14040-133](#)
[DPBS \(Dulbeccos Phosphate Buffered Saline\) 1x \[-\] calcium](#)
- [magnesium Corning Catalog #21-031-CV](#)
 **500 mL**

PEI Coating

- 1 At least one day before plating, make 2x borate buffer with  **4 mL** 20x borate buffer (Thermo, cat. no. 28341) in  **36 mL** water in a  **50 mL** Falcon tube.
- 2 Dilute  **800 µl** 5% PEI (poly-ethylenimine; stored at  **4 °C**) in  **40 mL** 2x borate buffer to make 0.1% PEI.
- 3 Filter-sterilize with SteriFlip: attach filter and tube to  **50 mL** tube with PEI solution; flip and attach to vacuum; when medium drains through, remove filter and original tube.
- 4 
Add  **150 µl** to each of the inner 60 wells of a 96 square-well plate (Brooks, cat. no. MGB096-1-2-LG-L);  **500 µl** for wells of 24-well plate (Corning, #3527).

Note: Poly-ornithine-laminin pre-coated plates (6-well, 24-well) can be substituted for this plating protocol-only 96-well must be PEI/laminin coated.

- 5 On 96-well plates, leave a border of empty wells to avoid edge effects; fill these with PBS to maintain humidity.



Label plate(s), wrap sides with parafilm and cover in plastic wrap; store at 4°C Overnight .

Laminin Coating

- 7 On the day of passaging, make laminin $5\text{ }\mu\text{g/ml}$ in 4°C PBS (with calcium and magnesium ions), ex. 1 mL PBS with $5\text{ }\mu\text{l}$ laminin (Sigma, cat. no. L2020, 1 mg/ml).



Wash PEI-coated plates twice with distilled water ($300\text{ }\mu\text{l}$ for 96-well plates).



Wash once with PBS.



Add $150\text{ }\mu\text{l}$ laminin to wells of 96-well plate, and 0.5 mL to wells of 12-well.



Incubate plates at least 2 hours at 37°C before re-plating neurons.