

VERSION 2

APR 06, 2023

OPEN ACCESS

**DOI:**  
[dx.doi.org/10.17504/protocols.io.q26g74r19gwz/v2](https://dx.doi.org/10.17504/protocols.io.q26g74r19gwz/v2)

**Protocol Citation:** Ying-Yu Hu, Zoe V. Finkel 2023. Estimate phospholipids from microalgae. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g74r19gwz/v2> Version created by Ying-Yu Hu

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
 We use this protocol and it's working

**Created:** Feb 16, 2023

**Last Modified:** Apr 06, 2023

**PROTOCOL integer ID:**  
 77080

**Keywords:** phospholipids, high temperature dry combustion

## Estimate phospholipids from microalgae V.2

Ying-Yu Hu<sup>1</sup>, Zoe V. Finkel<sup>1</sup>

<sup>1</sup>Dalhousie University

Marine Microbial Macroecology Lab  
 Tech. support email: [ruby.hu@dal.ca](mailto:ruby.hu@dal.ca)



Ying-Yu Hu  
 Dalhousie University

### ABSTRACT

Here we describe a protocol to estimate phospholipids from microalgae.

After extracting and measuring the total lipids from microalgae, the remaining lipid extract is dried using a nitrogen flow, followed by drying with magnesium sulfate at 90°C. However, it has been observed that traditional dry combustion at 500°C only decomposes approximately 50% of phospholipids (Hu et al., 2022). To achieve complete conversion of phospholipids to pyrophosphate, a temperature of around 800°C is required, but such high temperatures cannot be used with glassware. As the acid digestion method involves using only 500 µL of 0.2 M HCl, which must be placed in tightly capped glass vials to prevent concentration changes due to evaporation, combustion must be carried out using glassware instead of crucibles. It should be noted that the recovery rate of phospholipids is around 80% when combusted at 650°C, but this recovery rate is consistent, making the use of glass vials applicable. Therefore, we recommend using 650°C to combust phospholipids and using 80% to correct the final results.

The resulting ash is digested using 0.5 mL of 0.2 M HCl for 30 minutes at 90°C. After digestion, the resulting orthophosphate is detected by mixing the sample with a combination of molybdate and ascorbic acid to produce molybdenum blue, as described in Chen's work (1956).

### CITATION

P.S. Chen, T.Y. Toribara and Huber Warner. Microdetermination of Phosphorus. Anal. Chem..

LINK

<https://doi.org/10.1021/ac60119a033>


## CITATION

Ying-Yu Hu, Andrew J. Irwin, Zoe V. Finkel (2022). Improving quantification of particulate phosphorus. *Limnology and Oceanography: Methods*.

LINK

<https://doi.org/10.1002/lom3.10517>

## Prepare phospholipids sample

- 1 Dry remaining organic phase extract of total lipids at  37 °C under a stream of N<sub>2</sub> gas (<2 psi)



## Phosphate primary standard

2h

- 2 KH<sub>2</sub>PO<sub>4</sub> primary standard stock solution (≈ 1 mM)

 Potassium dihydrogen orthophosphate **ACP Chemicals Catalog #P-4550**

- 2.1 Transfer about 1 g KH<sub>2</sub>PO<sub>4</sub> into a beaker, cover the beaker with foil

- 2.2 Place the beaker into an oven, dry KH<sub>2</sub>PO<sub>4</sub> at  110 °C for at least  02:00:00

2h

- 2.3 Move KH<sub>2</sub>PO<sub>4</sub> into a vacuum desiccator, allow KH<sub>2</sub>PO<sub>4</sub> to cool to room temperature

- 2.4 Dissolve around  $\text{0.136 g}$  dried  $\text{KH}_2\text{PO}_4$  in  $\text{1 L}$  MilliQ water.
- Use 1 L volumetric flask
  - Take notes of the actual weight of  $\text{KH}_2\text{PO}_4$  for final concentration of standard stock solution
- 2.5 Transfer standard stock solution into a 1 L bottle and store in the fridge.

#### Note

This stock solution lasts quite a long time, unless there is evidence for growth of algae or other extraneous biotic material.

## High temperature dry combustion

9h

- 3 Use diamond pen to engrave the sample vials with numbers. Log number and sample code.

- 4  $\text{0.17 M}$   $\text{MgSO}_4$  reagent:

Dissolve  $\text{1.023 g}$   $\text{MgSO}_4$  in 50 mL MilliQ water

 Magnesium sulfate anhydrous **Fisher Scientific Catalog #M65500**

- 5 Add  $\text{200 } \mu\text{L}$   $\text{0.17 M}$   $\text{MgSO}_4$  to the dry extract.

#### Note

Sing-use pipet tip to avoid cross-contamination.

- 6 Cover the uncapped vials with foil and place in the oven at  $\text{90 } ^\circ\text{C}$  until samples are completely dry.

### Equipment

**Forced air oven**

NAME

VWR



BRAND

89511-410

SKU

### Note

Remove samples out of the oven as soon as they are dried. If muffle furnace is not available, keep samples in vacuum desiccator.

7 Combust dried samples at  650 °C for  09:00:00

9h

### Equipment

**Muffle furnace**

NAME

F30428C

TYPE

Thermo

BRAND

10-505-13

SKU

### Note

Only place glass vials in the muffle furnace.  650 °C turns foil into ash.

8 Allow samples to gradually cool down in the muffle furnace.

## Digestion

9 [M] 0.2 M HCl reagent:

In a reagent bottle, dissolve one part of [M] 12 N HCl in 59 parts of MilliQ water

⚗ 12 N Hydrochloric acid **Contributed by users**

### Note

Volume of HCl\_0.2M\_mL = (0.5\_mL) X (#Sample + #Blank)

10 Preheat oven to 🌡 90 °C

11 Add 🧪 0.5 mL [M] 0.2 M HCl to each vial.

12 Tightly cap the vial and vortex.

13 Place vials in the oven for ⌚ 00:30:00

30m

14 Cool samples down to 🌡 Room temperature

15 Preheat shaker/incubator to  37 °C

## Equipment

**SHAKING INCUBATOR**

NAME

71L

TYPE

Corning® LSE™

BRAND


6753

SKU

16 Standard working solutions and reagents can be prepared during sample digestion.

17 Standard working solution



	Standard	Primary (uL)	MilliQ (uL)
	S1	0	1000
	S2	5	995
	S3	10	990
	S4	20	980
	S5	50	950
	S6	100	900
	S7	150	850
	S8	200	800

18 Transfer  500 µL of each standard working solution to 2 mL microtube.

## Preparing working reagents



19 All reagents are freshly prepared before colorimetric measurement.

20 [M] 2.5 % ammonium molybdate reagent:

Weigh  0.25 g ammonium molybdate in a Falcon tube and top to  10 g with MilliQ water. Cap and shake until totally dissolved.

 Ammonium molybdate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #09878-100G**

21 [M] 10 % ascorbic acid reagent (avoid light exposure):

Weigh  1 g ascorbic acid in a Falcon tube and top to  10 g with MilliQ water; Cap and shake until all dissolved.

 Ascorbic acid **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A5960-100G**

22 [M] 6 N (3 M) sulfuric acid reagent:

Carefully add 1 part [M] 18 M concentrated sulfuric acid into 5 part MilliQ water

 18M sulfuric acid **Contributed by users**

23 Calculate the volume of molybdate-ascorbic reagent:


Total volume of reagent\_mL = (0.5 mL) X (#standard working solution + #samples + #blanks)

24 Mix the reagents into Falcon tube:

Reagent	Parts as in volume
MilliQ	2
6N sulphuric acid	1
2.5% ammonium molybdate	1
10% ascorbic acid	1

## Colorimetric measurement

3h

- 25 Add  500 µL reagent to each standard, sample (in the vial) and blank, starting from blanks, including blank for standards and blank for samples.

### Equipment

Finntip Stepper Tips

NAME

5 mL

TYPE

Thermo Scientific

BRAND

9404200

SKU

### Note

Before dispensing the reagent, wipe or dab the liquid drop on the outside of the tip, avoid wiping the open tip.

- 26 Vortex.

- 27 Incubate at  37 °C for  03:00:00 while shaking at 150 rpm

3h

- 28 Load microplate with 250 ul reactant from each tube, duplicate.



	<u>1</u>	<u>2</u>	<u>3</u>	<u>4</u>	<u>5</u>	<u>6</u>	<u>7</u>	<u>8</u>	<u>9</u>	<u>10</u>	<u>11</u>	<u>12</u>
<u>A</u>	S1	S1	Samples and sample blanks: 40 with duplicate									
<u>B</u>	S2	S2										
<u>C</u>	S3	S3										
<u>D</u>	S4	S4										
<u>E</u>	S5	S5										
<u>F</u>	S6	S6										
<u>G</u>	S7	S7										
<u>H</u>	S8	S8										

Example of loading the microplate

## 29 Read plate in microplate reader

A	B
Shake duration	00:00:05
Shaking type	Continuous
Shaking force	High
Shaking speed [rpm]	600
Wavelength [nm]	820
Use transmittance	No
Pathlength correction	No
Measurement Time [ms]	100

### Equipment

Varioskan LUX Multimode Microplate Reader

NAME

Thermo Fisher

BRAND

VL0L00D0


SKU

## Calculation

3h

## 30 Subtract the average absorbance at 820 nm of the blank standard replicates from the absorbance

at 820 nm of all other standard working solutions.

- 31** Subtract the average absorbance at 820 nm of the blank sample (i.e. blank filter) replicates from the absorbance at 820 nm of all other individual samples.
- 32** Prepare a standard curve by plotting the average blank-corrected 820 nm absorbance for each standard working solution versus its concentration in uM.  
Molar Mass of KH<sub>2</sub>PO<sub>4</sub>: 136.086 g/mol
- 33** Use the standard curve to determine the orthophosphate concentration of each unknown sample by using its blank-corrected 820 nm absorbance.
- 34**  $(P_{\text{measured}})_{\text{umol/sample}} = (\text{orthophosphate})_{\text{uM}} \times (V_{\text{HCl}})_{\text{mL}} \times (0.001)$   
 $(P_{\text{corrected}})_{\text{umol/sample}} = (P_{\text{measured}}) / 0.8$   
Where, 0.8 is the average recovery of phospholipids after a high temperature dry combustion at  650 °C .
- 35**  $(\text{Phospholipids})_{\text{ug/sample}} = (P_{\text{corrected}}) \times 30.97 / (0.01 \times 4.3)$