

Sep 08, 2020

# RNAPath Target Identification via Fluorescent Hybridization

Devin Willis<sup>1</sup><sup>1</sup>Scooter Willis

In Development

[dx.doi.org/10.17504/protocols.io.bk37kryn](https://dx.doi.org/10.17504/protocols.io.bk37kryn)

XPRIZE Rapid Covid Testing

RNAPath



Devin Willis

## EXTERNAL LINK

<http://rnath.com/>

## DOI

[dx.doi.org/10.17504/protocols.io.bk37kryn](https://dx.doi.org/10.17504/protocols.io.bk37kryn)

## PROTOCOL CITATION

Devin Willis 2020. RNAPath Target Identification via Fluorescent Hybridization. **protocols.io**<https://dx.doi.org/10.17504/protocols.io.bk37kryn>

## EXTERNAL LINK

<http://rnath.com/>

## KEYWORDS

Molecular Instruments, COVID, Microscopy, Hybridization, Fluorescence, RNA

## LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Sep 08, 2020

## LAST MODIFIED

Sep 08, 2020

## PROTOCOL INTEGER ID

41823

## MATERIALS

NAME	CATALOG #	VENDOR
SSC (20X), RNase-free	AM9770	Thermo Fisher
Molecular Instruments HCR COVID Probe Set		Molecular Instruments
Molecular Instruments Alexa Flour 488 B1 HCR Amplifier		Molecular Instruments

## STEPS MATERIALS

NAME	CATALOG #	VENDOR
Molecular Instruments HCR COVID Probe Set		Molecular Instruments
TE Buffer	12090015	Thermo Fisher


NAME	CATALOG #	VENDOR
SSC Buffer, 20X, 1L	V4261	Promega

#### BEFORE STARTING

Ensure you have cleaned the workspace with an appropriate RNase away solution to prevent RNA degradation

### Sample Preparation

- 1 Suspend RNA pellet in  **2 µl TE Buffer**




**TE Buffer**

by Thermo Fisher

Catalog #: 12090015

### Hybridization







- 2 Add  **2 µl Molecular Instruments COVID Probe Set** to the RNA suspension




**Molecular Instruments HCR COVID Probe Set**

by Molecular Instruments

[View](#)

- 3 Place Molecular Instruments Alexa Flour 488 B1 HCR Amplifier in  **90 °C** for  **00:02:00**
- 4 Remove Molecular Instruments Alexa Flour 488 B1 HCR Amplifiers and place at  **Room temperature** for  **00:30:00**
- 5 Add  **2 µl** of each Molecular Instruments Alexa Flour 488 B1 HCR Amplifier H1 and H2 to the RNA suspension
- 6 Dilute  **100 µl** of 20x SSC Buffer to 5x SSC Buffer



**SSC Buffer, 20X, 1L**

by Promega

Catalog #: V4261

7 Add  4  $\mu$ l of 5x SSC Buffer to the RNA Suspension

8 Incubate at  37 °C for  12:00:00

#### Imaging

9 Remove sample from incubation and pipette  1  $\mu$ l onto a glass slide with a coverslip

10 Image slide under a fluorescent microscope with filter cubes for FITC or Alexa Fluor 488