

NOV 08, 2022



WORKS FOR ME 1

Electroporation Transformation Protocol in S. Cerevisiae

This protocol is published without a DOI.

Stephanie Hood¹

¹University of Arkansas

Stephanie Hood

COMMENTS 0

ABSTRACT

Protocol adapted by Stephanie Hood, 09/2022, from:

An improved yeast transformation method for the generation of very large human antibody libraries

Overview:

This protocol calls for ~800-900 ul of cell pellet and is optimized for 1 ug of plasmid DNA in the lab strain (BY4742).

PROTOCOL CITATION

Stephanie Hood 2022. Electroporation Transformation Protocol in S. Cerevisiae. **protocols.io** https://protocols.io/view/electroporation-transformation-protocol-in-s-cerev-cizquf5w

LICENSE

This is an open access protocol distributed under the terms of the <u>Creative Commons</u>

<u>Attribution License</u>, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



CREATED

Nov 08, 2022

LAST MODIFIED

Nov 08, 2022

PROTOCOL INTEGER ID

72464

MATERIALS TEXT

Reagent List:

- 1X YPD
- 2X YPD
- Sterile Water
- 2M LiOAC
- 1 mM DTT
- 1M Sorbitol
- 2M Sorbitol
- 100 mM CaCl2
- SD Ura-liquid
- SD Ura plates

2-3 Days before starting experiment

2d

1

- Streak out yeast strains for electroporation protocol.
- Depending on how many cultures you are setting up you may need to streak more than 1 plate to get enough isolated colonies. Grow at 30°C for 2 days.

Day before experiment

16h

2

- Grow 5 ml culture(s) of *S. cerevisiae* cells overnight to stationary phase (OD600 to or about 3) in 1XYPD media on a platform shaker at 270 rpm and 30°C.
- Need enough cells to set up 100 ml culture(s).

4h 15m

Day of experiment: Preparation of electrocompetent cells

3 • Turn on spectrophotometer (cuvette)

15

- If you have multiple tubes from the same strain, pool the 5 mL O/N cultures into a sterile 50 mL conical tube
 - Vortex to mix
- 5 Spec the O/N cells
 - 1:50 dilution (20 ul cells: 980 ul 1X YPD)
- 6 Set up the culture(s) in 500 mL flasks.

Example:

100 ml 2X YPD culture - starting OD600 of 2.0

Use C1V1=C2V2 to calculate the amount of O/N cells you need to add to 100 ml for a starting OD600 of _____.

• Grow the inoculated cultures on a platform shaker at 30°C and 270 rpm for 4 hours.

While cultures are incubating: prepare media

- Put sterile water on ice and place at 4°C.
 - ~100 mL of sterile water needed per reaction.
- 9 Prepare **electroporation buffe**r and put in ice at 4°C.

500 mL Electroporation buffer: ~100 mL needed per reaction

- Add 250 mL 2M Sorbitol to a sterile 500 mL bottle
- Add 5 mL 100 mM CaCl2
- Bring volume to 500 mL with sterile water.

After 4 hour incubation

■ Collect yeast cells by centrifugation at **1500 g for 3 minutes** and remove the media.

• 2-50 ml conical tubes - 50 ml each

i protocols.io

4h

39m

■ Prepare conditioning buffer during spins- 0.1 M LiAc/10 mM DTT.

Conditioning Buffer: 20 mL needed per reaction

- Add 50 mL sterile water to a sterile 250 mL sterile bottle
- Add 5 mL 2M Lithium Acetate
- Add 1000 ul 1M DTT
- Bring volume to 100 mL (44 mL sterile water added with a serological pipette)

Put any remaining conditioning buffer in a well labeled 50 mL conical tube for disposal with EHS.

- Wash the cell pellets with 25 mL ice-cold water.
 - 50 ml conical tubes 25 ml each
 - Spin at 1500 g for 3 minutes and remove the media.
- Wash the cell pellets again with 25 mL ice-cold water.
 - 50 ml conical tubes 25 ml each
 - Spin at 1500 g for 3 minutes and remove the media.
- Wash the cell pellet with 25 mL of ice-cold electroporation buffer (1M Sorbitol / 1 mM CaCl2)

3m

- 50 ml conical tubes 25 ml each
- Spin at 1500 g for 3 minutes and remove the media.
- Condition the yeast cells by re-suspending the cell pellet in 10 mL 0.1 M LiAc/10 mM DTT.
 - 50 ml conical tubes 10 ml each
 - Pool the tubes culture into one flask for the outgrowth.
 - 10 ml of conditioning buffer for each tube, then pool in a 125 mL flask.
 - Shake at 270 rpm in a culture flask for 30 minutes at 30 °C.
- Pour cells into a 50 mL conical tube.
 - Collect the conditioned cells by centrifugation (1500g for 3 minutes).
- Wash the cells once with 50 mL ice-cold electroporation buffer.
 - Spin at 1500 g for 3 minutes and remove the media.
- Re-suspended the cell pellet in **100 to 200 μL** electroporation buffer to reach a final volume of 1 mL.
 - Record the volume needed to get to 1 mL

Example: 700 ul cell pellet + 300 ul electroporation buffer = 1 mL

2 electroporation reactions of 400 µl each.

You may need to add more electroporation buffer to reach 1 mL cell mixture.

If only doing 1 transformation, less volume can be added as only 400 ul is needed per transformation.

19 The cells are kept on ice until electroporation.

DNA preparation and electroporation and plating

- Pre-chill BioRad GenePulser cuvette(s) (0.2 cm electrode gap).
 - 1 cuvette per transformation
- Prepare enough outgrowth media for all your samples: 1:1 1M sorbitol: 2X YPD
 - 8 mL per sample aliquoted into 50 ml conical tubes.
- Use 1 ug of plasmid DNA per 400 ul of cells.
 - Add ____ uL of plasmid DNA to the 1 mL of cells = 1 ug per 400 uL.
 - 2.5 ug of plasmid DNA in 1 mL of cells.
- Gently mix 400 μl electrocompetent cells and DNA and transfer to a pre-chilled BioRad GenePulser cuvette (0.2 cm electrode gap).
 - Kept on ice for 5 minutes until electroporation.
- Electroporate the cells at 2.5 kV and 25 μF.
 - Typical time constant ranges from 3.0 to 4.5 milliseconds.
 - 200 Ohms
- Using a sterile pasteur pipette, transfer electroporated cells from each cuvette into 8 mL of 1:1 mix of 1 M sorbitol: 2X YPD media.
- Incubate the 50 mL conical tube containing the 8 mL of recovery media in the 30°C standing incubator for 1 hour.
 - Mix by gentle inversion (~4x) every 15 mins.
 - Use a rack to store 50 mL conical tubes in the incubator. Do not store conical tubes upside down as they tend to leak.

1h

■ Collect cells by centrifugation (1500g for 3 minutes) and resuspend in SD URA- + 2% glucose.

protocols.io

• Use 1 mL SD URA- + 2% glucose media for every 400 μl electroporated cells.

28

- Prepare 10-fold serial dilutions from the cell suspension and plate 100 μl of 10-2, 10-3, and 10-4 dilutions onto selection plates and incubate at 30 °C incubator for 2 days.
- If using an antibiotic marker for selection, let the cells grow on just YPD for 1 days and then replica print onto YPD + antibiotic the following day.

2d

2-3 days after experiment

29

• Count the number of colonies on each plate to determine the CFUs.