



AUG 26, 2023

OPEN ACCESS



DOI:
dx.doi.org/10.17504/protocols.io.kqdg3p7kzl25/v1

Protocol Citation: Carlos Guzman, Sascha Duttke, Camila De Arruda Saldanha, Christopher Benner, Sven Heinz 2023. TSS-MPRA Protocol. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kqdg3p7kzl25/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: May 11, 2022


Last Modified: Aug 26, 2023

TSS-MPRA Protocol

Carlos Guzman¹, Sascha Duttke^{1,2}, Camila De Arruda Saldanha¹, Christopher Benner¹, Sven Heinz¹

¹Department of Medicine, Division of Endocrinology, U.C. San Diego School of Medicine, La Jolla, California, 92093 USA;

²School of Molecular Biosciences, College of Veterinary Medicine, Washington State University, Pullman, WA, USA



Carlos Guzman

ABSTRACT

Cis-regulatory elements can be classified by the shapes of their transcription initiation patterns, which are indicative of distinct regulatory mechanisms. While massively parallel reporter assays (MPRAs) have enabled the functional study of sequence features within regulatory elements on an unprecedented scale, current MPRA approaches focus on quantifying transcript abundance, largely ignoring where transcription starts. This information however, could provide evidence that regulatory mechanisms in the context of the reporter assay resemble those active in the genome. Here we describe a transcription start site-capturing massively parallel reporter assay (TSS-MPRA) that simultaneously measures the location and frequency of transcription initiation. We characterize the degree to which plasmid-based MPRAs recapitulate endogenous initiation patterns ("TSS shapes") and transcription levels and evaluate the effects of increasing insert length and reporter chromatinization on plasmid-derived transcription initiation. Employing a new bioinformatic approach to compare TSS shapes, we find that shorter, episomal constructs most faithfully replicate endogenous initiation patterns and transcription levels. Finally, we illustrate how TSS-MPRA can be used to decode cis-regulatory grammar by assessing the effects of core promoter and transcription factor motif mutations and single nucleotide polymorphisms on transcription initiation. Taken together, TSS-MPRA reveals important caveats to consider when using MPRAs and enables high-resolution analysis of the sequence grammar underlying transcription initiation.

PROTOCOL MATERIALS

⊗ Chloroform:Isoamyl alcohol 24:1 Merck MilliporeSigma (Sigma-Aldrich) Catalog #C0549

In 2 steps

⊗ RNase A Qiagen Catalog #19101 Step 49

⊗ UltraFree MC Column (0.45 µM) Merck Millipore (EMD Millipore) Catalog #UFC30HVN

Step 57

⊗ Novex™ Hi-Density TBE Sample Buffer (5X) Thermo Fisher Catalog #LC6678

Step 50

⊗ SYBR Gold Nucleic Acid Gel Stain Contributed by users Catalog # S-11494

Step 52

⊗ TRIzol™ LS Reagent Thermo Fisher Catalog #10296010

In 2 steps

⊗ Quick CIP New England Biolabs Catalog #M0525 In 2 steps

⊗ Suprase-In RNase Inhibitor Thermofisher Catalog #AM2694 Step 16

⊗ 10% TBE Gel Invitrogen - Thermo Fisher Catalog #EC62752BOX Step 51

⊗ ChIP DNA Clean & Concentrator Zymo Research Catalog #D5205 Step 59

⊗ GlycoBlue Coprecipitant Thermo Fisher Scientific Catalog #AM9515 In 2 steps

⊗ RQ1 RNase-Free DNase, 1,000u Promega Catalog #M6101 Step 16

⊗ Sera-Mag SpeedBeads Carboxylate-Modified Magnetic Particles GE Healthcare Catalog #44152105050350

Step 50

⊗ CutSmart Buffer - 5.0 ml New England Biolabs Catalog #B7204S In 2 steps



⊗ 1.5 mL LoBind tubes Eppendorf Catalog #022431021 Step 54







BEFORE START INSTRUCTIONS







- * Make sure that cells are resuspended in media/PBS before TRIzol LS extraction (250 µL)
- * Wipe down work surface and pipettes with RNase Zap
- * Cool down centrifuge to 4°C











RNA Extraction

58m 20s

- 1 Add  750 μL of  TRIzol™ LS Reagent Emd
Millipore Catalog #10296010 per








 250 μL of sample volume (3x) and pipette up and down 5x (**Optional:** samples can be stored at  4 °C ON or at  -20 °C for up to a year).
- 2 Incubate at RT for  00:05:00 . 5m
- 3 Add  230 μL of  Chloroform:Isoamyl alcohol 24:1 Emd
Millipore Catalog #C0549 and 15s



shake vigorously by hand for  00:00:15 .
- 4 Incubate for  00:05:00 at RT. 5m
- 5 Centrifuge samples for  00:15:00 (13,000 x g, 4°C). 15m
- 6 Transfer the aqueous layer containing the RNA to a new 1.5 mL LoBind tube (~ 400 - 500 μL).
- 7 Add  1 μL of gly  GlycoBlue Coprecipitant Emd
Millipore Catalog #AM9515 to each sample.
- 8 Add 1/10th volume of 3M NaOAc ( 5.5) to each sample.



- 9 Vortex  00:00:05 and add 1x volume of isopropanol to each sample. Mix by inverting 10 times and then spin down briefly. 5s
- 10 Incubate ON at  -20 °C (optionally 20 minutes at -20°C).
- 11 Centrifuge samples for  00:30:00 (MAX x g, 4°C). 30m
- 12 Discard supernatant, wash with  1 mL 75% EtOH. Discard EtOH, quick-spin, remove rest of EtOH.
- 13 Air dry pellet at RT until it is translucent (~  00:03:00). (**Optional:** pellet can be frozen at  -80 °C indefinitely) 3m
- 14 Resuspend in  30 µL TE'T [ 0.05 % volume Tween-20,  0.1 millimolar (mM) EDTA,  10 millimolar (mM) Tris pH 7.5].



Capped MPRA 5' RNA-seq



1d 2h 45m 5s


- 15 Aliquot  15 µL of sample into 1.5 mL epi-tube and incubate at  75 °C for  00:02:00 , then chill on ice for  00:02:00 4m
- 16 Add  35 µL of CIP1 master-mix [ 25.25 µL of ddH2O + 0.05% Tween-20,  5 µL 1h




 CutSmart Buffer - 5.0 ml Emd
 Millipore Catalog #B7204S  0.75 µL




 Superase-In RNase Inhibitor Emd
 Millipore Catalog #AM2694  0.5 µL



 RQ1 RNase-Free DNase, 1,000u Emd
 Millipore Catalog #M6101  2 µL



 Quick CIP Emd
 Millipore Catalog #M0525] to each sample and incubate at  37 °C for


 01:00:00 .



17 Incubate at  75 °C for  00:01:00 , then chill on ice for  00:02:00 3m


18 Add  10 µL CIP 2 master-mix [ 8.5 µL of ddH2O + 0.05% Tween-20,  1 µL 30m



 CutSmart Buffer - 5.0 ml Emd
 Millipore Catalog #B7204S  0.5 µL


 Quick CIP Emd
 Millipore Catalog #M0525] to each sample and incubate at  37 °C for


 00:30:00 .

19 Add  500 µL  TRIzol™ LS Reagent Emd
 Millipore Catalog #10296010 , vortex, and 5m




incubate at RT for  00:05:00 .

20 Add  140 µL of TE'T and  140 µL of

 Chloroform:Isoamyl alcohol 24:1 Emd
 Millipore Catalog #C0549 to each sample.

21 Vortex well, then centrifuge samples for  00:10:00 at 12,000 x g (RT). 10m

22 Transfer the upper layer into new tubes.

- 23 Add  1 μ L of gly  GlycoBlue Coprecipitant Emd
Millipore Catalog #AM9515 to each sample.
- 24 Add 1/10th volume of 3M NaOAc ( 5.5) to each sample.
- 25 Vortex  00:00:05 and add 1x volume of isopropanol to each sample. Mix by inverting 10 times and then spin down briefly. 5s
- 26 Incubate  Overnight at  -20 °C (optionally 20 minutes at -20°C). 10m
- 27 Centrifuge samples for  00:30:00 (MAX x g, 4°C). 30m
- 28 Discard supernatant, wash with  1 mL 75% EtOH. Discard EtOH, quick-spin, remove rest of EtOH.
- 29 Air dry pellet at RT until it is translucent (~  00:03:00). (**Optional:** pellet can be frozen at  -80 °C indefinitely) 3m
- 30 Resuspend RNA in  6 μ L TET [ 0.05 % volume Tween-20,  1 millimolar (mM) EDTA,  10 millimolar (mM) Tris pH 7.5].

31 Denature RNA at 75 °C for 00:02:00 , then chill on ice for 00:02:00 . 4m

32 Add 9 µL of RppH master-mix to each sample, mix **really** well and incubate at 37 °C for 01:00:00 1h



Note

The PEG in the master-mix can be extremely viscous, so you need to be very careful when mixing both the master-mix together and when mixing it with the sample to ensure that it's properly mixed

Note

RppH master-mix (per sample):

3.25 µL ddH2O + 0.1% Tween-20

1.5 µL

T4 RNA Ligase Reaction Buffer - 3.0 ml New England
Biolabs Catalog #B0216L

3 µL 50% PEG8000 New England
Biolabs Catalog #B1004S

0.25 µL Superase-In RNase Inhibitor Thermofisher Catalog #AM2694

1 µL

RNA 5' Pyrophosphohydrolase (RppH) - 200 units New England
Biolabs Catalog #M0356S

33 Add 10 µL 5' Ligation master-mix, mix well and incubate at 21 °C for 02:00:00 or 18h
16 °C for 16:00:00 .

Note

5' Ligation master-mix (per sample):

1 μ L

T4 RNA Ligase Reaction Buffer - 3.0 ml New England
Biolabs Catalog #B0216L

2 μ L 10 mM ATP New England
Biolabs

1 μ L 10 micromolar (μ M) of denatured 5' SR Adapter from

NEBNext Multiplex Small RNA Library Prep Set for Illumina (1-12) - 96 rxns New England
Biolabs Catalog #E7300L

5 μ L 50% PEG8000 New England
Biolabs Catalog #B1004S

1 μ L

T4 RNA Ligase 1 (ssRNA Ligase) - 5,000 units New England
Biolabs Catalog #M0204L

Note

It is **CRITICAL** that you denature the 5' SR Adapter at 75 °C for 00:02:00 prior to creating the master-mix.

34 [go to step #19](#)

35 [go to step #20](#)

36 [go to step #21](#)

37  go to step #22

38  go to step #23

39  go to step #24

40  go to step #25

41  go to step #26

42  go to step #27

43  go to step #28

44  go to step #29

45 Resuspend pellet in 7 μL of Annealing master-mix.

Note

Annealing master-mix (per sample):

1 μL of 10 micromolar (μM) RS2 Primer (5'AGCGGATAACAATTTACACAGGA3')
2 μL of 700 millimolar (mM) KCl
4 μL TET

46 Denature RNA at 75 $^{\circ}\text{C}$ for 00:02:00, then incubate at 56 $^{\circ}\text{C}$ for 00:45:00, then cool on ice. 47m

47 Add 13 μL RT master-mix, incubate at 50 $^{\circ}\text{C}$ for 01:00:00. 1h

Note

RT master-mix (per sample):

7.5 μL ddH₂O + 0.05% Tween-20
2 μL RT_KCl_10x
2 μL 10X DTT New England
Biolabs Catalog #B1034A
1 μL 10 mM dNTPs Life
Technologies Catalog #10297-018
1 μL
ProtoScript II Reverse Transcriptase - 10,000 units New England
Biolabs Catalog #M0368L
0.5 μL
RNase Inhibitor, Murine - 15,000 units New England
Biolabs Catalog #M0314L



Note

PCR master-mix (per sample):

🧴 25 μ L

🧴 LongAmp Taq 2X Master Mix - 100 rxns New England
Biolabs Catalog #M0287S

🧴 2.8 μ L

🧴 Betaine 5M Sigma
Aldrich Catalog #B0300

🧴 0.2 μ L

🧴 100 micromolar (μ M) of blue cap SR Primer from

🧴 NEBNext Multiplex Small RNA Library Prep Set for Illumina (1-12) - 96 rxns New England
Biolabs Catalog #E7300L

🧴 2 μ L





















of 🧴 10 micromolar (μ M) 3' barcode primer (we use TruSeq HT Primers D70x-D7xx)


Note


The blue cap SR primer from the NEBNext Small RNA Library Prep kit comes in 10 μ M concentration. Therefore we order our own primer from IDT in order to use 100 μ M concentration primers. Decrease the amount of Betaine used if using 10 μ M concentrations.


A	B	C
PCR Program		
	94C	30 seconds
15x cycles	94C	15 seconds
	63C	30 seconds
	70C	18 seconds***
	70C	5 minutes
	4C	hold infinity

the extension time depends on the length of your inserts! LongAmp copies DNA at 1kb per 50 seconds



- 49 Add  1 μL  RNase A Emd
Millipore Catalog #19101, incubate for  00:15:00 at  37 °C. 15m
- 50 Add  47 μL 20% PEG / 2.5M NaCl +  3 μL  Sera-Mag SpeedBeads Carboxylate-Modified Magnetic Particles Emd
Millipore Catalog #44152105050350 5m
- to each sample. Incubate at RT for  0 μL . Wash 2x with 80% EtOH. Air-dry beads until cracked (~12 minutes). Resuspend in  15 μL of 1x  Novex™ Hi-Density TBE Sample Buffer (5X) Emd
Millipore Catalog #LC6678. Incubate at RT for  00:05:00.
- 51 Run the entire sample on  10% TBE Gel Emd
Millipore Catalog #EC62752BOX with 1h 30m
-  0.5 μL of any 25bp ladder (80V for  00:15:00, then 180V for  01:15:00). Perforate 0.5 mL qubit tubes 4x with 22g needle, and place qubit tubes inside 1.5 mL LoBind tube.
- 52 Stain the gel with  20 mL of 1x TBE. Buffer +  3 μL of 5m
-  SYBR Gold Nucleic Acid Gel Stain Emd
Millipore Catalog # S-11494 for  00:05:00 inside a RNase free container in the dark.
- 53 Visualize gel and cut gel depending on desired size (adapters are ~118 bp, add to size you expect).
- #### Note
- This is highly dependent on the adapters you chose to use and the design of your inserts. In our case our adapters add up to 118 bp and we expect an average transcript length of 60 bp, so we are looking for bands in the ~170-180 bp range.
- 54 Centrifuge samples for  00:05:00 (20,000 x g , RT). Ensure that all the gel is crushed and 5m

spun down into  1.5 mL LoBind tubes Emd
Millipore Catalog #022431021



55 Add  160 μL of TET to each.

56 Shake samples for  00:45:00 at RT (do not shake too quickly).



45m

57 Transfer slurry to  UltraFree MC Column (0.45 μM) Emd
Millipore Catalog #UFC30HVN and
spin for  00:02:00 at 1000 x g.

2m

58 Add 100 μL of TET + 300 mM NaCl to column and incubate for  00:30:00 at RT, then spin for  00:02:00 at 1000 x g.

32m

59 Clean up in  ChIP DNA Clean & Concentrator Emd
Millipore Catalog #D5205 columns
according to instructions. Elute in  10 μL EB.