




Oct 07, 2022

Electroporatic transformation of *Rhodobacter sphaeroides*

Jaya K Yakha¹, Deborah K Hanson¹, Rosemarie Wilton¹, laible¹¹Argonne National Lab

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.eq2ly7d6wlx9/v1 Jaya K Yakha

ABSTRACT

The expression host for *R. sphaeroides* *DrshI* was constructed by deleting Type II restriction enzyme *RshI* (locus tag RSP_3759; recognition sequence CGATCG) from strain DD11 is suitable for transformation by either chemical or electrical procedures.

DOI

dx.doi.org/10.17504/protocols.io.eq2ly7d6wlx9/v1

PROTOCOL CITATION

Jaya K Yakha, Deborah K Hanson, Rosemarie Wilton, laible 2022. Electroporatic transformation of *Rhodobacter sphaeroides*. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.eq2ly7d6wlx9/v1>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 12, 2022

LAST MODIFIED

Oct 07, 2022

PROTOCOL INTEGER ID

69895

MATERIALS TEXT

- 15 ml sterile tubes
- 1 ml microcentrifuge tubes
- Sterile glass beads
- Competent cells (Drsh1)
- Ice
- Electroporation tube
- Incubator at 33⁰C
- Media plate with appropriate Antibiotic selection
- ^GYCC liquid medium
- Centrifuge

1

In a microfuge tube on ice, mix  **40 µL** cells with 100 ng of plasmid DNA.


2

Turn on the electroporator and select the voltage to 2500 V with other setting of 25uF and 400 ohm.

3

Transfer the cells-DNA mixture to a chilled  **2 mm** electroporation cuvette on ice.




4

Transfer the cuvette to the pulse chamber and shock to desired number of pulses (only once)^{1m} (if doing multiple pulses, do first pulse and return cuvette to ice for one minute before doing next (and additional) pulses.  **00:01:00** e on ice between pulses.


5

Return the cuvette to ice and add 1ml ^GYCC **immediately**; record the time constant and the voltage achieved for each pulse.

6

Transfer cells/^GYCC from electroporation cuvette to a Falcon 2059 culture tube (15ml).^{4h}
Incubate at  **35 °C** with slow shaking ( **150 rpm**) for a  **04:00:00** outgrowth period.

7

Pellet the cells in microcentrifuge tube and spread onto ^GYCC agar containing appropriate antibiotic. Incubate plates at  **33 °C** for 2-4 days.