

VERSION 2

OCT 19, 2023

OPEN ACCESS



DOI:

[dx.doi.org/10.17504/protocols.io.6qpvr4bezgm/v2](https://dx.doi.org/10.17504/protocols.io.6qpvr4bezgm/v2)


**Protocol Citation:** Michael Meers Lab 2023. Multiple Targets Identified via Tagmentation (MulTI-Tag) v1.1. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.6qpvr4bezgm/v2> Version created by [Michael Meers Lab](#)

**MANUSCRIPT CITATION:** Meers MP, Llagas G, Janssens DH, Codomo CA, Henikoff S. Multifactorial profiling of epigenetic landscapes at single-cell resolution using MulTI-Tag. *Nat Biotechnol.* 2023 May;41(5):708-716. doi: 10.1038/s41587-022-01522-9. Epub 2022 Oct 31. PMID: 36316484; PMCID: PMC10188359.

## Multiple Targets Identified via Tagmentation (MulTI-Tag) v1.1 V.2

Michael Meers  
Lab<sup>1</sup>

<sup>1</sup>Department of Genetics, Washington University School of Medicine in St. Louis



Michael Meers Lab

### ABSTRACT

We introduce a public protocol for Multiple Targets Identified via Tagmentation (MulTI-Tag), a chromatin profiling approach that ascertains the genomic enrichment of multiple chromatin protein targets in the same cellular context. Multi-Tag uses unique, antibody-specific barcodes incorporated into pA-Tn5-loaded adapters to sequentially tagment targets and read out protein identity via barcode sequencing. MulTI-Tag attains profiles of similar quality to CUT&Tag, and crucially it attains exemplary specificity of barcode-derived target identification. MulTI-Tag can be used to profile multiple chromatin targets in individual cells via plate-based combinatorial indexing approaches, and in principle can be used with microfluidic platforms.

### MATERIALS

- Chilling device (e.g. metal heat blocks on ice or cold packs in an ice cooler)
- Pipettors (e.g. Rainin Classic Pipette 1 mL, 200 µL, 20 µL, and 10 µL)
- Disposable tips (e.g. Rainin 1 mL, 200 µL, 20 µL)
- Disposable centrifuge tubes for reagents (15 mL or 50 mL)
- Standard 1.5 ml microfuge tubes
- 0.5 ml maximum recovery PCR tubes (e.g. Fisher cat. no. 14-222-294)
- Frozen native or lightly cross-linked whole cells or nuclei suspension (e.g. human K562 or H1 cells) prepared as described in the Appendix (Steps 52-64, cite: [Kaya-Okur et al. \(2020\) Nature Protocols](#)).
- Concanavalin A (ConA)-coated magnetic beads (Bangs Laboratories, ca. no. BP531)
- Strong magnet stand (e.g. Miltenyi Macsimag separator, cat. no. 130-092-168)
- Vortex mixer (e.g. VWR Vortex Genie)
- Mini-centrifuge (e.g. VWR Model V)
- PCR thermocycler (e.g. BioRad/MJ PTC-200)
- Distilled, deionized or RNase-free H<sub>2</sub>O (dH<sub>2</sub>O e.g., Promega, cat. no. P1197)
- 1 M Hydroxyethyl piperazineethanesulfonic acid pH 7.9 (HEPES (K+); Sigma-Aldrich, cat. no. H3375)

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Oct 19, 2023

**Last Modified:** Oct 19, 2023

**PROTOCOL integer ID:**  
89627

- 1 M Manganese Chloride (MnCl<sub>2</sub>; Sigma-Aldrich, cat. no. 203734)
- 1 M Calcium Chloride (CaCl<sub>2</sub>; Fisher, cat. no. BP510)
- 1 M Potassium Chloride (KCl; Sigma-Aldrich, cat. no. P3911)
- Roche Complete Protease Inhibitor EDTA-Free tablets (Sigma-Aldrich, cat. no. 5056489001)
- 1 M Hydroxyethyl piperazineethanesulfonic acid pH 7.5 (HEPES (Na<sup>+</sup>); Sigma-Aldrich, cat. no. H3375)
- 5 M Sodium chloride (NaCl; Sigma-Aldrich, cat. no. S5150-1L)
- 2 M Spermidine (Sigma-Aldrich, cat. no. S0266)
- 10% Triton X-100 (Sigma-Aldrich, cat. no. X-100)
- 0.5 M Ethylenediaminetetraacetic acid (EDTA; Research Organics, cat. no. 3002E)
- Antibody to an epitope of interest. Because in situ binding conditions are more like those for immunofluorescence (IF) than those for ChIP, we suggest choosing IF-tested antibodies if CUT&RUN/Tag-tested antibodies are not available
- CUTAC control antibody to RNA Polymerase II Phospho-Rpb1 CTD Serine-5 phosphate (PolII S5P) or histone H3K4me<sub>2</sub>. We have obtained excellent results with these rabbit monoclonal antibodies: Cell Signalling Technology Phospho-Rpb1 CTD (Ser5), CST #13523 (D9N5I) and Epicypher H3K4me<sub>2</sub> #13-0027.
- Secondary antibody, e.g. guinea pig α-rabbit antibody (Antibodies online cat. no. ABIN101961) or rabbit α-mouse antibody (Abcam cat. no. ab46540)
- Protein A/G–Tn5 (pAG–Tn5) fusion protein loaded with double-stranded adapters with 19mer Tn5 mosaic ends (Epicypher cat. no. 15-1117).
- 1 M Magnesium Chloride (MgCl<sub>2</sub>; Sigma-Aldrich, cat. no. M8266-100G)
- 1 M [tris(hydroxymethyl)methylamino]propanesulfonic acid (TAPS) pH 8.5 (with NaOH)
- 1,6-hexanediol (Sigma-Aldrich cat. no. 240117-50G)
- N,N-dimethylformamide (Sigma-Aldrich cat. no. D-8654-250mL)
- NEBNext 2X PCR Master mix (New England Biolabs cat. no. ME541L)
- PCR primers: 10 μM stock solutions of i5 and i7 primers with unique barcodes [Buenrostro, J.D. et al. Nature 523:486 (2015)] in 10 mM Tris pH 8. Standard salt-free primers may be used. We do not recommend Nextera or NEBNext primers.
- 10% Sodium dodecyl sulfate (SDS; Sigma-Aldrich, cat. no. L4509)
- Thermolabile Proteinase K (New England Biolabs, cat. no. P8111S)
- SPRI paramagnetic beads (e.g. HighPrep PCR Cleanup Magbio Genomics cat. no. AC-60500)
- 10 mM Tris-HCl pH 8.0
- Ethanol (Decon Labs, cat. no. 2716)

## Reagent setup

- 1 **Bead binding buffer** Mix 200 μL 1M HEPES-KOH pH 7.9\*, 100 μL 1M KCl, 10 μL 1M CaCl<sub>2</sub> and 10 μL 1M MnCl<sub>2</sub>, and bring the final volume to 10 mL with dH<sub>2</sub>O. Store the buffer at 4 °C for up to several months. \*HEPES-NaOH pH 7.5 is OK.

**Wash buffer** Mix 1 mL 1 M HEPES pH 7.5, 1.5 mL 5 M NaCl, 12.5 µL 2 M spermidine and 250 µL 10% Triton-X100, bring the final volume to 50 mL with dH<sub>2</sub>O, and add 1 Roche Complete Protease Inhibitor EDTA-Free tablet. Store the buffer at 4 °C for up to 2 days.

**Wash 300 buffer** Add 750 µL 5M NaCl to 24.3 ml Triton-Wash buffer. Store at 4 °C for up to 2 days.

**TAPS wash buffer** Mix 1 mL dH<sub>2</sub>O, 10 µL 1 M TAPS pH 8.5, 0.4 µL 0.5 M EDTA (10 mM TAPS, 0.2 mM EDTA).

**TAPS-SDS buffer** (For 32 samples) Mix 20 µL 10% SDS and 2 µL 1 M TAPS pH 8.5 in 178 µL dH<sub>2</sub>O. Just before use add 20 µL Thermolabile Proteinase K (NEB cat. no. P8111S).

**Triton neutralization solution** Mix 30 µL 12.5% Triton X-100 with 370 µL ddH<sub>2</sub>O.

## Conjugate generation (4 hours)

- 2 Resuspend dehydrated, 5'-aminated (NH<sub>2</sub>) P5\_i5 oligo in 1xPBS at a concentration of 200 µM.
- 3 Resuspend dehydrated Tn5MErev oligo (5'-[phos]CTGTCTCTTATACACATCT-3') in 1xPBS at a concentration of 200 µM.
- 4 Anneal P5\_i5 oligo with an equal volume of Tn5MErev oligo to generate 100 µM P5\_i5 adapter by incubating the mixture at 95°C for 2 minutes on a thermal cycler, then ramping down the temperature by 0.5°C every 30 seconds until it reaches 25°C.
- 5 Use annealed P5\_i5 adapter, along with 100 µg primary antibody purified in 1xPBS in the absence of BSA, glycerol, and sodium azide, to generate antibody-adapter conjugates at a molar ratio of 1:10 according to manufacturer's protocols (Abcam ab218260).

## Conjugate assembly (1 hour)

- 6 Resuspend free P5\_i5 oligo in 1xTE (10 mM Tris-HCl pH 8.0, 0.1 mM EDTA) at a concentration of 200  $\mu$ M.
- 7 For bulk experiments: Resuspend free P7\_i7 oligo in 1xTE at a concentration of 200  $\mu$ M. For single cell experiments: Resuspend 96 uniquely barcoded P7\_i7 oligos in 1xTE at a concentration of 200  $\mu$ M.
- 8 Resuspend dehydrated Tn5Merev oligo (5'-[phos]CTGTCTCTTATACACATCT-3') in 1xTE at a concentration of 200  $\mu$ M.
- 9 Anneal P5\_i5 oligo with an equal volume of Tn5Merev oligo to generate free 100  $\mu$ M P5\_i5 adapter as described in step 4.
- 10 **For bulk experiments:** Anneal P7\_i7 oligo with an equal volume of Tn5Merev oligo to generate free P7\_i7 adapter as described in step 4. **For single cell experiments:** In each well of a 96 well plate, anneal one of 96 uniquely barcoded P7\_i7 oligos with an equal volume of Tn5Merev oligo.
- 11 For each conjugate to be used, combine 0.5  $\mu$ g antibody-i5 adapter conjugate with 1  $\mu$ L 5  $\mu$ M pA-Tn5 and 16 pmol free P5\_i5 adapter in a minimal volume and incubate 30-50 minutes at RT.
- 12 **For bulk experiments:** Combine 1  $\mu$ L pA-Tn5 with 32 pmol P7\_i7 adapter in a minimal volume and incubate for 30-50 minutes at RT. **For single cell experiments:** In each well of a 96 well plate, combine 10  $\mu$ L pA-Tn5 with 320 pmol of one of 96 uniquely barcoded P7\_i7 adapters in a minimal volume and incubate for 30-50 minutes at RT. Seal the plate, store at 4°C, and use for several future experiments.
- 13 Use assembled conjugates within 24 hours—when incubating overnight, conjugates should be assembled directly before incubation.

## MulTI-Tag (1+ days)

- 14**     **For bulk experiments:** Wash 10  $\mu$ L paramagnetic Concanavalin A beads (Bangs Laboratories) 2x with 1 mL Bead Binding Buffer, then resuspend in 10  $\mu$ L Bead Binding Buffer. **For single cell experiments:** Wash 5  $\mu$ L paramagnetic Wheat Germ Agglutinin (WGA)-MyOneC1 Dynabeads 2x with 1 mL Bead Binding Buffer, then resuspend in 5  $\mu$ L Bead Binding Buffer. While this quantity of beads is sufficient for 50000 nuclei, combinatorial indexing experiments should prepare to use an input of 2 million nuclei that will be distributed equally across a 96 well plate (~20833 nuclei/well) For example: you might prepare 8 samples of 250000 nuclei and 25  $\mu$ L beads each, which will then be processed normally until step 26, at which point each sample will be aliquoted equally across 12 wells in the 96 well plate.
- 15**     Add resuspended beads to 50000 nuclei while gently vortexing to immediately mix.
- 16**     Incubate beads and nuclei with rotation for 10 minutes at RT.
- 17**     Clear supernatant from nuclei-bead mix of magnet stand and resuspend in 50  $\mu$ L Wash 300 buffer plus 2 mM EDTA and assembled conjugate of interest.
- 18**     Incubate 1 hour at RT or overnight at 4°C with nutation. Upon addition of last conjugate, go to Step 24.
- 19**     Clear supernatant and wash 2x with 200  $\mu$ L Wash 300 buffer.
- 20**     Clear supernatant and add 50  $\mu$ L Wash 300 buffer plus 10 mM MgCl<sub>2</sub>.

- 21 Incubate 1 hour at 37°C.
- 22 Clear supernatant on magnet stand and wash 2x with 200 µL Wash 300 buffer plus 2 mM EDTA.
- 23 Repeat steps 17-22 for all conjugates to be added; after addition of the last conjugate, skip from Step 17 to Step 24.
- 24 Clear supernatant and add 100 µL Wash Buffer with 1:100 Secondary antibody.
- 25 Incubate 30 minutes-1 hour at RT with nutation.
- 26 Clear supernatant and Wash 1x with 200 µL Wash Buffer.
- 27 **For bulk experiments:** Clear supernatant and add 100 µL Wash 300 Buffer plus assembled P7\_i7 adapter pA-Tn5. **For single cell experiments:** Clear supernatant and add Wash 300 buffer to samples such that the sum total volume of all samples is 26 µL \* 96 (2.469 mL). Distribute 25 µL to each well in a Lo Bind 96 well plate, and add 2 µL uniquely barcoded P7\_i7 adapter pA-Tn5 mix to each well using an 8-channel p20 micropipette.
- 28 Incubate 1 hour at RT with nutation.
- 29 **For bulk experiments:** Clear supernatant and wash 2x with 200 µL Wash 300 buffer. **For single**

**cell experiments:** Clear supernatant and wash each well in 96 well plate 2x with 100  $\mu$ L Wash 300 buffer.

**30** Clear supernatant and add 50  $\mu$ L Wash 300 buffer plus 10 mM MgCl<sub>2</sub>.

**31** Incubate 1 hour at 37°C.

**32** Clear supernatant and add 100  $\mu$ L TAPS wash buffer.

### **Bulk Multi-Tag DNA preparation (1.5 hours)**

**33** Clear supernatant and add 5  $\mu$ L 10 mM TAPS-SDS buffer.

**34** Incubate 1 hour at 58°C.

**35** Add 16  $\mu$ L Triton neutralization solution and mix well by pipetting.

**36** Add 2  $\mu$ L each of forward and reverse barcoded PCR primers.

**37** Add 25  $\mu$ L NEBNext 2x PCR Master Mix and mix well with pipetting.

## Single cell Multi-Tag DNA preparation (1.5+ hours)

- 38** Clear supernatant from columns 1-6 and, using wide-bore pipette tips and an 8-channel p200 micropipette, add 100  $\mu$ L TAPS wash buffer to column 1.
- 39** Pipette up and down to dislodge all nuclei from the walls, then remove TAPS-nuclei mix and pipette into column 2.
- 40** Repeat step 39 for all columns 2 through 6 to ensure that all nuclei are resuspended and consolidated in column 6.
- 41** Repeat steps 38-40 to clear any remaining nuclei, resulting in 200  $\mu$ L TAPS wash buffer per well in column 6.
- 42** Repeat steps 38-41 for columns 7-12.
- 43** Clear supernatant from column 6 and, using a wide-bore pipette tip, add 100  $\mu$ L TAPS wash buffer to row A column 6.
- 44** Pipette up and down to dislodge all nuclei from the walls, then remove TAPS-nuclei mix and pipette into row B column 6.
- 45** Repeat step 44 for all rows B through H in column 6 to ensure that all nuclei are resuspended and consolidated in Row H column 6.



- 46** Repeat steps 43-45 to clear any remaining nuclei, resulting in 200  $\mu$ L TAPS wash buffer per well in row H column 6.
- 47** Repeat steps 43-46 for column 12.
- 48** Transfer 200  $\mu$ L nuclei mix from each of row H column 6 and row H column 12 into fresh, low-bind 1.7 mL tubes, then chill on 0°C cold block on ice for 2 minutes.
- 49** To gently digest extracellular DNA to reduce nuclei/bead clumping, proceed to step 50. Otherwise, proceed to step 53.
- 50** Prepare RQ1 RNase-free DNase solution (Promega) by combining 10  $\mu$ L 10x DNase buffer, 80  $\mu$ L ddH<sub>2</sub>O, and 10  $\mu$ L RQ1 DNase.
- 51** Clear supernatant and transfer empty tube back to 0°C cold block. Squirt entire 100  $\mu$ L of DNase solution onto beads, flick tube gently to mix, and incubate on cold block for 5-10 minutes.
- 52** Add 100  $\mu$ L DNase stop solution and flick tube gently to mix.
- 53** Immediately transfer to magnet stand, clear supernatant and add 100  $\mu$ L 10 mM TAPS.

- 54** Bring nuclei-bead solution up to 650  $\mu$ L and add to 20-micron filter affixed to fresh Eppendorf tube; spin 5 minutes at 100xg and collect supernatant.
- 55** Prepare nuclei for nanowell/microfluidic platform of choice according to manufacturer's instructions.

### **Amplification and cleanup (1.5 hours)**

- 56** Amplify DNA using the following conditions:
1. 5 minutes at 58°C
  2. 5 minutes at 72°C
  3. 45 seconds at 98°C
  4. 15 seconds at 98°C
  5. 10 seconds at 60°C
  6. Repeat steps 4-5 13 times (14 total cycles)
  7. 5 minutes at 72°C
  8. Hold at 8°C
- 57** Purify DNA once using 1.1x Ampure XP magnetic beads and resuspend in 25  $\mu$ L 10 mM Tris-HCl pH 8.0 plus 1 mM EDTA.