

Aug 22, 2024 Version 2

Splitting Adherent Cell Lines V.2

DOI

dx.doi.org/10.17504/protocols.io.x54v9p51qg3e/v2

Carolina Lopez¹

¹Washington University

Lopez Lab



Carolina Lopez

Washington University

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.x54v9p51qg3e/v2

Protocol Citation: Carolina Lopez 2024. Splitting Adherent Cell Lines. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.x54v9p51qg3e/v2> Version created by **Sydney Faber**

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 17, 2023

Last Modified: August 22, 2024

Protocol Integer ID: 106114

Disclaimer

Timing for trypsin treatment and following splits may need to be adjusted based on the different cell types considering how adhesive they are and their growth rate.

Abstract

This protocol describes how to split and maintain adherent cell lines in culture. Examples of these cells are: A549 cells, LLCMK2 cells, and MDCK cells.

Materials

TISSUE CULTURE MEDIUM (TCM): Filter through 0.2 µm filter

Component	Amount	Conc. Supp.	Product information
DMEM	500 mL		Gibco Cat. # 11965092 (#11965118-cs)
Gentamicin	500 µL	50µg/mL	Gibco Cat #15750060 (#15750078-pk) – 50 mg/mL
Sodium Pyruvate	5.0 mL	1mM	Corning Cat #25-000-CI – 100mM
L-Glutamine	5.5 mL	2 mM	Sigma Aldrich Cat #G7513 – 200mM
FBS	50 mL	10%	



Protocol

1 **Splitting Adherent cells**

1. Wash flask 2x with sterile PBS.
2. Add 2mL of trypsin/T75. Incubate at 37°C for 2-3 mins or until cells are detached from the flask.
3. Add 8mL of cell culture media and pipette up and down. Transfer all the media to a 15mL tube.
4. Centrifuge at 1200 rpm for 5 min.
5. Discard media and resuspend cells in 10mL Tissue Culture Medium
6. Add 2mL of cell suspension to T75 flask, then add 8mL Tissue Culture Medium. The cell will be ready for next split two days later.