

Sep 06, 2022


🌐 Thawing of mouse embryonic fibroblasts (MEFs) for hPSC cultures

 In 1 collection

Hanqin Li¹, Oriol Busquets², Steven Poser², Dirk Hockemeyer¹, Frank Soldner²

¹University of California, Berkeley; ²Albert Einstein College of Medicine

1 Works for me

 Share

dx.doi.org/10.17504/protocols.io.b4n4qvgw



Devin E Snyder

ABSTRACT

This protocol describes the thawing of mouse embryonic fibroblasts (MEFs) as feeder cells for human pluripotent stem cell (hPSC) culture.

General notes

1. Throughout this protocol, the term hPSC is used to collectively refer to both hiPSCs and hESCs. All described procedures have been tested and work equally well for hiPSCs and hESCs.
2. MEFs were obtained as described in Manipulating the Mouse Embryo: A Laboratory Manual, Third Edition (ISBN: 0879695919)

Andras Nagy, Marina Gertsenstein, Kristina Vintersten, & Richard Behringer. Manipulating the Mouse Embryo: A Laboratory Manual, 3rd ed.. Cold Spring Harbor Laboratory Press.

DOI

dx.doi.org/10.17504/protocols.io.b4n4qvgw

PROTOCOL CITATION

Hanqin Li, Oriol Busquets, Steven Poser, Dirk Hockemeyer, Frank Soldner 2022. Thawing of mouse embryonic fibroblasts (MEFs) for hPSC cultures. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.b4n4qvgw>



FUNDERS ACKNOWLEDGEMENT

Aligning Science Across Parkinson's

Grant ID: ASAP-000486

COLLECTIONS ⓘ



Maintenance and inactivation of mouse embryonic fibroblasts (MEFs) as feeder cells for human pluripotent stem cell culture

KEYWORDS

ASAPCRN

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Feb 03, 2022

LAST MODIFIED

Sep 06, 2022

PROTOCOL INTEGER ID

57788

PARENT PROTOCOLS

Part of collection

[Maintenance and inactivation of mouse embryonic fibroblasts \(MEFs\) as feeder cells for human pluripotent stem cell culture](#)

MATERIALS TEXT

Item	Vendor	Catalog #
DMEM	Corning	10-013-CV
FB Essence	Avantor	10803-034
FBS	Gibco	10437028
200mM L-Glutamine	Sigma	G8540
Penicillin & Streptomycin	Gibco	15140-122
MEM Non-Essential Amino Acids	Gibco	11140-050
Gelatin powder	Sigma	G2625

BEFORE STARTING

All cell culture plates which are used as feeders to maintain hPSCs are coated for at least 1 hour with autoclaved 0.2% gelatin solution at room temperature. Remove gelatin solution immediately before plating MEF cells.

0.2% Gelatin Solution

A	B
Sterile H ₂ O	1L
Gelatin powder	2g

After preparation, the gelatin solution should be autoclaved. Final volume: 1L

- 1 To recover frozen stocks for MEF expansion, P0 MEF tubes will be thawed in a water bath at **37 °C** by gently shaking
- 2 Thawed cells are transferred into a 15 ml conical tube containing 9 ml pre-warmed MEF medium. Centrifuge the tube at **250 x g, 00:05:00** 5m

2.1 MEF medium

A	B
DMEM	435 ml
FB Essence/FBS*	75 ml
200mM L-Glutamine	5 ml
Penicillin & Streptomycin (100x)	5 ml
MEM Non-Essential Amino Acids	5 ml

*We have successfully used either FB Essence or FBS and have not observed an obvious difference. Final volume: 500ml

- 3 Resuspended MEFs (10x10⁶ cells/10 cm plate) are plated in fresh MEF medium and maintained in a humidified incubator (37°C; 5% CO₂)
- 4 3-4 days after thawing, the cells should be ready for passaging (90-95% confluent)