

Jul 24, 2024

MSD + S + AS

This protocol is a draft, published without a DOI.

Eleanor Harman¹

¹UGA



Eleanor Harman

UGA

OPEN  ACCESS



Protocol Citation: Eleanor Harman 2024. MSD + S + AS. protocols.io <https://protocols.io/view/msd-s-as-dhrj354n>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 23, 2024

Last Modified: July 24, 2024

Protocol Integer ID: 103947

Abstract

.



Media

- 1 MS Media: 4.4g/L 2.2g/500mL
Sucrose: 30g/L 15g/500mL
Sorbitol(5%): 50g/L 25g/500mL
Thymidine: 300mg/L 150mg/500mL
2,4D: 2mg/L 1mg/500mL
Phytigel(1.4%): 14g/L 7g/500mL
- 2 Add 80% final volume DI water to beaker
- 3 Add MS Media, Sucrose, Sorbitol, Thymidine, and 2,4D until fully dissolved
- 4 Adjust pH to 5.8
- 5 Add DI water until final volume
- 6 Autoclave on 20 minute liquid cycle