

FEB 27, 2023

## Plant Extraction and Fractionation

Matin Mahmood<sup>1</sup>, Enas Jawad Kadhim<sup>2</sup>,  
Abdulkareem Hameed Abd<sup>3</sup>

<sup>1</sup>Department of Pharmacology, College of Pharmacy, Alkitab University, Altun Kopre, Kirkuk, Iraq;

<sup>2</sup>Department of Pharmacognosy and Medicinal Plants, College of Pharmacy, University of Baghdad, Baghdad, Iraq;

<sup>3</sup>Department of Pharmacology, College of Medicine, Al-Nahrain University, Kadhimiya, Baghdad, Iraq

OPEN ACCESS

DOI:  
[dx.doi.org/10.17504/protocols.io.q26g7y6bqgwz/v1](https://dx.doi.org/10.17504/protocols.io.q26g7y6bqgwz/v1)

**Protocol Citation:** Matin Mahmood, Enas Jawad Kadhim, Abdulkareem Hameed Abd 2023. Plant Extraction and Fractionation. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g7y6bqgwz/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited


**Protocol status:** Working  
We use this protocol and it's working

**Created:** Feb 18, 2023

**Last Modified:** Feb 27, 2023

**PROTOCOL integer ID:**  
77255

**Keywords:** Plant extraction, Soxhlet, Hot method, Fractionation



Matin Mahmood

Department of Pharmacology, College of Pharmacy, Alkitab Uni...

### ABSTRACT

Plant extraction is **a process that aims to extract certain components present in plants**. It is a solid/liquid separation operation: a solid object (the plant) is placed in contact with a fluid (the solvent). The plant components of interest are then solubilized and contained within the solvent

### MATERIALS

500g of granular powdered plant

[hexane](#) (BDH chemicals, England cat-no. BDH24575.100E) (in a ratio of (1:3) W/V

[Soxhlet apparatus](#)(BOECO, Germany)

[80% ethanol](#) (Sigma-Aldrich, Germany cat-no1070172511)

IKA RV 10 Rotary Evaporator (Germany)

5% [HCl](#) (Sigma-Aldrich cat-no. 1009861000)

[ethyl acetate](#) (Sigma-Aldrich, Germany cat-no. 319902-1L)

[NaOH](#) (Honeywell, USA cat-no. 30620)

[chloroform](#)(Honeywell, USA cat-no. C2432)

separatory funnel


[methanol](#) (Biochem Chemopharma, France cat-no. 213032500)

[petroleum ether](#) (Sigma-Aldrich, Germany cat-no. 32299)

- 1 A (500gm) of shade-dried for 12 days.
- 2 Coarsely powdered plant materials
- 3 Powdered plant materials were defatted with hexane in ratio 1:3 W/V for 24 hr. ( 500 gm of plant material with 1500 mL of hexane)
- 4 Allowed to dry at room temperature
- 5 The defatted plant materials were extracted with (2 Liters) of 80% ethanol in soxhlet apparatus until complete exhaustion.
- 6 The alcoholic extract was evaporated to dryness, under reduced pressure at a temperature not exceeding 40 °C by using IKA RV 10 Rotary Evaporator to give a dark greenish-yellow residue designated as a crude fraction (A)
- 7 Crude fraction was acidified with 300ml of 5% hydrochloric acid to pH 2
- 8 Crude fraction then partitioned (three times) with equal volume of ethyl acetate to get two layers (aqueous acidic and ethyl acetate layer), this step is necessary to get rid from any basic compound found in the crude extract.
- 9 The ethyl acetate layer of the original alcoholic extract (crude fraction) was evaporated to dryness under reduced pressure

- 10 Basified with 300ml of 5% sodium hydroxide to pH 10 and extracted with chloroform in the separatory funnel to get two layers
- 11 The aqueous basic layer which was separated, evaporated to dryness and acidified with 5% HCL to pH 2 then extracted with ethyl acetate to get fraction designated as fraction 3 (F-B)
- 12 The second layer-chloroform layer which was also separated and evaporated to dryness under reduced pressure then partitioned with 80% methanol and petroleum ether to get two layers petroleum ether fraction (C) and 80% methanol fraction which designated as fraction 4 (F-D)

## Protocol Scheme

- 13  scheme plant extraction.pdf