



Nov 02, 2021

# Live, in vivo, imaging (GCaMP6F), Video Processing, and Analysis

Bryan Yoo<sup>1</sup>, Jessica Griffiths<sup>1</sup>, Sarkis Mazmanian<sup>1</sup><sup>1</sup>California Institute of Technology

1

[dx.doi.org/10.17504/protocols.io.bzp7p5rn](https://dx.doi.org/10.17504/protocols.io.bzp7p5rn)

Mazmanian Lab

Jessica Griffiths

Protocol for GCaMP6F imaging and analysis in gut neurons used in Yoo et al 2021

DOI

[dx.doi.org/10.17504/protocols.io.bzp7p5rn](https://dx.doi.org/10.17504/protocols.io.bzp7p5rn)

Bryan Yoo, Jessica Griffiths, Sarkis Mazmanian 2021. Live, in vivo, imaging (GCaMP6F), Video Processing, and Analysis. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bzp7p5rn>



Aligning Science Across Parkinson's

Grant ID: ASAP-000375

Department of Defense

Grant ID: PD160030

National Institutes of Health

Grant ID: GM007616 and DK078938

Center for Environmental and Microbial Interactions (CEMI)

Grant ID: n/a

Heritage Medical Research Institute

Grant ID: n/a

Emerald Foundation

Grant ID: n/a

protocol ,

Nov 02, 2021

## Preparation


- 1 AAV-PHP.S-CAG-GCaMP6F ( $10^{12}$ VGs) was delivered systemically to WT C57.
- 2 3-4 weeks after infection, mice were anesthetized with 2% isoflurane on a heating pad (Kent Scientific, Torrington, CT-DCT15) with plastic sleeve covers (Kent Scientific, Torrington, CT-DCT1520P).
- 3 The abdominal cavity was surgically opened to expose the intestines.

## Imaging

- 4 The proximal colon was identified, and this portion was placed on top of a stack of 4-6 glass microscopy slides (depending on size of animal) (VWR, Radnor PA-Cat. No. 48300-026).
- 5 Tissue was secured onto glass slide with a biorthogonal silicon elastomer (Kwik-Sil) (World Precision Instruments, Sarasota, FL-KWIK-SIL), and a glass coverslip.
- 6 Elastomer stiffens within 1 minute of application and coverslip must be held steadily to ensure a flat imaging surface. Anesthetized mouse is placed under an upright confocal microscope (Zeiss LSM 880).
- 7 Using a 10X objective, GCaMP6F fluorescence was taken at 5Hz (1 image every 200ms).
- 8 Periods of movement of tissue and luminal contents is normal during live imaging.

## Processing & Analysis

- 9 Cell tracking was performed in 2D by using TrackMate ImageJ plugin (<https://imagej.net/TrackMate>) and fluorescence intensity was recorded from cells.

- 
- 10 Average background was determined by taking the fluorescence of a region of interest that did not contain a cell, over the duration of the video. Background was subtracted from cell fluorescence intensities.