

AUG 08, 2023

Mitophagy induction using Oligomycin/Antimycin A

Louise Uoselis¹

¹Lazarou Lab, WEHI



Louise Uoselis WEHI

ABSTRACT

Mitophagy induction in HeLa cells using Oligomycin/Antimycin A.

OPEN ACCESS



Protocol Citation: Louise Uoselis 2023. Mitophagy induction using Oligomycin/Antimycin A. protocols.io https://protocols.io/view/mito phagy-induction-usingoligomycin-antimycin-acybixske

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: Aug 08, 2023

Last Modified: Aug 08,

2023

PROTOCOL integer ID:

86090

Keywords: ASAPCRN

Dav	1

1 Seed cells, aiming for a confluency of 80-90% at the time of treatment the next day.

Day 2

- 2 Feed cells for 1 h in an appropriate volume of standard growth media.
- To start the treatment, replace the media in each well with standard growth media that contains 10 uM Oligomycin, 4 uM Antimycin A, and 10 uM QVD.
- 4 Harvest the samples after the desired treatment times.