

May 14, 2024

Staining of Dissociated Lung Cells for Cellular Indexing of Transcriptomes and Epitopes by Sequencing (CITE sequencing)

DOI

dx.doi.org/10.17504/protocols.io.e6nvw1jk7lmk/v1

Gautam Bandyopadhyay¹, Jeffrey Malik², Heidie Huyck³, Ravi Misra⁴, John_Ashton⁵, Gloria S Pryhuber⁵

¹Department of Pediatrics, University of Rochester Medical Center;

²Genomic Research Center, University of Rochester Medical Center; ³University of Rochester; ⁴URMC;

⁵University of Rochester Medical Center

URMC Pryhuber Lab

Gautam Bandyopadhyay



Gautam Bandyopadhyay

Department of Pediatrics, University of Rochester Medical Ce...

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.e6nvw1jk7lmk/v1**

Protocol Citation: Gautam Bandyopadhyay, Jeffrey Malik, Heidie Huyck, Ravi Misra, John_Ashton, Gloria S Pryhuber 2024. Staining of Dissociated Lung Cells for Cellular Indexing of Transcriptomes and Epitopes by Sequencing (CITE sequencing). **protocols.io** **<https://dx.doi.org/10.17504/protocols.io.e6nvw1jk7lmk/v1>**

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 14, 2024

Last Modified: May 14, 2024

Protocol Integer ID: 99762



Funders Acknowledgement:

**NHLBI Molecular Atlas of
Lung Development Program
Human Tissue Core Grant
Grant ID: U01HL14886**

Abstract


This article describes step by step protocol for antibody-staining frozen lung cells for CITE sequencing.



Thawing of frozen dissociated human lung cells

- 1 Frozen cells were thawed in 37°C water bath and soon after the freezing medium melted cells were transferred into a 15 ml conical tube and the tubes were filled with staining buffer [1% BSA (w/v, Millipore) in DPBS (Biowhittaker)] to dilute the DMSO in freezing medium.
- 2 Cells were centrifuged at 800g, 4°C for 10 minutes.
- 3 After centrifugation supernatants were discarded and cells were resuspended in 1.0 ml staining buffer.

Fc receptor blocking and antibody staining

- 4 Viable cell counts were determined by trypan blue exclusion assay
- 5 Cells were then incubated in staining buffer containing 2% v/v normal mouse serum (Millipore Sigma) for 10 minutes over ice for Fc-receptor blocking.
- 6 Cells were then washed with staining buffer and incubated at a final concentration of 1×10^6 cells per 100 μ l of staining cocktail containing analyte-specific antibodies (**attached table**) for 30 min at 4°C.  Covid Supp CITE Seq Staining Antib... 14KB
- 7 After antibody-staining, cells were washed three times (800g, 4°C for 10 minutes) and resuspended in 100 μ l staining buffer and passed through a 100- μ m strainer (Falcon, Corning) and sent Genomic Research Center, University of Rochester Medical Center for further processing for CITE sequencing.