

## Antioxidant activity by DPPH assay: in vitro protocol

Adrieli Sachett<sup>1</sup>, Matheus Gallas-Lopes<sup>1</sup>, Greicy M M Conterato<sup>2</sup>, Ana P Herrmann<sup>1</sup>, Angelo Piato<sup>1</sup>

<sup>1</sup>Universidade Federal do Rio Grande do Sul; <sup>2</sup>Universidade Federal de Santa Catarina





dx.doi.org/10.17504/protocols.io.btbpnimn



Angelo Piato Universidade Federal do Rio Grande do Sul

## ABSTRACT

Considering the role of oxidative stress in the pathology of several diseases and the use of antioxidants as treatment and/or adjuvants in these conditions. Here we propose a protocol to evaluate the antioxidant capacity of compounds by the DPPH method through the scavenging capacity of free radicals by reducing the DPPH radical. This protocol was standardized at LAPCOM (Psychopharmacology and Behavior Laboratory at UFRGS) to assess biochemical parameters *in vitro*.

DOI

dx.doi.org/10.17504/protocols.io.btbpnimn

PROTOCOL CITATION

Adrieli Sachett, Matheus Gallas-Lopes, Greicy M M Conterato, Ana P Herrmann, Angelo Piato 2021. Antioxidant activity by DPPH assay: in vitro protocol. **protocols.io** https://dx.doi.org/10.17504/protocols.io.btbpnimn

LICENSE

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Mar 13, 2021

LAST MODIFIED

Jun 15, 2021

OWNERSHIP HISTORY

Mar 13, 2021



Matheus Gallas-Lopes

Universidade Federal do Rio Grande do Sul

Jun 15, 2021



ngelo Piato Universidade Federal do Rio Grande do Sul

PROTOCOL INTEGER ID

48207

protocol. <a href="https://dx.doi.org/10.17504/protocols.io.btbpnimn">https://dx.doi.org/10.17504/protocols.io.btbpnimn</a>

Citation: Adrieli Sachett, Matheus Gallas-Lopes, Greicy M M Conterato, Ana P Herrmann, Angelo Piato (06/15/2021). Antioxidant activity by DPPH assay: in vitro

```
⊗Gloves Contributed by users
⊗96 well plate Contributed by users

⊗1.5 mL Eppendorf tubes Contributed by users
⊗Surgical mask Contributed by users

⊗Micropipette (0.5 - 10 μL) Contributed by users
⊗Micropipette (100 - 1000 μL) Contributed by users

⊗pH meter Contributed by users
⊗Synergy™ HTX Multi-Mode Microplate Reader Contributed by users

⊗Multichannel pipette (5 μL; 30- 300 μL) Contributed by users

⊗22-Diphenyl-1-picrylhydrazyl Sigma-aldrich Catalog #D9132
Step 1.1

⊗Ethanol Merck

Millipore Catalog #100983
Step 1.1
```

SAFETY WARNINGS

Use personal protective equipment (including lab coat, masks, and gloves) whenever manipulating chemical and biological samples. Make sure to read all Safety Data Sheets for the reagents.

## Preparing the reagents

- 1 The first step is to prepare the reagents to be used in this protocol;
  - 1.1 2,2-Diphenyl-1-picrylhydrazyl (DPPH) [M]0.24 mg/mL:

1.1.1 Weigh **0.024** g of DPPH in a piece of aluminum foil;

**⊠**22-Diphenyl-1-picrylhydrazyl **Sigma-**

aldrich Catalog #D9132

1.1.2 Add **50 mL** of absolute ethanol to the beaker to dissolve the salt;

⊠ Ethanol Merck

Millipore Catalog #100983

- 1.1.3 Transfer your solution to a **□100 mL** volumetric flask;
- 1.1.4 Using absolute ethanol, complete the solution's volume to reach **100 mL**;
- 1.1.5 Read the absorbance of the solution at **517 nm** in a microplate reader;

If the absorbance of the solution is above 1.12, dilute the solution by adding absolute ethanol; If the absorbance of the solution is lower than 1.08, add DPPH;

Incubation of the samples

1d

2 🔲 🎾

To optimize the reaction, an incubation step is needed.

2.1 Prepare 1.5 mL microtubes, to be used to store the samples, with the correct information. Wrap the microtubes in aluminum foil. The number of microtubes depends on the number of samples. You should provide at least five replicates (n = 5) of each sample with at least one control tube per sample. Perform the test at least two times to ensure the results are as correct as possible.

 For each sample, fill the plastic microtubes as described below. Using a micropipette fill the tubes in this order: Absolute ethanol + DPPH + Sample. Mix the solution with the pipette tip to homogenize the content;

Α	В	С	D
Microtubes	EtOH (µL)	DPPH (µL)	Sample (µL)
Control	20	280	-
Sample	10	280	10

2.3 Incubate all your samples in the dark at § 25 °C for © 24:00:00;

1d

06/15/2021

Reading your samples



Prepare to read the absorbance of your samples in a microplate reader;

Use a conventional 96-well microplate to run your samples. Before start pipetting, each well of the microplate should be marked for sample identification.

3.2

■250 µL of absolute ethanol to work as the negative control for the absorbance reading of all of your samples;

Read the absorbance of the samples at **517 nm** in a microplate reader;

Calculating data and determining results

- Prepare to analyze the results obtained after reading the absorbance of the samples;
  - Calculate the mean absorbance of your replicates;
  - Determine the percentage of inhibition of the DPPH radical: 4.2

% of inhibition of the DPPH radical= 
$$\left[\frac{Abs_{control} - (Abs_{sample} - Abs_{blank})}{Abs_{control}}\right] \times 100$$

 $Abs_{control}$  = The absorbance of DPPH

Abs<sub>sample</sub> = The absorbance of your sample

Abs<sub>blank</sub> = The absorbance of ethanol negative control

mprotocols.io

protocol. <a href="https://dx.doi.org/10.17504/protocols.io.btbpnimn">https://dx.doi.org/10.17504/protocols.io.btbpnimn</a>

- 4.3 Results should be expressed as % of inhibition of the DPPH radical.
- Determine the median effective concentration (EC50) for each sample, which expresses the minimum amount of the extract capable of reducing the initial concentration of DPPH radical by 50%. Calculated in a nonlinear regression model using GraphPad Prism.

Prism 8.0 © by GraphPad

- 5.1 On an XY graph, plot the concentrations of the samples in column X and the results of the DPPH assay (% inhibition) in column Y.
- 5.2 Transform the values to logarithmic scale:

Analyze -> Transform -> Ok -> transform X values using = X = Log(x) -> Ok.

5.3 Perform a nonlinear regression analysis:

Analyze -> XY analyzes -> Nonlinear regression -> Ok -> Dose responds stimulation— Log (agonist) VS. Normalized response -> variable slope -> Ok.

5.4

Results should be expressed as EC50.