



Mar 29, 2021

# Single-cell suspension preparation from Human bronchial biopsies to perform scRNA-sequencing using 10x chromium

Martijn Nawijn<sup>1</sup>, Leonie Apperloo<sup>1</sup><sup>1</sup>Department of Pathology and Medical Biology, GRIAC Research Institute, University of Groningen, University Medical Center Groningen, Groningen, Netherlands

1

Works for me

This protocol is published without a DOI.

University medical center Groningen

m.khan

SUBMIT TO PLOS ONE

## ABSTRACT

This protocol enables the dissociation of human bronchial biopsies into a single-cell suspension starting from fresh materials. These single cells are suitable for use in 10x chromium kits for single-cell RNA sequencing.

## PROTOCOL CITATION

Martijn Nawijn, Leonie Apperloo 2021. Single-cell suspension preparation from Human bronchial biopsies to perform scRNA-sequencing using 10x chromium. **protocols.io**  
<https://protocols.io/view/single-cell-suspension-preparation-from-human-bron-btf2njqe>

## KEYWORDS

Human bronchial biopsies, Single-cell suspension, 10x chromium, scRNA sequencing, discovAIR

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Mar 18, 2021

## LAST MODIFIED

Mar 29, 2021

## PROTOCOL INTEGER ID

48346

**Solution/buffers required:****1. Tissue dissociation solution per sample**

Reagents	Required amount
*HBSS +1% P/S	3ml
DNase I (Roche Diagnostics GmbH, REF #: 10104159001)	0,0003 grams
Collagenase D (Roche Diagnostics GmbH, REF #: 11088866001)	0,003 grams

\*Hanks' Balanced Salt solution (HBSS with calcium, magnesium and phenol red, 500mL, Lonza, Cat #: BE10-508F), supplemented with 1% P/S (Penicillin-Streptomycin, Gibco, Cat #:150700630).

**2. Resuspension buffer**

Reagents	Required amount
PBS pH 7.4 (Gibco, REF #:10010-15, LoT #: 2288905)	1x
BSA (7.5% in DPBS, SIGMA, Cat #: A8412)	0.04% (400µg/ml)

To prepare resuspension buffer in 50mL volume, take (266.66µl) from the 7.5% BSA stock solution and fill the tube up to 50mL with 1x PBS solution.

**3. Red blood cells lysis buffer (Ammonium Chloride Kalium, ACK 10x), pH 7.4**

Reagents (conc)	Required amount
155mM Ammoniumchloride (NH <sub>4</sub> Cl) (Sigma-Aldrich, CAS #: A9434)	41.45 grams
10mM Kaliumhydrogencarbonat (KHCO <sub>3</sub> ) (Merck, CAS #:237205)	5.0 grams
0.1mM EDTA (Invitrogen, Cat #: AM9912)	186 mg
Deionized water	500 mL

\* dilute the buffer to 1x before use.

**4. Cell counting**

- Trypan Blue Solution, (0.4% Sigma, Lot #:RNBj6324)
- Hemocytometer (BRAND™, 719520; Bürker Türk Counting Chamber)

**Before you start**

- 1 All steps should be performed in a sterile environment.  
Prepare all the buffers and reagents at room temperature as described in the "Materials" section.
- 2 Freshly prepare the tissue dissociation solution.  
  
Dilute the 10x RBCs lysis buffer to 1x.
- 3 Collect the tissue in a coulter counter cup filled with 15ml cold PBS (1x). Keep the tissues on ice before starting.

**Preparation of single-cell suspension**

- 4 Directly proceed by transferring the biopsies to a petri dish using sterile tweezers.
- 5 Cut the biopsies using a scalpel into smaller pieces at room temperature.
- 6 Transfer the pieces to the FACS tube containing 1ml of tissue dissociation solution (HBSS +1%P/S + Collagenase D +

DNaseI). Rinse the petri dish with tissue dissociation solution to collect all the materials.

- 7 Incubate for 1 hour @37°C, shake regularly. After this incubation keep the samples on ice.
- 8 Strain the mixture through a 70µm filter and collect the filtrate in a 50ml tube while keeping the samples on ice.  
Push the sample through the filter using the back of a plunger.  
Rinse the FACS tube with 1 ml cold HBSS+1%P/S and add it to the filter.  
Wash the filter with 1 ml HBSS+1%P/S. Collect any remaining cell suspension materials from the back of the filter using a pipette.
- 9 Centrifuge @590xg, 5 min at 4°C .
- 10 Resuspend the pellet in 500 µL resuspension buffer (PBS+0,04% BSA) and transfer the sample to a 1.5ml Low-Bind tube. Add 1 ml RBCs lysis buffer (1x) and incubate for ~2-5 minutes on ice.
- 11 Centrifuge @590xg, 5 min at 4°C.
- 12 Resuspend the pellet in 50µL resuspension buffer (PBS+0,04% BSA).
- 13 Count the cells using Trypan blue staining (1:1 ratio of Trypan blue to sample) and a manual hemocytometer. While keeping the rest of the cells on ice.

13.1 Proceed with 10x chromium protocol if 80-90% of the cells are alive.

#### 10x single-cell RNA-seq kit instructions

- 14 Continue with 10x Chromium (Next GEM Single Cell 3' reagent kit) and load 10,000 cells on a single lane.