

Sep 05, 2024

## Probe qPCR assay for Charybdis japonicum

This protocol is a draft, published without a DOI.

Jackson Treece<sup>1</sup>, Gert-Jan Jeunen<sup>1</sup>, Jo-Ann Stanton<sup>1</sup>, Neil Gemmell<sup>1</sup>

<sup>1</sup>University of Otago



Jackson Treece

University of Otago

OPEN  ACCESS



External link: <https://www.biosecurity-toolbox.org.nz/>

**Protocol Citation:** Jackson Treece, Gert-Jan Jeunen, Jo-Ann Stanton, Neil Gemmell 2024. Probe qPCR assay for Charybdis japonicum . protocols.io <https://protocols.io/view/probe-qpcr-assay-for-charybdis-japonicum-dki64uhe>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** September 02, 2024

**Last Modified:** September 05, 2024

**Protocol Integer ID:** 106814

**Keywords:** charybdis, charybdis japonicum, invasive crustacean, invasive crab, charybdis qpcr, charybdis japonicum qpcr, biosecurity toolbox, marine biosecurity, marine biosecurity toolbox

**Funders Acknowledgement:**

**New Zealand Ministry for  
Business, Innovation and  
Employment**

**Grant ID:** CAWX1904



## Abstract

*Charybdis japonica* or the Asian Paddle Crab is a non-native crab species to Aotearoa, New Zealand. It was first detected in 2000 and has since spread along the eastern Northland coast where it is a threat to local native species. This test was developed to aid in detecting the invasive crustacean's presence using environmental DNA sampling.

The probe based qPCR assay targets a 208bp portion of the COX-1 gene with an indicative limit of detection of 1.0E-4 ng/ $\mu$ l genomic DNA input or 92 molecules/ $\mu$ l with a synthetic standard (IDT gBlock). Species specificity was successfully challenged against a panel of 14 New Zealand native crab species. Further challenge with local species is recommended.



## Materials

- 1 Meridian Bioscience SensiFAST Probe Lo-ROX Mix (2x)  
PCR quality water  
Template DNA  
Oligonucleotides (10uM concentration):  
Forward primer: Cj\_1295F ATGGCATGCCTCGTCGTTAT  
Reverse primer: Cj\_1498R CTGCAGGTGGGTAAGAGTGG  
FAM Probe: Cj\_1323P CCCCATGCATACACTACGT

## Master Mix

- 2  
  
10.0 µl SensiFAST Probe 2x Mix  
0.8 µl 10µM Forward primer  
0.8 µl 10µM Reverse primer  
0.2 µl 10µM FAM Probe  
6.2 µl PCR Water  
2.0 µl Template  
**20µl total reaction volume**

## Thermal Conditions

- 3 95°C - 3 minutes  
**95°C - 15 seconds**  
**60°C - 30 seconds x45 cycles**

QuantStudio 3 (ThermoFisher Scientific) used during assay development.

## Control Sequence

- 4 Synthesized positive control sequence (IDT gBlocks)

tattactgccattcttactcctctcgcttctgtgtagctggagctattactatattattaacagaccgaaactaaatacttcattct  
tcgacctgctggggggtggggatcctgttcttcaacacttggttgattcttggacatcctgaagtttatatttaacttacctgcc  
ttcggaataattcacatattgtcagacaagaatctggtaaaaaggaatcatttggtactttgggtataatttatgctataatggctatt  
ggcattttaggctttattgtatgagctcatcacatgtttacagttggaatagatgtagatactcgagcatattttacatcagctactata  
attattgctgtgcctactggaattaaaatttttagatgacttagaactcttcacggcactcaaattaactatagaccttctatattatgg  
gccctaggattattttcttattcacagtaggaggtcttacaggagtagtttagcaaactcttcattgatattattcttcacgacacat  
actatgtggtagccacttcattacgttcttctataggagctgtatttggaattttcgccggcattgcccactgattcccctattttac  
cggactttcattaaatcctaaatgaataaaaaattcattttactattatatttttggtgtaaatattactttcttctcaacacttttagga  
ctaaatggcatgcctcgctgttattctgactacccgatgcatacactacgtgaaatattgtctcttctataggatcaatagatcact

tattgtataactaattttataattgtaattgagaagcacttatctcaaactcgacctgttctatcttctccttctccttctctattgaat  
gaaaccactcttaccacactgcagatcactcatatagaaatccccctaactcttaacttctaaaatggcagatagatgtatagg

## Species for Assay Development and Specificity Testing

- 5 Many thanks to the National Institute of Water and Atmospheric Research (NIWA), the Marine Invasives Taxonomic Service (MITS) and Deakin University for contributions from their invertebrate tissue collection. Donation ID: 2023-032 D

Species	Accession
Austrohelices crassa	MITS_74324
Carcinus maenas	N/A (Deakin)
Charybdis japonica	MITS_70207
Elamena producta	NIWA_89937
Halicarcinus cookii	NIWA_81032
Halicarcinus whitei	NIWA_147705
Heterozius rotundifrons	NIWA_163347
Liocarcinus corrugatus	NIWA_73086
Nectocarcinus antarcticus	NIWA_74132
Notomithrax minor	NIWA_92580
Notomithrax peronii	NIWA_72993
Ovalipes catharus	NIWA_60579
Pilumnus lumpinus	NIWA_75667
Pilumnus novaezealandiae	NIWA_92423
Pinnotheres atrinicola	NIWA_74906