




DEC 14, 2022

 Eosin Y**This protocol is published without a DOI.**joshdm¹¹UC Davis


joshdm

COMMENTS 0

WORKS FOR ME

1

ABSTRACT

Preparation of stock and dilution of  Eosin Y powder **Fisher Scientific Catalog #E511-25**

PROTOCOL CITATION

joshdm 2022. Eosin Y. **protocols.io**
<https://protocols.io/view/eosin-y-ckhqut5w>



LICENSE

_____ This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 14, 2022

LAST MODIFIED

Dec 14, 2022

PROTOCOL INTEGER ID

74000

PARENT PROTOCOLS




In steps of
[H&E staining](#)

MATERIALS TEXT

 Eosin Y powder **Fisher Scientific Catalog #E511-25****Prepare 1% Eosin Y stock**

1 Weigh  2 g of  Eosin Y powder **Fisher Scientific Catalog #E511-25** and add it to a 1000mL flask

Note

- 2 Add  40 mL of ddwater and swirl to mix
- 3 add  160 mL of 95% ethanol and swirl to mix
- 4 add  1 mL of Glacial Acetic acid

Preparing working dilution of Eosin Y

- 5 Use 80% ethanol to dilute stock to the desired concentration of eosin.

Note