



NOV 08, 2023

OPEN ACCESS



DOI:
dx.doi.org/10.17504/protocols.io.j8nlko6d5v5r/v1

Protocol Citation: Ching-Chieh Chou, Judith Frydman 2023. Neuronal transdifferentiation from human primary adult fibroblasts. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.j8nlko6d5v5r/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited


Protocol status: Working
 We use this protocol and it's working

Created: Nov 08, 2023

Neuronal transdifferentiation from human primary adult fibroblasts

Ching-Chieh Chou¹, Judith Frydman¹

¹Department of Biology, Stanford University



Ching-Chieh Chou
 Stanford University

ABSTRACT

This is a protocol for the direct conversion of human primary fibroblasts into neurons using a combination of transcription factor- and small molecule-based approach. The majority of converted neurons showing characteristics of cortical glutamatergic neurons.

Last Modified: Nov 08,
2023

PROTOCOL integer ID:
90649

Keywords: ASAPCRN,
Neuron, Transdifferentiation

Funders

Acknowledgement:

Aligning Science Across
Parkinson's
Grant ID: ASAP-000282

Preparation of culture medium and reagents

- 1
 1. Fibroblast medium @4°C lasts for 1 month. A total of 500 mL:
 - 435 ml DMEM (High glucose, GlutaMAX) (Thermo Fisher Scientific)
 - 50 ml Fetal bovine serum (FBS) (Thermo Fisher Scientific)
 - 5 ml Penicillin-Streptomycin (Thermo Fisher Scientific)
 - 5 ml MEM NEAA (100X; Thermo Fisher Scientific)
 - 5 ml Sodium pyruvate (100 mM; Thermo Fisher Scientific)
 - 0.5 ml beta-mercaptoethanol (55mM; Thermo Fisher Scientific)
 - Sterilized by filtration through a 0.22 µm filter
 2. Neuronal reprogramming medium @4°C lasts for 1 month. A total of 500 mL:
 - 240 mL DMEM/F-12 (Thermo Fisher Scientific)
 - 240 mL Neurobasal Medium (Thermo Fisher Scientific)
 - 10 mL B-27 Supplement (50X) (Thermo Fisher Scientific)
 - 5 mL N-2 Supplement (100X) (Thermo Fisher Scientific)
 - 1.25 mL GlutaMAX Supplement (Thermo Fisher Scientific)
 - 5 ml Penicillin-Streptomycin (Thermo Fisher Scientific)
 - Sterilized by filtration through a 0.22 µm filter
 3. Neuronal maturation medium @4°C lasts for 1 month. A total of 500 mL:
 - 480 mL BrainPhys Neuronal Medium (STEMCELL Technologies)
 - 10 mL B-27 Supplement (50X) (Thermo Fisher Scientific)
 - 5 mL N-2 Supplement (100X) (Thermo Fisher Scientific)
 - 1.25 mL GlutaMAX Supplement (Thermo Fisher Scientific)
 - 5 ml Penicillin-Streptomycin (Thermo Fisher Scientific)
 - Sterilized by filtration through a 0.22 µm filter
 4. Hexadimethrine Bromide (Polybrene) (Sigma-Aldrich)
 5. Forskolin (Sigma-Aldrich)

6. Dorsomorphin (Tocris)
7. SB 431542 (Tocris)
8. XAV939 (Stemgent)
9. Doxycycline (Cayman)
10. Puromycin (Thermo Fisher Scientific)
11. BDNF Recombinant Protein (Peprotech)
12. NT-3 Recombinant Protein (Peprotech)
13. Poly-L-ornithine hydrobromide (Sigma-Aldrich)
14. Human rhLaminin-521 (Thermo Fisher Scientific)
15. Human Vitronectin (Thermo Fisher Scientific)

Transduction/induction/selection of human fibroblasts

2 **Day -1:** Plate human fibroblasts at 200K cells per well of a poly-L-ornithine coated 6-well plate.

2.1 Dilute sterile-filtered 0.01% poly-L-ornithine in PBS at a 1-to-10 ratio.

2.2 Coat the plate and incubate in a 37°C incubator for 2 hr.

2.3 Wash the plate with H₂O for 5 min on a shaker and repeat for three times. Airdry the plate before plating the cells.

2.4 Rinse the cells with sterile 1x PBS to remove all complete medium.

- 2.5** Add 0.05% Trypsin-EDTA to cover the bottom of the dish. Incubate at 37°C with 5% CO₂ for ~5 minutes. Cells should round up and become dislodged.
- 2.6** Neutralize trypsin-EDTA activity by adding the complete media.
- 2.7** Gently pipette to resuspend the cells and transfer the cell-containing medium to a 15 mL conical tube.
- 2.8** Centrifuge at 200xG for 5 min to pellet the cells.
- 2.9** Remove the supernatant. Cell numbers are determined for each sample by a hemocytometer.
- 2.10** Add fresh culture medium to adjust the cell concentration to 200K cells per well and each well contain 2 mL of culture medium.
- 3** **Day 0:** Infect cells with 2 mL of viral supernatants per well of a 6-well plate.
- 3.1** Make virus dilutions by adding each viral supernatant to cold Fibroblast medium. 100 to 400 µL of each viral supernatant (M2rtTA, Ngn2-puro, Ascl1, Brn2, Myt1l) depending on viral titers.

- 3.2** Combine these viral supernatants and add cold Fibroblast medium to reach a total of 12 mL for one 6-well plate.
- 3.3** Add 4 µg/mL Polybrene to the virus dilutions.
- 3.4** Remove old Fibroblast medium and incubate cells with virus dilutions for 24 hr.
- 4** **Day 1:** Remove virus-containing medium and add fresh Fibroblast medium plus 1 µg/mL doxycycline.
- 5** **Day 2:** Change to fresh Fibroblast medium plus doxycycline and puromycin.
- 6** **Day 4:** Coat glass coverslips or plates with Vitronectin/rLaminin-521 and replate the transduced cells.
- Note**
- The conversion efficiency and cell survival are much better after replating as compared with no replating.
- 6.1** Coat with vitronectin (1:100 dilution of stock concentration 0.5 mg/mL in PBS) for 1 hr at room temperature. No need to rinse for vitronectin-coating coverslips or plates.
- 6.2** Coat with rhLaminin-521 (1:100 dilution of stock concentration 100 µg/mL in DPBS with Ca²⁺ and Mg²⁺) for 2 hr at 37°C.

- 6.3** In parallel, trypsinize cells by 0.05% Trypsin-EDTA for 5 min at 37°C and perform PSA-NCAM selection using magnetic-based cell sorting.
- 6.4** Replate the PSA-NCAM+ cells in Fibroblast medium plus doxycycline. Generally ~50K cells per cm².
- 7** **Day 5:** Switch to Neuronal reprogramming medium plus small molecules.
- 5 µM Forskolin
 - 2 µM Dorsomorphin
 - 10 µM SB 431542
 - 2 µM XAV939
 - 1 µg/mL Doxycycline
- 8** Change half of the culture medium every 2 to 3 days.
- 9** **Day 12,** Switch to Neuronal reprogramming medium supplemented with small molecules and neurotrophic factors.
- 5 µM Forskolin
 - 2 µM Dorsomorphin
 - 10 µM SB 431542
 - 2 µM XAV939
 - 1 µg/mL Doxycycline
 - 10 ng/mL BDNF
 - 10 ng/mL NT-3
- 10** Change half of the culture medium every 2 to 3 days.
- 11** **Day 20,** switch to Neuronal maturation medium supplemented with small molecules and neurotrophic factors.
- 5 µM Forskolin
 - 2 µM Dorsomorphin

- 10 μ M SB 431542
- 1 μ g/mL Doxycycline
- 10 ng/mL BDNF
- 10 ng/mL NT-3

Note

Suggest half medium change to gradually replace the neuronal reprogramming medium with maturation medium. This reduces the exposure of neurons to the air.

- 12 Change half of the culture medium every 3 to 4 days.
- 13 Neurons will mature on and after 5 weeks in culture and be ready for experiments
- 14 Neurons can be cultured for about 9 weeks.