



May 02, 2022

# Measurement of pancreatic islet beta-cell proliferation by flow cytometry

Julien Ghislain<sup>1</sup>, Vincent Poitout<sup>2</sup><sup>1</sup>Montreal Diabetes Research Center, CRCHUM, Montréal, QC, Canada.;<sup>2</sup>Montreal Diabetes Research Center, CRCHUM, and Department of Medicine, Université de Montréal, Montréal, QC, Canada.

1

[dx.doi.org/10.17504/protocols.io.n92ldz8d9v5b/v1](https://dx.doi.org/10.17504/protocols.io.n92ldz8d9v5b/v1)

Laboratory of Vincent Poitout

Tech. support email: [julien.ghislain.chum@ssss.gouv.qc.ca](mailto:julien.ghislain.chum@ssss.gouv.qc.ca)

Julien Ghislain

This protocol describes the steps to assess pancreatic beta-cell proliferation by flow cytometry. It is suitable for islets isolated from both rodents and humans. We routinely apply this protocol to quantify total and proliferating beta-cells but, with the appropriate antibody-fluorophore conjugates, can be easily adapted to assess other islet cells types (eg. alpha and delta) and markers (eg. apoptosis). Briefly, islets labelled with the proliferation marker EdU are dispersed, dead cells labeled and fixed. EdU is detected using the Click-iT Plus EdU Flow Cytometry Assay Kit and beta-cells immunostained with fluorophore-coupled anti-insulin antibodies. Beta-cell proliferation is then determined by flow cytometry by calculating the percentage of double-positive cells for EdU and insulin over the corresponding total insulin-positive cell population.

DOI

[dx.doi.org/10.17504/protocols.io.n92ldz8d9v5b/v1](https://dx.doi.org/10.17504/protocols.io.n92ldz8d9v5b/v1)

Julien Ghislain, Vincent Poitout 2022. Measurement of pancreatic islet beta-cell proliferation by flow cytometry. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.n92ldz8d9v5b/v1>

CIHR

Grant ID: 0406014428

NIH

Grant ID: 5R01DK058096-15

protocol

Pronounced proliferation of non-beta cells in response to beta-cell mitogens in isolated human islets of Langerhans. Maachi H, Ghislain J, Tremblay C, Poitout V. Sci Rep. 2021 May 28;11(1):11283. doi: 10.1038/s41598-021-90643-3. PMID: 34050242.

flow cytometry, pancreatic islet, beta-cell, proliferation

\_\_\_\_\_ protocol ,

Pronounced proliferation of non-beta cells in response to beta-cell mitogens in isolated human islets of Langerhans. Maachi H, Ghislain J, Tremblay C, Poitout V. Sci Rep. 2021 May 28;11(1):11283. doi: 10.1038/s41598-021-90643-3. PMID: 34050242.

Mar 23, 2022

May 02, 2022

59811

 [Phosphate Buffered Saline 10x \(solution\)](#) **Bio Basic**

**Inc. Catalog #PD8117**

 [EDTA \(0.5 M\), pH 8.0](#) **Life**

**Technologies Catalog #AM9260G**

 [StemPro™ Accutase™ Cell Dissociation Reagent](#) **Thermo Fisher**

**Scientific Catalog #A1110501**

 [RPMI 1640 Medium](#) **Thermo Fisher**

**Scientific Catalog #11875093**

 [FBS Invitrogen - Thermo Fisher](#)

 [LIVE/DEAD™ Fixable Aqua Dead Cell Stain Kit, for 405 nm excitation](#) **Thermo**

**Fisher Catalog #L34965**

 [Click-iT™ Plus EdU Alexa Fluor™ 488 Flow Cytometry Assay Kit](#) **Thermo**

**Fisher Catalog #C10632**

 [BSA fraction V](#) **Gibco - Thermo**

**Fischer Catalog #15260037**

 [Alexa Fluor® 647 Mouse Anti-Insulin](#) **BD**

**Biosciences Catalog #565689**

Islet dispersion and fixation

3m

17m

## 1 Dispersing islets into single cells

Wash approximately 100-200 islets two times with ice cold PBS containing **12 millimolar (mM)** EDTA (**pH8**) in a 1.5 ml microcentrifuge tube using a microcentrifuge at **339 x g, 4°C, 00:03:00**.

Keeping the tubes on ice add **200 µL - 300 µL** of ice cold Accutase solution.

Transfer the tubes to a water bath at **37 °C** and gently mix the islet suspension using a 1 ml pipette tip during a maximum of **00:10:00**.

Return the tubes to ice and add an equal volume of islet media (eg.10% FBS in RPMI 1640), mix and centrifuge at **339 x g, 4°C, 00:04:00**.

1h 9m

## 2 Live/Dead staining and fixation

Discard the supernatant and wash the cells once with PBS and suspend in **1 mL** of PBS on **ice**. Add **2 µL** (for human islets) or **1 µL** (for rat and mouse islets) of the reconstituted Aqua LIVE/DEAD fluorescent reactive dye, mix and incubate **ice** for **00:30:00** (**00:20:00** for rat or mouse islets), protected from light.

Centrifuge at **339 x g, 4°C, 00:04:00** and wash once with **1 % (v/v)** BSA in PBS **ice**.

Suspend the cell pellet and add **100 µL** of Click-iT™ plus EdU kit fixative (Component D), mix and incubate for **00:15:00** at **Room temperature** protected from light.

Add **1 mL** of **1 % (v/v)** BSA in PBS **ice**, centrifuge as above and wash with **1 mL** Click-iT™ plus EdU kit perm/wash reagent. (To store for up to one week at **4 °C**, perform instead 2 washes with **1 % (v/v)** BSA in PBS **ice**.)

Staining and flow cytometry

55m

55m

## 3 Click-iT™ plus EdU reaction

Centrifuge as above and suspend the cell pellet in **50 µL** Click-iT™ plus EdU kit perm/wash reagent and mix. Incubate for **00:10:00** **37 °C**. During the incubation step prepare the Click-iT™ plus EdU kit reaction cocktail.

Start by preparing a 1X Click-iT™ plus EdU kit buffer additive by diluting the 10X stock solution in water.

Then, prepare the Click-iT™ plus EdU kit reaction cocktail at **Room temperature** and use within

**00:15:00**.

### 2 rxn

**875 µL**

**20 µL**

**5 µL**

**100 µL**

### 5 rxn

**2.188 mL**

**50 µL**

**12.5 µL**

**250 µL**

PBS

Component F

488-dye azide

1X Buffer additive

Add **450 µL** of Click-iT™ plus EdU kit reaction cocktail to each tube, mix and incubate for **00:30:00** at **Room temperature**, protected from light.

Centrifuge and wash the cells once with **1 mL** of Click-iT™ plus EdU kit perm/wash reagent.

20m

## 4 Antibody staining

Centrifuge and remove the supernatant. Dislodge the cell pellet and resuspend in **50 µL** of Click-iT™ plus EdU kit perm/wash reagent containing **2 µL** of Alexa Fluor® 647 Mouse Anti-Insulin and incubate for **00:20:00** at **Room temperature** in the dark.

Wash twice with **1 mL** Click-iT™ plus EdU kit perm/wash reagent and suspend in **200 µL** of **1 % (v/v)** BSA in PBS.

## 5 Flow cytometry

Dead-cell stain, EdU and insulin labelled cells are detected using 405-, 488-, 640-nm lasers coupled with 525/50-, 530/30-, 670/14-nm BP filters, respectively, or similar. Proliferation is calculated as the

percentage of double-positive  $\text{EdU}^+$  and  $\text{Ins}^+$  cells over the total  $\text{Ins}^+$  cell population.