



SEP 23, 2023

OPEN ACCESS



DOI:
dx.doi.org/10.17504/protocols.io.ewov1qyepgr2/v1

Protocol Citation: Elias Adriaenssens 2023. In vitro phosphatase assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.ewov1qyepgr2/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jun 27, 2023

Last Modified: Sep 23, 2023

PROTOCOL integer ID:
84093

Keywords: In vitro phosphatase assay

In vitro phosphatase assay

Elias Adriaenssens¹

¹Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna (AT)

ASAP Collaborative Research Network

Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

ABSTRACT

This protocol describes *in vitro* phosphatase assay.

ATTACHMENTS

[756-1924.pdf](#)

MATERIALS

Materials

- 6 well plates
- MnCl₂

Pierce & Warriner; Detergent Compatible Bradford Assay Kit **Thermo**
Fisher Catalog #23246

Lambda Protein Phosphatase - 100,000 units **New England**
Biolabs Catalog #P0753L


Lysis buffer

A	B
HEPES pH 7.4	50 mM
NaCl	150 mM
MgCl ₂	2.5 mM
DTT	2 mM
NP-40	0.50%
protease inhibitor cocktail	

In vitro phosphatase assay

5m

- 1 Seed HAP1 wild-type or FIP200 knockout cells in 6 well plates and grow until confluency.




2 Collect cells by trypsinization and pellet by centrifugation at  300 x g, 4°C, 00:05:00 .

5m



3 After a PBS wash to remove the remaining cell medium, resuspend the cell pellets in lysis buffer.








4 Lyse the samples for  00:20:00  On ice . Clear the cell lysates by centrifugation at  20000 x g, 4°C, 00:10:00 .



30m



5 Determine the protein concentrations of the cleared protein lysates with the Pierce Detergent Compatible Bradford Assay Kit (23246, Thermo Fisher).

6 For both samples, wild-type and FIP200 knockout lysates, incubate  100 µg of cell lysate with  5 µL of 10x NEBuffer for Protein MetalloPhosphatases (P0753, New England Biolabs) and  5 µL of  10 millimolar (mM) of MnCl₂ to make a total reaction volume of  50 µL .



7 Add  1 µL of Lambda Protein Phosphatase (NEB) to the reaction and incubate the samples at  30 °C for the indicated time.



8 Terminate the phosphatase reactions by the addition of 6x Protein Loading dye and heat inactivation at  95 °C for  00:05:00 .

5m



9 Analyze the samples by western blot analysis as described above.

