

Version 2

Sep 10, 2020

Rosbash/Janelia StickLAMP Protocol V.2

Albert Yu¹, Tim Brown², Jasmine Quynh Le¹, Kristina Galatsis¹, Michael Rosbash¹

¹Brandeis University; ²Janelia Research Campus

1

Works for me

dx.doi.org/10.17504/protocols.io.bk89kzz6

XPRIZE Rapid Covid Testing

Albert Yu

ABSTRACT

A protocol for the detection of SARS-CoV-2 from saliva samples featuring a rapid purification step and a high-contrast colorimetric readout. Saliva is first inactivated using a 100x inactivation reagent consisting of 2.5M TCEP, 100 mM EDTA, 1.2N NaOH solution diluted to approximately 1x final concentration and heated to 95C for 5 minutes. RNA is rapidly purified and concentrated with magnetic beads in a PEG/NaCl-based buffer using a 3D-printed magnetic stick that enables selective separation of beads without carryover of saliva contaminants. Beads are eluted directly into an RT-LAMP reaction mix, which uses a novel high contrast dye that turns from purple to clear when acidified by nucleic acid amplification products that enables unambiguous identification of successful amplification. This protocol is sensitive down to 1 copy/μl of SARS-CoV-2 in 300 μl of saliva. This degree of sensitivity enables faithful detection of SARS-CoV-2 even in pooled samples.

DOI

dx.doi.org/10.17504/protocols.io.bk89kzz6

PROTOCOL CITATION

Albert Yu, Tim Brown, Jasmine Quynh Le, Kristina Galatsis, Michael Rosbash 2020. Rosbash/Janelia StickLAMP Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bk89kzz6>
 Version created by Albert Yu

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 10, 2020

LAST MODIFIED

Sep 10, 2020

PROTOCOL INTEGER ID

41985

MATERIALS

NAME	CATALOG #	VENDOR
NaCl	53014	Sigma Aldrich
Twist synthetic SARS-CoV-2 RNA control	Mt007544.1	Twist Bioscience
SARS-CoV-2 Master Mix		
Actin Master Mix		
100x Inactivation Reagent		
Bead Mix		

NAME	CATALOG #	VENDOR
Magnetic Tips		
Heat Block at 65C		
Heat Block at 95C		
Magnetic Stick		

STEPS MATERIALS

NAME	CATALOG #	VENDOR
100x Inactivation Reagent		
Bead Mix		
Magnetic Tips		
NaCl	53014	Sigma Aldrich
SARS-CoV-2 Master Mix		
Actin Master Mix		
Water		
Twist synthetic SARS-CoV-2 RNA control	Mt007544.1	Twist Bioscience

MATERIALS TEXT

100x Inactivation Reagent
2.5M TCEP
150mM EDTA
1.2N NaOH

SARS-CoV-2/Actin Master Mix

12.5µl SARS-CoV-2/Actin Buffer/Dye/Primer Mix (Currently only available from us)
0.5µl WarmStart RTx NEB M0380L
1µl Bst2.0 NEB M0537L
11µl H2O

Primers used

See https://docs.google.com/spreadsheets/d/11n-9754VqtsXszTC2tUxFq-_gKlGgjPL-KtKqXVevH4/edit#gid=0

Bead Mix

See https://ethanomics.files.wordpress.com/2012/08/serapure_v2-2.pdf with 300µl beads instead of 1000µl

EQUIPMENT

NAME	CATALOG #	VENDOR
ThermoMixer	5382000023	
Magnetic Stick	None	

SAFETY WARNINGS

Do not open up PCR tubes after amplification.

BEFORE STARTING

Prepare:

Saliva collection kit (2.0ml Tube and funnel provided by us, or your own saliva collection device from standard labware, such as 1.5ml, 5ml, 15ml, or 50ml tubes. Saliva samples >1ml will likely have to be subsampled)

Magnetic stick

1 magnetic tip per sample

Bead mix: Let bead mix come to room temperature for 20 minutes prior to use, and ensure beads are suspended in solution by vortexing or pipetting up and down

130mM NaCl

Saliva Collection

- 1 Instruct patient to avoid food, drink, toothbrushing, and nasal sprays for a minimum of 🕒 **00:30:00** prior to sample collection
- 2 Begin pooling saliva in your mouth. Saliva production can be stimulated by thinking about food, or about the saliva collection itself.
- 3 Gently expel saliva into the funnel, tapping to collect in the tube, until amount of saliva is approximately flush with the base of the funnel 📏 **750 µl Approximately**

Inactivation

5m

- 4 Add inactivation reagent to approximately 1x final concentration. Reaction is tolerant of between 0.7x to 2x final concentration. 📏 **7.5 µl Approximately**



100x Inactivation Reagent

- 5 Invert 40 times to mix
- 6 Heat tube to approximately 🔥 **95 °C** for 🕒 **00:05:00** . Viral RNA release is similar between 93-98C. Use tube clip to prevent popping. ^{5m}







ThermoMixer
Benchtop Incubator

Eppendorf 5382000023






Any heat block will suffice



- 7 Remove tube from heat and let rest at  **Room temperature** for at least  **00:03:00** OR  **On ice** for at least  **00:00:30**.

Quick Purification and LAMP Reaction

1h

- 8 While tube is resting, aliquot  **25 µl** SARS-CoV-2 mastermix and  **25 µl** Actin mastermix to separate wells of PCR strip tube, 96-well plate, or 1.5ml tube per sample  **On ice**.

Per run, prepare two additional  **25 µl** SARS-CoV-2 mastermixes for positive and negative controls.



SARS-CoV-2 Master Mix



Actin Master Mix

Step 8 includes a Step case.

If pooling

step case

If pooling

Prepare one 25ul SARS-CoV-2 reaction and one 25ul Actin reaction per 5 samples

- 9 Add approximately 0.7x volumes of bead mix -  **525 µl Approximately**. Sample is tolerant of between 0.7x-1.2x volumes of bead mix. Pipette up and down to mix.



Bead Mix
View

Step 9 includes a Step case.

If pooling

If you would like to preserve some sample

step case

If pooling

Remove 60ul of inactivated saliva from 5 samples and add to a single tube, for a total of 300µl. Add 210µl bead mix to pooled tube.

10 Let stand at **Room temperature** for **00:03:00** 10m

11 Cap magnetic stick with a clean tip and dip in bead/sample mix for **00:02:00** , dipping up and down 5 times every 2m
 00:00:30 . Meanwhile, prepare **500 µl** 130mM NaCl in a separate 1.5ml or 2ml tube.



Magnetic Stick

Rosbash/Brown

None

Magnetic stick used for bead purifications



Magnetic Tips



NaCl

by Sigma Aldrich


Catalog #: 53014


12 Remove magnetic stick from sample and swirl in clean 130mM NaCl solution for **00:00:05** . Discard NaCl solution. 5s

13 Remove magnetic stick from wash sample and place in SARS-CoV-2 mix for **00:00:30** 30s

14 Remove magnetic stick from SARS-CoV-2 mix and place in Actin mix for **00:00:30** 30s

15 Add **5 µl** water to additional SARS-CoV-2 Mix (negative control) and **5 µl** synthetic Twist SARS-CoV-2 positive RNA control to additional SARS-CoV-2 Mix, prepared in Step 8.


Water


Twist synthetic SARS-CoV-2 RNA control
by Twist Bioscience
Catalog #: Mt007544.1

16 Cap tubes and place on  **65 °C** heating apparatus for  **00:40:00** .

40m

If using a thermal cycler, run with the following program:

65C for 40 minutes

4C indefinitely


ThermoMixer
Benchtop Incubator
Eppendorf 5382000023 
Any heat block will suffice



17 Remove tubes from heating apparatus and examine color change.

Positive		Negative		Inconclusive			
SARS-CoV-2	Actin	SARS-CoV-2	Actin	SARS-CoV-2	Actin		
							
SARS-CoV-2	Actin			Negative Control			
							
				Positive Control			
							

18 If a positive sample is found when pooling, re-test pooled samples individually.

