



MAR 07, 2024

## 🌐 Immunofluorescent Imaging and Analysis

📁 In 1 collection

Dan Tudorica<sup>1</sup>

<sup>1</sup>University of California, Berkeley



Dan Tudorica  
Hurley Lab, QB3, UC Berkeley

### ABSTRACT

Following treatment with antibodies/genetic methods of introducing fluorescent tags

OPEN  ACCESS



#### DOI:

[dx.doi.org/10.17504/protocols.io.6qpvr37r2vmk/v1](https://dx.doi.org/10.17504/protocols.io.6qpvr37r2vmk/v1)

**Protocol Citation:** Dan Tudorica  
2024. Immunofluorescent  
Imaging and Analysis.  
**protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.6qpvr37r2vmk/v1>

**License:** This is an open  
access protocol distributed under  
the terms of the [Creative  
Commons Attribution License](#),  
which permits unrestricted use,  
distribution, and reproduction in  
any medium, provided the original  
author and source are credited

**Protocol status:** Working

**Created:** Nov 09, 2023

**Last Modified:** Mar 07, 2024

## Preparation of chambered slides

- 1 Seed cells on 8-well chambered slide to be 80% confluent on day of analysis.
- 2 Wash cells once with PBS, then fix for 15 minutes using 4% PFA diluted in PBS. Wash cells once more with PBS.
- 3 Permeabilize cells for an additional 15 minutes using a 1x stock of SLO prepared in PBS + 10 mM TCEP. Prepare 1x stock according to the protocol in this article (<https://doi.org/10.7554/eLife.20378>)
- 4 Remove SLO, wash cells with PBS once
- 5 Incubate permeabilized cells in 2% BSA in PBS for 1 H at RT to block
- 6 Replace blocking solution with blocking solution + IFA competent primary antibody. For this study, used a 1:200 dilution of ab302494 anti pS72 Rab7 Rb antibody. Incubate at RT for 1 H
- 7 Wash cells 3 x 5 min with PBS to remove primary antibody

- 8 Replace PBS with a solution of fluorescently labelled secondary antibody diluted in 2% BSA in PBS. In this study, used ThermoFisher antibody A11008, Goat anti Rabbit Alexa 488. Incubate for 1 H in the dark.
- 9 Wash cells with 3 x 5 min of PBS, then leave cells sitting under PBS. Image immediately.

## Imaging

- 10 Prepare Confocal microscope with 60x immersion oil objective. Place slides on stage, and find focus.
- 11 Adjust laser power, pinhole, and gain for optimal brightness and resolution.
- 12 Capture each fluorescent channel individually in order to minimize risk of bleedthrough
- 13 In order to test for bleedthrough, check to make sure that when excitatory laser for a particular channel is dialed down to an intensity of 0, image is completely dark.
- 14 Acquire images and export to FIJI for analysis

## Colocalization analysis

- 15 Merge fluorescent channels into a single stack for colocalization analysis
- 16 For each channel, measure the max brightness of a known dark region of the image. For this study, I used the nucleus of a cell as an expected dark. Use FIJI's *math>subtract* command to subtract this value from each pixel in the image.
- 17 In each image, define ROIs around target cells. If a transient transfection was performed, cells should be expressing a moderate amount of the target protein.
- 18 Use the *Coloc 2i* plugin to calculate the Mander's colocalization coefficient with a Costes' threshold regression.
- 19 Perform this analysis on at least 20 different cells imaged from different parts of the well. The average of these cells constitute a single biological replicate. Repeat this experiment from the start several more times, and then perform statistics on aggregated experimental averages.