

Mar 17,
2020

SARS-CoV-2 virus plaque assays [Biosafety Level 3]

tenOever Lab¹

¹Icahn School of Medicine at Mount Sinai

1 Works for me dx.doi.org/10.17504/protocols.io.bdt6ni6me

Coronavirus Method Development Community



Emma Ganley
Protocols.io



ABSTRACT

Method to plaque the SARS-CoV-2 virus

Image attribution:

Y tambe / CC BY-SA (<http://creativecommons.org/licenses/by-sa/3.0/>)

Downloaded from: https://commons.wikimedia.org/wiki/File:Plaque_assay_macro.jpg

EXTERNAL LINK











<https://bit.ly/3d6OUny>

SAFETY WARNINGS








Biosafety Level 3

Media Recipes

- Media recipe for modified DMEM: [DMEM with 2% FBS]:
 - 475 ml DMEM (w/ L-Glutamate, Sodium pyruvate)
 - 5 ml 100x P/S
 - 5 ml 100x Non-Essential Amino Acids (NEAA) Solution
 - 5 ml 1M HEPES
 - 10 ml FBS

- 2 Media recipe for 2XMEM: Sterile filter with 22 µm filtration unit:
 -  **100 ml** 10x MEM
 -  **10 ml** 100x P/S
 -  **10 ml** 100x L-Glutamine
 -  **6 ml** 35%BSA
 -  **340 ml** sterile ddH2O
 -  **10 ml** 1M HEPES
 -  **24 ml** 5% NaHCO₃
- 3 Overlay Media (Ensure this is just warm to the touch when added to cells)
 -  **1 ml** PBS
 -  **31.5 ml** 2X MEM
 -  **17.5 ml** Oxoid Agar

SARS-CoV2 virus plaque assays (Biosafety level 3)

- 4 Seed Vero E6 cells in 6-well plates overnight at  **37 °C** (10⁶ cells/well)
- 5 Dilute the SARS-CoV-2 virus logarithmically in infection media
- 6 Replace media with modified DMEM and add diluted virus in each well (200 µl/well)
- 7 Allow virus adsorption for  **01:00:00** , agitating the plate every  **00:10:00** .
- 8 Heat the overlay media/oxoid agar mixture to ensure homogenous consistency and add it to each well (2ml/well), incubate at  **37 °C** incubator for  **72:00:00** .
- 9 Fix the plate with 5% formaldehyde (w/ methanol), (1.5 ml/per well)  **Overnight** to ensure virus inactivation and Stain with crystal violet for  **01:00:00** to determine viral titers (PFU/ml)



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited