



Aug 17, 2022

Diagnostic Restriction Digest (Instructor Protocol)

Brian Teague¹¹University of Wisconsin - Stout

1 Works for me

Share

This protocol is published without a DOI.

Yeast ORFans CURE

Brian Teague
University of Wisconsin - Stout

ABSTRACT

This is the instructor protocol for



Diagnostic Restriction Digest
by **Brian Teague**,
University of Wisconsin - Stout

PREVIEW

RUN

PROTOCOL CITATION

Brian Teague 2022. Diagnostic Restriction Digest (Instructor Protocol).
protocols.io
<https://protocols.io/view/diagnostic-restriction-digest-instructor-protocol-cffftjjn>



LICENSE

_____ This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 17, 2022

LAST MODIFIED

Aug 17, 2022

PROTOCOL INTEGER ID

68807

MATERIALS TEXT

-  **700 µL** or  **1.7 mL** tubes for aliquots

 **CutSmart® Buffer New England**

- **Biolabs Catalog #B7204S** Step 1

 **PvuII-HF - 5,000 units New England**

- **Biolabs Catalog #R3151S** Step 2

 **Diluent B - 5.0 ml New England**

- **Biolabs Catalog #B8002S** Step 2

SAFETY WARNINGS

None of the materials are hazardous.

HOWEVER, we are shedding nucleases -- enzymes that degrade DNA -- all the time. Wear lab coats and gloves to keep your samples nuclease-free.

Setup

1

 **CutSmart® Buffer New England**

Aliquot the **Biolabs Catalog #B7204S**

:  **20 µL** ul

aliquots, 1 per 4 students

2

 **PvuII-HF - 5,000 units New England**

Aliquot the **Biolabs Catalog #R3151S**

enzyme:

 **4 µL** of enzyme in  **16 µL** of

 **Diluent B - 5.0 ml New England**

Biolabs Catalog #B8002S

, 1 per 4 students.

Instructor Tips & Common Student Errors

3 Instructor Tips

- In a more advanced course, I would let students select which enzyme to use (chosen from the ones in the freezer.) In an intro course, I choose the enzyme for them. PvuII is a good choice because it has two cut sites in the backbone *and one in the GFP insert*. This makes it easy to distinguish between the correct plasmid and one with a GFP still present (maybe not glowing because of a mutation? Or some other error?)
- PvuII is a little on the expensive side -- but if you're giving the digest a full hour, diluting down to 2 units per ul gives two advantages: it decreases the amount used, and increases the volume of the pipetting step. Make sure you use the correct diluent, though.

- I used to include a heat-inactivation step. I don't any more because the SDS in the purple loading dye denatures the enzyme.
- Sometimes a student's miniprep isn't concentrated enough to get a full microgram of DNA into the digest. As long as they can get at least 500 ng in the digest, that should be enough to see on a gel. It's more acceptable to decrease the mass of DNA than it is to increase the volume -- contaminants in the miniprep often get in the way of the digest, particularly leftover ethanol.
- I don't have a positive control in this experiment as written. Maybe add one?

4 Common Student Errors

- Not mixing the reaction well.
- Not loading the entire reaction onto the gel. (Because the first one they did was a PCR, it's easy to assume that you just need 1-2 ul without reading the protocol carefully.)