



DEC 19, 2022

WORKS FOR ME

1

UPitt TriState SenNet TMC Tissue Collection

COMMENTS 0

DOI

dx.doi.org/10.17504/protocols.io.yxmvm294og3p/v1

Marta Bueno¹, Lanping Guo¹, Ricardo Pineda¹,
John Sembrat¹, [koenigshoffm](#)¹,
[Oliver Eickelberg](#)¹

¹Division of Pulmonary and Critical Care Medicine.
Department of Medicine. School of Medicine

TriState SenNet

Cellular Senescence Network (SenNet) Method
Development Community

Marta Bueno

[Division of Pulmonary and Critical Care Medicine. Department...](#)

ABSTRACT

This document outlines the tissue sampling of lung and heart specimens at the TriState SenNet TMC Biospecimen Core at the University of Pittsburgh, as part of the Cellular Senescence Network Program (SenNet).

DOI

dx.doi.org/10.17504/protocols.io.yxmvm294og3p/v1

PROTOCOL CITATION

Marta Bueno, Lanping Guo, Ricardo Pineda, John Sembrat, [koenigshoffm](#), [Oliver Eickelberg](#) 2022. UPi TriState SenNet TMC Tissue Collection. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.yxmvm294og3p/v1>

FUNDERS ACKNOWLEDGEMENT

TriState SenNET (Lung and Heart) Tissue Map and Atlas consortium – NIA
Grant ID: U54AG075931

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 18, 2022

LAST MODIFIED

Dec 19, 2022

PROTOCOL INTEGER ID

74184

Preparation

- 1 Remove specimen from shipping or storage container.
- 2 Place specimen inside the hood over dissection table and photograph prior to sampling.
- 3 Select the sites from which samples will be collected, annotate and mark (if serial sampling is done)
 - 3.1 For lung:
 1. Pleura (A)
 2. Broncho-vascular bundle (B)
 3. Parenchyma (C)
 - 3.2 For heart:
 1. Left Ventricle (D)
 2. Right Ventricle (E)
 3. Left Atrium (F)
 4. Right Atrium (G)
 - 3.3 (OPTIONAL - If needed)
Wash the full-thickness core in cold PBS in a petri dish.

Flash freezing and storing

- 4 Cut each of the selected regions from the previous step into 9 small chunks. Chunks should be small enough to just cover the bottom of a microfuge tube (~5mm x 5mm x 5mm).
- 5 Dab the small chunks on a sterile gauze pad and place each one into a labeled microfuge tube.

- 6 Place 9 small chunks into the empty tubes to be snap frozen.
Drop tubes into liquid Nitrogen.
Leave for ~2-10 minutes.
Remove and store in -80C freezer.

Storing in formalin

- 7 Place one piece of the selected regions into a tubes containing 20x volume of 10% formalin (to be processed into FFPE blocks).
Store in 4C fridge.

Note: Fixation occurs at a rate of 1 mm per hour. Over-fixation causes tight cross-linking and under-fixation can compromise the tissue.

- 7.1 Make a notation of the time the sample went into formalin.
- 7.2 Based on the size of the tissue, fix the required length of time (4 to 24 hours).
- 7.3 After the fixation period, remove the cassette and place in 70% ETOH until scheduled processing.

Storing for fresh processing

- 8 Place one piece of the selected regions into a tube with cold DMEM with anti/anti.
Store in 4C fridge.

Process as soon as possible.