



JUL 13, 2023

OPEN ACCESS

**DOI:**  
[dx.doi.org/10.17504/protocols.io.q26g74x79gwz/v1](https://dx.doi.org/10.17504/protocols.io.q26g74x79gwz/v1)

**Protocol Citation:** Mark Ellisman, Marianne Martone, Thomas Deerinck 2023. 0.3M Sodium Cacodylate Buffer pH 7.4 Stock Solution Recipe.

**protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.q26g74x79gwz/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
 We use this protocol and it's working

**Created:** Apr 26, 2022

**Last Modified:** Jul 13, 2023

**PROTOCOL integer ID:**  
 61479

**Keywords:** cacodylate ,  
 buffer, stock, solution

## 0.3M Sodium Cacodylate Buffer pH 7.4 Stock Solution Recipe

Mark Ellisman<sup>1,2</sup>, Marianne Martone<sup>1</sup>, Thomas Deerinck<sup>1,2</sup>

<sup>1</sup>National Center for Microscopy and Imaging Research, University of California, San Diego, La Jolla CA;

<sup>2</sup>Department of Neurosciences, University of California, San Diego, La Jolla CA

NCMIR@UCSD



NCMIR  
 University of California, San Diego

### ABSTRACT

Standard preparation protocol for stock 0.3M sodium cacodylate buffer.

### GUIDELINES

Cacodylate acid, sodium salt, trihydrate (MW = 214.03)

Vendor: Ted Pella

Product number: 18851

Calculation:

$M \text{ (moles/liter)} = g / MW \text{ (g/mol)} \times V \text{ (liter)}$

$0.3M = g / 214.03 \text{ g/moles} \times 1.0L$

$g = M \times V \times MW$

$g = 0.3 \text{ moles/liter} \times 214.03 \text{ g/moles} \times 1.0 \text{ liter}$

**$g = 64.209 \text{ grams for } 1.0 \text{ liter}$**

For 3L =  $3 \times 64.209 = 192.627 \text{ grams}$

## MATERIALS

- Sodium Cacodylate trihydrate salt [TED PELLA, INC. (Catalog no. 18851)]
- Double distilled H<sub>2</sub>O (ddH<sub>2</sub>O)
- Appropriately sized container
- Concentrated hydrochloric acid

## SAFETY WARNINGS



Always wear appropriate PPE.  
Cacodylate is always toxic.

- 1** Add 900 ml of DDH<sub>2</sub>O to 64.209 grams of sodium cacodylate in a large 1.0L or greater container, such as a capped glass bottle. Shake and wait until powder is fully dissolved.
- 2** Maximum HCL concentration is 12.18mole/liter (36-38%). Start with 1.3mL of 12.18 M HCL per 900 ml of sodium cacodylate dissolved solution. Make sure solution is mixed thoroughly.
- 3** Check pH. Keep adding 50-100µl of HCL, mix, and check pH again until pH 7.4 has been reached. (Approximately 1.4 ml of conc. HCL total).
- 4** Top off to 1000ml with DDH<sub>2</sub>O.
- 5** Do a final pH check and keep stock solution refrigerated at 4C.