



MAR 17, 2023

Economic and easy bacterial and yeast colony PCR with 2x premix Gotaq-Green (Promega) or DreamTaq-Green (Thermo Fisher)

Darek Abramczyk¹¹University of Edinburgh

Darek Abramczyk

ABSTRACT

Slightly modified method Gotaq [Promega](#) or DreamTaq ([Thermo](#)), economical, quick and ready-to-load. Perfect for yeast genotyping and screening

OPEN ACCESS

Protocol Citation: Darek Abramczyk 2023. Economic and easy bacterial and yeast colony PCR with 2x premix Gotaq-Green (Promega) or DreamTaq-Green (Thermo Fisher) . **protocols.io** <https://protocols.io/view/economic-and-easy-bacterial-and-yeast-colony-pcr-w-cqv7vw9n>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Mar 10, 2023

Last Modified: Mar 17, 2023

PROTOCOL integer ID:
78495

1 Bacterial colony PCR

 GoTaq Green Master Mix **Promega Catalog #M7122**

PCR reactions using ProFlex PCR system, Applied Biosystems. [Thermal Cycler](#)

1.1





A	B	C
Reagent Name	Initial Conc.	Volume (μL)
gotaq green 2x	2X	5
oligo 1	10uM	0.6
oligo 2	10 uM	0.6
cell suspension in water	cells suspension	0.8
ddH2O	N/A	3
TOTAL Volume:	10 uL	

A	B	C	D
Number of Cycles	Step Name	Temperature	Duration
1	Initial Denaturation	98oC	10min
	Denaturation	98oC	10 sec
30	Annealing	50-55oC	30 sec
	1min (depends on the length)	72oC	1min per 1kb
1	Hold	4-10oC	hold

Note

Maximum (good quality) PCR product up to ~3kb. The best yield is 0.2kb-2kb


Note

pick a colony from agar plate and swirl in sterile water  50-100 μL . From liquid culture, dilute saturated culture  5-10 μL in sterile  50-100 μL water. This is a source for reaction (taking  0.8 μL).





Keep the rest in ice, for re-culture of positively verified clones.

1.2 Load on 1-1.2% agarose gel  8-10 µL and visualise under GelDoc or another UV source.

2 Yeast colony PCR (*S.cerevisiae*, *Pichia*, *Candida*)

2.1 Alkaline yeast cells lysis treatment. Pick a yeast colony from agar plate and swirl in sterile 20mM NaOH  50-100 µL . 15m

Optionally, a liquid culture or a agar-plate pick colony resuspend in the sterile water or YPD, take small volume (~20uL) and mix with the same volume of 40mM NaOH.

20mM NaOH-cells resuspension incubate  95 °C  00:10:00 and immediately transfer to ice/  4 °C for at least  00:05:00

2.2

A	B	C
Reagent Name	Initial Conc.	Volume (µL)
gotaq green 2x	2X	5
oligo 1	10uM	0.6
oligo 2	10 uM	0.6
cell suspension in 20mM	yeast suspension	0.8
ddH2O	N/A	3
TOTAL Volume:	10 uL	

A	B	C	D
A	B	C	D
Number of Cycles	Step Name	Temperatur e	Duration

A	B	C	D
1	Initial Denaturation	98oC	2min
30-33	Denaturation	98oC	10 sec
	Annealing	52-55oC	1min 30sec
	1min (depends on the length)	72oC	1min per 1kb
1	Hold	4-10oC	hold

Note

Maximum PCR product up to ~3kb. The best yield is 0.2kb-2kb

Note

Important - annealing 1min30sec

2.3 Load on 1-1.2% agarose gel  8-10 µL and visualise under GelDoc or another UV source.