



Oct 01, 2024

Cloning with Golden Gate

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.yxmvme22bg3p/v1

Raining Wang¹, Melinda Wheelock¹

¹University of Washington

IGVF



Sarah Heidl

University of Washington, Starita Lab

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.yxmvme22bg3p/v1**

Protocol Citation: Raining Wang, Melinda Wheelock 2024. Cloning with Golden Gate. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.yxmvme22bg3p/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: May 09, 2024

Last Modified: October 01, 2024

Protocol Integer ID: 99784

Abstract

This protocol details the cloning and transformation.

Materials

PCR mix:

A	B	C
Component	Final conc.	1Rxn (uL)
H2O		Up to 25ul
T4 DNA ligase buffer	1X	2.5
T4 DNA ligase (2,000U/ul)		0.5
Plasmid (100ng/ul)	2.5ng	2.5
SapI (10U/ul)	20U	2
Insert (ul)	1:2 mol ratio	1.09
Total		25

Digestion mix:

A	B	C	D	E	F
		Total samples #	2	Err	1.05
Component	Final conc.	1Rxn (uL)	Total		
H2O		3	6.3		
rCutSmart Buffer		1	2.1		
SapI (10U/ul)	10U	1	2.1		
Total		5	#REF!		

Wild-type target assembly

1d

1 Golden Gate assembly:

- If the wild-type insert is obtained from PCR amplification, ensure the PCR product is cleaned.


A	B	C
Component	Final conc.	1Rxn (uL)
H2O		Up to 25ul
T4 DNA ligase buffer	1X	2.5
T4 DNA ligase (2,000U/ul)		0.5
Plasmid (100ng/ul)	2.5ng	2.5
SapI (10U/ul)	20U	2
Insert (ul)	1:2 mol ratio	1.09
Total		25

1.1 A negative control reaction with H2O instead of insert is recommended.

1.2 Set the thermal cycler for the program below:


A	B	C
Step	Temperature	Time (hh:mm:ss)
1	37°C	00:05:00
2	16°C	00:05:00
4	Go to step 1 for 10-15 total cycles	
5	60°C	00:15:00
6	85°C	00:15:00
7	4°C	hold

Note

SapI enzyme is less thermostable. Consider adding  1 µL of SapI enzyme after 2 hours.

2 Made additional digestion mix.

		Total samples #	2	Err	1.05
Component	Final conc.	1Rxn (uL)	Total		
H2O		3	6.3		
rCutSmart Buffer		1	2.1		
SapI (10U/uL)	10U	1	2.1		
Total		5	#REF!		



2.1 Aliquot  5 µL of digestion mix to each of the golden gate reactions.

2.2 Set the thermal cycler for the program below:

A	B	C
Step	Temperature	Time (hh:mm:ss)
1	55°C	At least 2 hours
2	80°C	00:15:00
4	Hold at 4C	

3 Purify each reaction with Zymo Clean & Concentrate -5 kit.

4 Transformation:

Add  1 µL -  2 µL of DNA from sample reaction and negative control to chemical competent E. coli cells by following NEB stable competent cell transformation protocol.

Note

Quick transformation without an outgrowth period is adequate. It will reduce transformation efficiency.

4.1 Consider to add Puc 19 and H₂O as controls to the transformation process.




5 Incubate overnight at  30 °C for  24:00:00 .

1d



6 **Colony selection and mini prep:**

Check negative control plate to ensure the background colony is less than 10% of the colony counts from the sample plate.

7 Pick 2-3 colonies for an overnight growth in  3 mL of LB + amp (1:1000) per colony for mini prep.

8 Sequence confirmation by nanopore sequencing.