


PBMC Stimulation with
Peptide Pools and
Fluorospot Assay

Jul 19, 2021

PBMC Stimulation with Peptide Pools and Fluorospot Assay

cecilia¹, Alessandro Sette¹¹La Jolla Institute for Immunology

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.bphjmj4n Yaqian Xu

ABSTRACT

The purpose of this protocol is to stimulate PBMCs with peptide pools for 14 days followed by Fluorospot assay.

ATTACHMENTS

[cpvbbg4xf.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.bphjmj4n

PROTOCOL CITATION

cecilia , Alessandro Sette 2021. PBMC Stimulation with Peptide Pools and Fluorospot Assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bphjmj4n>

KEYWORDS

PBMC, peptide pools, Fluorospot assay

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited


CREATED

Nov 07, 2020

LAST MODIFIED

Jul 19, 2021

OWNERSHIP HISTORY

Nov 07, 2020  Megan FreundNov 11, 2020  Yaqian Xu

PROTOCOL INTEGER ID

44299

MATERIALS TEXT

- Peptide pools
- Cryopreserved PBMC
- Media (complete RPMI with 5% human serum)
- 50 mL tubes
- Water bath
- Centrifuge

 **Benzonase Sigma Catalog #70664-3**

 **Recombinant Human IL-2 Protein R&D**

Systems Catalog #202-IL










 **Trypan Blue Sigma**

Aldrich Catalog #T8154

SAFETY WARNINGS

For safety wans hazard warnings, please see the Safety Data Sheet (SDS).

Thawing of Cryopreserved PBMC

- 1 Add  **10 mL media** (complete RPMI with 5% human serum) and  **20 µl Benzonase** per vial of cells to a 50 ml tube.
- 2 Place PBMC vials in  **37 °C** water bath for approximately 60-90s.  **00:01:30** 1m 30s
- 3 Before cells are completely thawed, add them to the prepared 50 ml tube.
- 4 
Centrifuge at  **1200 rpm, 00:07:00** .
- 5 Pour off media and resuspend in media.
- 6 
Count cells and record viability (trypan blue).
- 7 
Centrifuge at  **1200 rpm, 00:07:00** .
- 8 Pour off media and resuspend according to cell count to a final concentration of 4×10^6 cells/ml.

9 Record peptide pool information:

Donor	<input type="text"/>
Location	<input type="text"/>
Vials	<input type="text"/>
Cells/ml ($\times 10^6$)	<input type="text"/>
Volume (mL)	<input type="text"/>
Total # PBMCs ($\times 10^6$)	<input type="text"/>
Volume to reach 4×10^6 /ml (mL)	<input type="text"/>
% viability	<input type="text"/>

In vitro Restimulation for 14 days (Plate cells (2×10^6 per well) in 24-well plate)

10 Resuspend thawed PBMC at 4×10^6 /ml.

11 Add mL cells per well and mL stimulus ($5 \mu\text{g}/\text{ml}$ final concentration) .

Plate 1	1	2	3	4	5	6
A						
B						
C						
D						

Plate Layout

12 Add media supplemented with 10 U/mL (final concentration) of rIL-2 every 3-4 days.

12.1 Day 4 – add mL media supplemented with IL-2.

12.2 Day 7 and 11 – replace  1 mL media supplemented with IL-2 without disturbing cells.

13