



APR 02, 2024

🌐 Multiplex Labeling with Tyramide Fluorophores (Free-Floating Tissues)-Killinger Lab 2024

Bryan Killinger¹, Solji_G_Choi¹, Tyler_Tittle¹

¹Rush University

ASAP Collaborative Research Network

Killinger



Bryan Killinger
Rush University

OPEN  ACCESS



ABSTRACT

This protocol details the multiplex labeling of free-floating tissues using tyramide fluorophores in the killinger lab (2024).

ATTACHMENTS

[2024 Dual Labeling Using Tyramide Fluorophores.docx](#)

DOI:

dx.doi.org/10.17504/protocols.io.yxmvme7zng3p/v1

Protocol Citation: Bryan Killinger, Solji_G_Choi, Tyler_Tittle 2024. Multiplex Labeling with Tyramide Fluorophores (Free-Floating Tissues)-Killinger Lab 2024.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.yxmvme7zng3p/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Mar 18, 2024

Last Modified: Apr 02, 2024

PROTOCOL integer ID: 97660

Funders Acknowledgement:
NIH-NINDS
Grant ID: 1R01NS128467
Michael J Fox Foundation
Grant ID: ASAP-024442

MATERIALS

Sodium Citrate Buffer, (1L):

A	B
2.94 g	Sodium citrate-Trisodium salt(Dihydrate)in 1000 mL DI water.
6.0 pH	Adjust pH to
0.5 mL	Tween-20 (Mix well)

Blocking buffer:

A	B
100 mL	Dilution media
3 mL	Normal serum
2 g	Bovine serum albumin
0.4 mL	Triton x100 (Mix well so the Triton is completely dissolved)

[M] 0.05 Molarity (M) **Borate buffer** 

A	B
300mL	DI H2O
5.72 g	Sodium tetraborate decahydrate (P17, big bottle)

Day 1:

1h 50m

1 Wash free-floating tissue (3 x 10 minutes) in dilution media (DM).



1.1 Wash free-floating tissue for  00:10:00 in dilution media (DM) (1/3).

10m



1.2 Wash free-floating tissue for 00:10:00 in dilution media (DM) (2/3).

10m



1.3 Wash free-floating tissue for 00:10:00 in dilution media (DM) (3/3).

10m



2 Heat water bath 01:30:00 before the antigen retrieval step.

1h 30m



1. Human samples: 90 °C - 95 °C

2. Mouse samples: 80 °C - 85 °C

3 Place the dish containing sodium citrate buffer in the water bath and heat it for 00:10:00 .

10m



a. **Sodium Citrate Buffer, pH 6 (1L):**

- 2.94 g Sodium citrate-Trisodium salt (Dihydrate) in 1000 mL DI water.
- Adjust pH to 6.0.
- 0.5 mL Tween-20. Mix well.

4 Wash the tissues in sodium citrate buffer for 00:05:00 .

5m



5 Incubate the tissues in the heated sodium citrate buffer for 00:30:00 .

30m



6 Cool down the tissues by placing container in an ice bucket for 00:15:00 .

15m

7 Wash in DM for 10 minutes x 2 times.



7.1 Wash in DM for 00:10:00 (1/2).

10m



7.2 Wash in DM for 00:10:00 (2/2).

10m



8 Endogenous peroxidase inhibition and serum blocking step (01:00:00 incubation): 0.3% H₂O₂+0.1% Sodium Azide in blocking buffer.

1h



a. Blocking buffer:

A	B
Dilution media	100 mL
Normal serum	3 mL
Bovine serum albumin	2 g
Triton x100 (Mix well so the Triton is completely dissolved)	0.4 mL

b. In 50 mL blocking buffer, add 0.5 mL 30% H₂O₂ + 0.5 mL 10% Sodium Azide.

9 Dilute primary antibody in blocking buffer. Incubate Overnight at 4 °C .


10m



Day 2:

10 Wash (3 x 10 minutes) in dilution media.



10.1 Wash for  00:10:00 in dilution media (1/3).


10m



10.2 Wash for  00:10:00 in dilution media (2/3).


10m



10.3 Wash for  00:10:00 in dilution media (3/3).

10m



11 HRP-Secondary antibody incubation 1:1000 dilution ( 01:00:00).


1h



▪ Solvent is  100 mL DM/  1 mL normal serum/  1 g BSA.

12 Wash (2 x 10 minutes) in dilution media.



12.1 Wash for  00:10:00 in dilution media (1/2).

10m



12.2 Wash for 00:10:00 in dilution media (2/2).

10m



13 Wash in borate buffer for 00:10:00 .

10m



a. [M] 0.05 Molarity (M) **Borate buffer** pH 8.5

A	B
300mL	DI H2O
5.72 g	Sodium tetraborate decahydrate (P17, big bottle)

1. Mix well to dissolve completely.

2. Adjust to pH 8.5 .

14 Incubate with tyramide fluorophore (TF) for 00:30:00 while blocking light.

30m



a. 10 mL Borate buffer + 1 µL H₂O₂ + 5 µL TF.

15 View under the microscope to confirm successful staining.

16 Store in PBS and leave at 4 °C . It can be stored for up to 2 weeks. Otherwise, proceed with the antigen retrieval step.



Day 3:

2h 50m

17 Heat water bath 01:30:00 before the antigen retrieval step.

1h 30m



a. Human samples: 90 °C - 95 °C

b. Mouse samples: 80 °C - 85 °C

18 Place the dish containing sodium citrate buffer in the water bath and heat it for 00:10:00 . 10m

a. **Sodium Citrate Buffer, pH 6 (1L):**

- 2.94 g Sodium citrate-Trisodium salt (Dihydrate) in 1000 mL DI water.
- Adjust pH to 6.0
- 0.5 mL Tween-20. Mix well.

19 Wash the tissues in sodium citrate buffer for 00:05:00 . 5m



20 Incubate the tissues in the heated sodium citrate buffer for 00:30:00 . 30m



21 Cool down the tissues by placing container in an ice bucket for 00:15:00 . 15m

22 Wash in DM for 10 minutes x 2 times.



22.1 Wash in DM for 00:10:00 (1/2). 10m



22.2 Wash in DM for 00:10:00 (2/2). 10m





23

Endogenous peroxidase inhibition and serum blocking step (🕒 00:10:00 incubation): 0.3% H_2O_2 +0.1% Sodium Azide in blocking buffer. 10m



a. Blocking buffer:

A	B
Dilution media	100 mL
Normal serum	3 mL
Bovine serum albumin	2 g
Triton x100 (Mix well so the Triton is completely dissolved)	0.4 mL

b. In 50 mL blocking buffer, add 🧪 0.5 mL 30% H_2O_2 + 🧪 0.5 mL 10% Sodium Azide.

24

Dilute primary antibody in blocking buffer. Incubate 🕒 Overnight at 🌡️ 4 °C . 10m



Day 4:

2h 20m

25

Wash (3 x 10 minutes) in dilution media.



25.1

Wash for 🕒 00:10:00 in dilution media (1/3). 10m



25.2

Wash for 🕒 00:10:00 in dilution media (2/3). 10m



25.3

Wash for 00:10:00 in dilution media (3/3).

10m



26

HRP-Secondary antibody incubation 1:1000 dilution (01:00:00).

1h



a. Solvent is 100 mL DM/1 mL normal serum/1g BSA.

27

Wash (2 x 10 minutes) in dilution media.

27.1

Wash for 00:10:00 in dilution media (1/2).

10m

27.2

Wash for 00:10:00 in dilution media (2/2).

10m

28

Wash in borate buffer for 00:10:00 .


10m




a. 0.05 Molarity (M) Borate buffer pH 8.5



A	B
DI H2O	300mL
Sodium tetraborate decahydrate (P17, big bottle)	5.72 g

1. It takes a while to dissolve completely.

2. Adjust to  8.5 .

29 Incubate with tyramide fluorophore (TF) on a shaker for  00:30:00 while blocking light.

30m

a.  10 mL Borate buffer +  1 μ L H_2O_2 +  5 μ L TF.


30 View under the microscope to confirm successful staining.

31 DAPI staining ( 00:20:00)

20m

a. 1:2000 dilution PBS. Block the light.

32 Mount the tissues on a slide, cover the slide with Fluoroshield, and coverslip. Seal with nail polish on all sides of the coverslip.

33 When the nail polish is completely dried, view under the microscope. Always protect the slides from light. Slides can be stored at  4 $^{\circ}$ C .