



DEC 06, 2022

WORKS FOR ME

1

CTAB-based DNA extraction for citrus

DOI

dx.doi.org/10.17504/protocols.io.dm6gpj8jpgzp/v1Shingo Goto¹¹Division of Citrus Research, Institute of Fruit Tree and Tea Science, NARO

Shingo Goto

[Division of Citrus Research, Institute of Fruit Tree and Tea...](#)

COMMENTS 0

ABSTRACT

CTAB-based protocols are used for genomic DNA extraction from many kinds of plant species. However, the protocols can't necessarily completely remove contamination of polysaccharide and RNA in extracted genomic DNA solution. Especially, citrus leaves generally contain high polysaccharide. This protocol is a simple and efficient method for extracting genomic DNA for citrus without contamination of polysaccharide and RNA.

DOI

dx.doi.org/10.17504/protocols.io.dm6gpj8jpgzp/v1

PROTOCOL CITATION

Shingo Goto 2022. CTAB-based DNA extraction for citrus. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.dm6gpj8jpgzp/v1>

KEYWORDS

CTAB, DNA extraction, citrus, High-salt precipitation solution

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 12, 2022

LAST MODIFIED

Dec 06, 2022

PROTOCOL INTEGER ID

69837

Buffer preparation

1












- 2×CTAB solution: 2% (w/v) CTAB, 100mM Tris-HCl pH8.0, 1.4M NaCl, 20mM EDTA pH8.0 [1]
- High-salt precipitation solution: 1.2M NaCl, 0.8M Sodium citrate [2]

- 10 mg/ml RNase (Nippon gene)
- Chloroform
- Isopropanol
- 70% Ethanol
- TE buffer: 10 mM Tris-HCl, 1 mM EDTA, pH 8.0




Reference

1. Allen GC, Flores-Vergara MA, Krasynanski S et al. A modified protocol for rapid DNA isolation from plant tissues using cetyltrimethylammonium bromide. Nat Protoc 2006;1:2320–5.
2. Chomczynski P, Mackey K. Modification of the TRI Reagent procedure for isolation of RNA from polysaccharide- and proteoglycan-rich sources. Biotechniques 1995;19:492–5.

Homogenization and cell lysis

- 2 Preheat 2×CTAB solution to  60 °C in water bath. Add  2 % (v/v) of 2-Mercaptoethanol to the 2×CTAB solution just before use.
- 3 Homogenize  100 mg of fresh leaf in liquid Nitrogen. Add  800 µL of 2×CTAB solution and completely suspend homogenized leaf. Transfer the suspended solution to 2 ml tube.
- 4 Add  4 µL of  10 mg/mL RNase, mix by inversion, and incubate at  37 °C for  00:15:00 (In order to prevent RNase contamination in laboratory, RNase treatment is conducted before denaturing proteins by chloroform).
- 5 Incubate at  56 °C for  00:30:00 inverting the tube once every  00:10:00 .



Chloroform extraction

- 6 Add  300 µL of Chloroform and mix gently with tube rotator for  00:15:00 .
- 7  13000 rpm, 25°C, 00:10:00

8 Transfer supernatant carefully to new 2 ml tube.

9 Repeat step 6 and 7


Precipitation and wash of DNA pellet

10 Transfer  600 µL of supernatant to new 1.5 ml tube, add  300 µL of High-salt precipitation solution, and mix by inversion.

11 Add  300 µL of Isopropanol and mix by inversion.

12  15000 rpm, 25°C, 00:10:00

13 Discard and remove supernatant (DNA pellet is often transparent).

14 Add  1000 µL of 70% ethanol and mix by inversion 10 times to wash salts.


15  15000 rpm, 25°C, 00:05:00

16 Completely remove supernatant and dry up DNA pellet in air.

17 Dissolve pellet in \sphericaltriangle 100 μ L of TE buffer.

Result of gel electrophoresis

18





0.8% agarose gel electrophoresis

Lane 1-3: Genomic DNA extracted by isopropanol precipitation .

Lane 4: size marker

Lane 4-6: Genomic DNA extracted by isopropanol precipitation with high salt precipitation solution.

All genomic DNA is a citrus cultivar, 'Kiyomi'.