



APR 12, 2024

## 🌐 CODEX Antibody Conjugation Protocol

Bei Wei<sup>1</sup>, Joanna Bi<sup>1</sup>

<sup>1</sup>Stanford University

Human BioMolecular Atlas Program (HuBMAP) Method Development Community

Tech. support email: [Jeff.spraggins@vanderbilt.edu](mailto:Jeff.spraggins@vanderbilt.edu)



Joanna Bi  
Stanford University

OPEN  ACCESS



### ABSTRACT

Used by the Snyder Lab at Stanford University for CODEX for human intestine.

#### DOI:

[dx.doi.org/10.17504/protocols.io.4r3l2q8q4l1y/v1](https://dx.doi.org/10.17504/protocols.io.4r3l2q8q4l1y/v1)

**Protocol Citation:** Bei Wei, Joanna Bi 2024. CODEX Antibody Conjugation Protocol. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.4r3l2q8q4l1y/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Apr 12, 2024

Last Modified: Apr 12, 2024

PROTOCOL integer ID: 98091

- 1 Block non-specific absorption of antibody onto 50k MWCO spin column with PBS-tween (0.2% tween solution) by adding 500 uL into the top of each column (rinse all areas of filter) and spinning down at 12,000g for 2 minutes
- 2 Remove all liquid from the top of the column and discard flow-through. Label each column and filter.
- 3 Spin down the antibodies
- 4 Optional: measure the concentration of your antibodies
  - 4.1 Nanodrop, protein A280, blank with PBS, IgG
- 5 Add 100 ug of antibody to the top of each pre-blocked 50k MWCO filter columns. If the volume is less than 100 uL, add 200 uL of PBS to the column before adding the antibody
- 6 Spin column down at 12,000g for 8 minutes. Discard flow-through.

**6.1** If the volume of antibody needed did not fit, repeat this step until all of the antibody is contained in the filter unit.

**7** Optional: If antibody needs to be washed (Azide/BSA) add 400 uL more of PBS and spin down again at 12,000g for 8 minutes. Repeat as necessary

**8** Create the following solution and vortex to mix: 5 uL 500mM TCEP, 5 uL 500mM EDTA, 990 uL 1X PBS

**9** Add 360 uL of this solution to the top of the MWCO filter unit. Pipette mix then briefly vortex the tube

**10** Incubate at RT for 30 minutes (do not go longer than 30 minutes)

**11** Spin the solution at 12,000g for 8 minutes. Discard flow-through

**12** Add 400 uL Buffer C and spin the solution at 12,000g for 8 minutes. Discard flow-through.

**12.1** Remove oligo aliquots from the freezer at this point

- 13**      Make high salt buffer C: 380 uL Buffer C, 20 uL 5M NaCl
  
- 14**      Solubilize oligos. Oligo:antibody ratio is 2:1. For 100 ug antibody use 200 ug oligo. Do one tube at a time. Once oligo is solubilized the next step should follow immediately. Add high salt buffer C to the tube containing the top oligo. Mix with pipette until oligo is completely solubilized. There should be no visible cloudiness to the solution.
  - 14.1**      PCR tubes (100 ug/tube): use 2 PCR tubes, add 200 uL high salt buffer C to each
  
  - 14.2**      Large tube (lyophilized) (200 ug/tube): add 400 uL high salt buffer C
  
  - 14.3**      Large tube (prediluted) (200 ug/tube): add high salt buffer C to a total of 400 uL
  
- 15**      Add entire oligo solution to the top of the filter column containing the the reduced antibody solution.
  
- 16**      Mix well with the pipette. Do not vortex.
  
- 17**      Incubate the oligo-antibody mix for 2h at room temperature (can go longer than 2 hours)

- 18** Spin the column down at 12,00g for 8 minutes. Discard flow-through.
- 19** Add 450 uL of high salt PBS to the top of the column and spin down at 12,000g for 8 minutes. Discard flow-through.
  - 19.1** High salt PBS (1M NaCl PBS): 5 M NaCl diluted to 1M in PBS
- 20** Add another 450 uL of high salt PBS to the top of the column and spin down at 12,000g for 8 minutes. Discard flow-through.
- 21** Repeat step #20 one more time.
- 22** Add CODEX antibody stabilizer solution to the top of the column in a 2:1 ratio (2 uL solution : 1 ug antibody)
  - 22.1** CODEX antibody stabilizer (can make stock and store at 4C for a few weeks): 900 uL antibody stabilizer (BOCA scientific, + 0.2% sodium azide), 10 uL 0.5M EDTA, 100 uL 5M NaCl
- 23** Invert into new collection tube. Spin down at 3,000g for 2 minutes. KEEP COLLECTED SOLUTION

- 24** Combine total volume of conjugated antibody for a given antibody into a single screw-top tube and store at 4C – DO NOT FREEZE