



SEP 20, 2023

🌐 Dual In Situ Hybridization/Immunofluorescence


📁 In 1 collection

Michael Henderson¹

¹Van Andel Institute

ASAP Collaborative Research Network

Team Biederer



maria.matos

OPEN  ACCESS



ABSTRACT

This protocol details about the Dual In Situ Hybridization/Immunofluorescence in tissue.

ATTACHMENTS

[338-741.pdf](#)

GUIDELINES

ACD protocol notes that tissues should be fixed in 10% NBF for 16-32 hours, and embedded in paraffin. Then, sectioned and dried overnight at RT. They suggest sectioned tissue be used in less than a year (4°C) or less than 3 months at room temperature.

Protocol Citation: Michael Henderson 2023. Dual In Situ Hybridization/Immunofluorescence. **protocols.io** <https://protocols.io/view/dual-in-situ-hybridization-immunofluorescence-b3g4qjyw>

MANUSCRIPT CITATION:
Adapted from ACD Standard Protocol/Cheadle/Otero-Garcia Protocols

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jan 05, 2022

Last Modified: Sep 20, 2023

MATERIALS

Solutions

A	B	C
Needed (mL)	Stock Solution	Final Concentration
5 L	dH2O	
485 g	Tris base	0.5 M
240 mL	Concentrated HCl	
pH to 7.6		
To 8L	dH2O	

0.5 M Tris (8 L)

Reagents

A	B	C	D	E
Vendor	Catalog #	Qty	Unit Price	Description
RNAscope® Multiplex Fluorescent Reagent Kit V2	323100		1330	Contains H2O2, protease reagents, target retrieval reagent, wash buffer, HRP reagents
RNAscope® 3-plex	320861		100	Polr2a (C1 channel) and PPIB (C2 channel), UBC (C3 channel)
Positive Control Human Sigma	199664-25G	1	66.6	Sudan Black B
Vector Laboratories	H-4000	1	120	ImmEdge Hydrophobic Barrier Pen
Southern Biotech	0100-01	1	45.14	DAPI Fluoromount-G



Sudan black B Sigma
Aldrich Catalog #199664



ImmEdge hydrophobic barrier pap pen Vector
Laboratories Catalog #H-4000



Fluoromount-G Southern
Biotech Catalog #0100-01



Preparing Tissue (Day 1): Prepare Tissue

1h 19m

- 1 Bake slides in a dry oven for  01:00:00 at  60 °C . Use slides within a week. 1h
- 2 De-paraffinize slides in fresh xylenes, then in 100% ethanol.
 - 2.1 De-paraffinize slides for  00:05:00 in fresh xylenes. (1/4) 5m
 - 2.2 De-paraffinize slides for  00:05:00 in fresh xylenes. (2/4) 5m
 - 2.3 De-paraffinize slides for  00:02:00 in 100% ethanol. (3/4) 2m
 - 2.4 De-paraffinize slides for  00:02:00 in 100% ethanol. (4/4) 2m
- 3 Place slides on absorbent paper and dry in the oven from  00:05:00 at  60 °C or until dry. 5m

Preparing Tissue (Day 1): Hydrogen Peroxide Treatment

1h 19m

- 4 Place slide horizontally in an incubation tray. Add ~5-8 drops of RNAscope Hydrogen Peroxide to cover each section. Incubate for  00:10:00 at  Room temperature .



10m



- 5 Dab solution off and move to a rack in distilled water. Move up and down 5 times. Repeat with a fresh boat of distilled water.

Preparing Tissue (Day 1): Target Retrieval


1h 19m

- 6 Dilute Target Retrieval Regents (RNAscope) 1:10 in dH₂O ( 25 mL /  225 mL dH₂O/boat). Mix well.






- 7 Place in microwave for  00:15:00 at  95 °C . 15m

- 8 Transfer slides to a slide boat with  200 mL distilled water for  00:00:15 . 15s

- 9 Transfer the slides to 100% ethanol for  00:03:00 . 3m



- 10 Dry the slides in a  60 °C incubator (or  Room temperature) for  00:05:00 . 5m

- 11 Draw a hydrophobic barrier onto slides with ImmEdge pen. Do NOT due for fluorescent slides. Let the barrier dry for  00:05:00 . OPTIONAL PAUSE POINT  Overnight at  Room temperature .



RNAscope Multiplex Fluorescent v2 Assay (Day 2): Protease T...

1h 15m

- 12 Place a wet Humidifying Paper in an incubation tray and warm for  00:30:00 at  40 °C (TC incubator). Keep the tray in the incubator when not in use. Insert the slides into the



incubation tray.

- 13** Add ~5 drops RNAScope Protease Plus (Protease III-Cheadle) to cover each section and place tray into the incubator at 40 °C for 00:30:00 (standard) (00:15:00 -Otero-Garcia). 45m



Note

Prepare RNAScope assay reagents during this step.

- 14** Wash slides with 200 mL+ distilled water and slight agitation.

- 14.1** Wash slides with 200 mL + distilled water and slight agitation. (1/2)



- 14.2** Wash slides with 200 mL + distilled water and slight agitation. (2/2)



RNAScope Multiplex Fluorescent v2 Assay (Day 2): Preparation 1h 15m

- 15** **Wash Buffer:** Warm 50x Buffer to 40 °C for 00:10:00 to 00:20:00. Add 980 mL distilled water to 20 mL of RNAScope Wash Buffer in a 1 L bottle. May need 1 L to 2 L per run. Mix well. Can be stored for up to one month. 30m





- 16** **Probes:** Prepare only those probes needed. 10m











Note

If you are only using C2 and C3, dilute in probe diluent instead of C1.

Warm probes for 00:10:00 at 40 °C, then let cool to Room temperature. Add 1 volume C2 and volume C3 probes to 50 volumes C1 probe in a tube (e.g. 200 µL C1 +

 4 μL C2). Invert to mix. Store at  4 °C for up to 6 months.

17 **Reagents:** Warm AMP1-3, HRP-C1-3 and HRP blocks at  Room temperature .

18 (Optional) **Saline Sodium Citrate:**  175.3 g NaCl +  88.2 g sodium citrate in  800 mL distilled water. Adjust to  7.0 with  1 Molarity (M) HCl. Add water to a final volume of  1 L . Sterilize by autoclaving and store at  Room temperature for up to 2 months.

RNAscope Multiplex Fluorescent v2 Assay (Day 2): Hybridize P.. 1h 15m

19 Remove liquid from slides. Add 4-6 drops (6 drops=  180 μL) of the probe mix to slides. Incubate in incubator for  02:00:00 at  40 °C . 2h





20 Wash slides with Wash Buffer.

20.1 Wash slides with Wash Buffer  00:02:00 at  Room temperature . (1/2) 2m



20.2 Wash slides with Wash Buffer  00:02:00 at  Room temperature . (2/2) 2m





21 **OPTIONAL PAUSE POINT:** Store slides in 5x SSC  Overnight at  Room temperature . 2m



RNAscope Multiplex Fluorescent v2 Assay (Day 2): Hybridize A..

1h 15m

- 22** Remove liquid from slides. Add 4-6 drops RNAscope Multiplex FL v2 Amp 1 to each slide. Incubate in incubator for  00:30:00 at  40 °C .

30m



- 23** Wash slides with Wash Buffer.

- 23.1** Wash slides with Wash Buffer  00:02:00 at  Room temperature . (1/2)



2m



- 23.2** Wash slides with Wash Buffer  00:02:00 at  Room temperature . (2/2)

2m





- 24** Repeat steps 22 and 23 for Amp 2 and Amp 3. Amp 3 only requires  00:15:00 at  40 °C .

15m

- 25** During this incubation, dilute necessary Opal Dye fluorophores in TSA Buffer (1:1500 standard).

RNAscope Multiplex Fluorescent v2 Assay (Day 2): Develop HR..

1h 15m

- 26** Remove liquid from slides. Add 4-6 drops RNAscope Multiplex FL v2 HRP-C1 to each slide. Incubate in incubator for  00:15:00 at  40 °C .

15m



- 27** Wash slides with Wash Buffer.

27.1 Wash slides with Wash Buffer  00:02:00 at  Room temperature . (1/2)



2m



27.2 Wash slides with Wash Buffer  00:02:00 at  Room temperature . (2/2)

2m



28 Remove liquid from slides. Add  200 μ L Opal 520 to each slide. Incubate in HybEZ Oven for  00:30:00 at  40 °C .

30m



29 Wash slides with Wash Buffer.

29.1 Wash slides with Wash Buffer  00:02:00 at  Room temperature . (1/2)



2m



29.2 Wash slides with Wash Buffer  00:02:00 at  Room temperature . (2/2)

2m



30 Remove liquid from slides. Add 4-6 drops RNAScope Multiplex FL v2 HRP Blocker to each slide. Incubate in incubator for  00:15:00 at  40 °C .

15m



31 Wash slides with Wash Buffer.

31.1 Wash slides with Wash Buffer  00:02:00 at  Room temperature . (1/2)

2m





31.2 Wash slides with Wash Buffer  00:02:00 at  Room temperature . (2/2)

2m



32 STOP HERE IF USING JUST C1 PROBE. Continue to **Immunofluorescence**.

33 Repeat steps 26-32 with HRP-C2 and Opal 570, and again with HRP-C3 and Opal 690.

Note

*Note that after additional of fluorophores, slides should be kept out of the light as much as possible.


33.1 4% PFA fix for  00:15:00 at  4 °C .

15m

33.2 Then, wash with PBS-Otero-Garcia for  00:04:00 . (1/2)

4m




33.3 Wash with PBS-Otero-Garcia for  00:04:00 . (2/2)

4m



Day 2: Immunofluorescence

1h 15m

34 Wash in  0.1 Molarity (M) Tris buffer,  7.6  00:05:00 . Discard all Tris washes.

5m



35 Block in **1mM 0.1 Molarity (M)** Tris/2% FBS (Tris/FBS) **00:30:00** +. Keep blocking solution for up to 2 weeks @ **4 °C**. **30m**

36 Dilute primary antibodies in Tris/FBS), and prepare humidified chamber(s) by soaking towel in the middle of the slide chamber(s).

37 Wipe excess fluid off back of slides and from around tissue and apply **200 µL** of primary antibody to slides.

38 Incubate at **4 °C** in humidified chamber **00:45:00** to **02:00:00** at **Room temperature** or **Overnight** at **4 °C**. **4h 45m**



Day 3

1h 15m

39 Rinse off antibody from tissue using Tris.

Note

Carefully direct spray from wash bottle around tissue, NOT directly on it.

40 Wash in Tris **00:05:00**. **5m**



41 Block in Tris/FBS **00:05:00**. **5m**

42 Dilute fluorophore-conjugated secondary antibody 1:500 in Tris/FBS and apply **200 µL** to **4h**





wiped slides. Incubate at Room temperature for 02:00:00 or Overnight at 4 °C .

43 Rinse off slides with Tris.

44 Wash in running tap H₂O for 00:05:00 .



5m

45 Wash in Tris for 00:05:00 in green boats.



5m

46 Coverslip using non-photobleaching reagent (Prolong Gold with DAPI or FluorMount with DAPI). Allow to dry completely before imaging on scanner.