

FEB 07, 2024

🌐 Anti-condensation agent

John Bergqvist¹

¹Behavioural Phenomics, MRC Laboratory of Medical Sciences

Behavioural Genomics



John Bergqvist

MRC Laboratory of Medical Sciences, Imperial College London

ABSTRACT

Preventing plate lids from fogging up by applying an anti-condensation agent. This allows longer recordings of e.g. worm behaviour.

MATERIALS

50 mL test tube
Sterile water
99% ethanol
Triton X-100

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.e6nvwdrq2lmk/v1

Protocol Citation: John Bergqvist 2024. Anti-condensation agent. [protocols.io](https://dx.doi.org/10.17504/protocols.io.e6nvwdrq2lmk/v1)
<https://dx.doi.org/10.17504/protocols.io.e6nvwdrq2lmk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Feb 07, 2024

PROTOCOL integer ID: 94838

Keywords: condensation, 96-well plates, 24-well plates, 6-well plates



Triton X-100

Pictograms:

- GHS05, GHS07, GHS09

Hazard statements:

- H302 Harmful if swallowed
- H318 Causes serious eye damage
- H411 Toxic to aquatic life with long lasting effects

Ethanol

Pictograms:

- GHS02, GHS07

Hazard statements:

- H226 Flammable liquid and vapour.
- H319 Causes serious eye irritation.

BEFORE START INSTRUCTIONS

Make sure to treat lids that are for use at least a day after treatment.

Making anti-condensation agent

- 1 Make a mixture of 0.05% Triton X-100 and 20% ethanol in a 50 mL test tube.

For 20 mL, mix:

🧪 4 mL 99% ethanol




🧪 0.1 mL Triton X-100

🧪 16 mL sterile water

1h 0m 30s

Treating plate lids

- 2 Treating lids with the anti-condensation agent **should be done at least a day before using the plates** to study worm behaviour.

3 In a Class II Microbiological Safety Hood, pour ca.  3-4 mL of the anti-condensation agent on the inside  30s of the lids and leave for  00:00:30 seconds.

Ensure all of the lid is covered by the anti-condensation agent.

4 Pour the liquid back into the 50 mL tube for later use, or onto new lids to be treated.

5 Leave lids to dry for ca  01:00:00 hour  1h.

6 After drying, you will most likely see rings of anti-condensation agent residue on the plate. Wipe these off the lids with Kimberly-Clark Kimwipes until there are no signs of these rings.