




Sep 12, 2022

MagAttract + Metapolyzyme metagenomic gDNA extraction from skin swabs

Natalie Ring¹¹Roslin Institute & Royal (Dick) School of Veterinary Studies, University of Edinburgh

1 Works for me

 Share

This protocol is published without a DOI.

Dogstails Natalie Ring

ABSTRACT

A protocol for the metagenomic extraction of bacterial DNA from skin swab samples (optimised using canine swabs), for use in a rapid diagnostics pipeline. At the end of the protocol, the DNA is cleaned up and ready for rapid barcoding (SQK-RBK004) library preparation for nanopore sequencing (or whatever other application you want to do).

Unless otherwise stated, all reagents should be included in the listed kits.

PROTOCOL CITATION

Natalie Ring 2022. MagAttract + Metapolyzyme metagenomic gDNA extraction from skin swabs. **protocols.io**
<https://protocols.io/view/magattract-metapolyzyme-metagenomic-gdna-extractio-cggxttn>



FUNDERS ACKNOWLEDGEMENT

Dogs Trust
Grant ID: Dogstails

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Coming soon.

LICENSE

_____ This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 12, 2022

LAST MODIFIED

Sep 12, 2022

PROTOCOL INTEGER ID

69879

GUIDELINES

This protocol, an adaptation of Qiagen's MagAttract HMW DNA kit, was developed by Natalie Ring and Alison Low for the Dogstails project, a collaboration between the Roslin Institute and the Royal (Dick) School of Veterinary Studies funded by the Dogs Trust. We are grateful to the dogs (and their owners) who donated samples to the R(D)SVS's Hospital for Small Animals, many of which were used in the development of this protocol.

Please follow on Twitter for latest updates, papers and results:

@NatalieAnneRing

MATERIALS TEXT

Kits

 [MagAttract HMW DNA](#)

[kit Qiagen Catalog #67563](#)


 [ProNex Size-Selective Purification](#)

[System Promega Catalog #NG2001](#)

 [Qubit® dsDNA HS Assay Kit Thermo Fisher](#)

[Scientific Catalog #Q32854](#) In 2 steps

Other reagents

- 50 mM Tris, 10 mM EDTA, pH8.0 ("buffer P1")
-  [1X PBS \(Phosphate-buffered saline\)](#)

[MetaPolyzyme](#) Sigma

- Aldrich Catalog #MAC4L-5MG
- [Nuclease-free Water](#) Contributed by users
- [Distilled Water](#) Contributed by users

Equipment

[Swabs in tubes \(no liquid\)](#) VWR

Scientific Catalog # 710-0928

Any swab tube with no transport/storage medium is fine.

DNA LoBind tubes, 1.5 mL

Tubes

Eppendorf 022431021 [Link](#)

1.5 mL

OR

SafeSeal reaction tube, 1.5 ml, PP, PCR

Performance Tested, Low DNA-binding

Tubes

Sarstedt 72.706.700 [Link](#)

1.5 mL

Magnetic Stand

Magnetic Stand

Thermo Scientific

MR02



Any magnetic rack that fits your tubes will suffice.



Centrifuge

Benchtop Centrifuge

Eppendorf

5405000441



Any benchtop centrifuge will suffice



ThermoMixer

Benchtop Incubator

Eppendorf

5382000023



Any heat block will suffice



Mini-centrifuge Centrifuge

Fisher S67601B 

Any standard mini centrifuge with adapters for different tube sizes will suffice



Vortex Mixer Vortex Mixer

VWR 97043-562 



BEFORE STARTING

- **"Buffer P1"** is required for the metapolyzyme lysis incubation: 50 mM Tris, 10 mM EDTA, pH 8.0
- **Metapolyzyme** is used here at a concentration of 3.3 mg/ml (resuspend 5 mg lyophilized powder in 1.5 ml PBS pH 7.5)
- We recommend using **low DNA-binding tubes** throughout, but definitely for the elution/storage of DNA

Extended pre-lysis spin down

- 1 Bathe swab tip in 3 ml PBS in the swab tube for 10 minutes, with occasional vortexing. 10m
Remove swab from tube, squeezing the sides as you do.

 **3 mL PBS**

 **00:10:00**

Bathing the swab tip in PBS overnight yields much more DNA, if you have time.

- 2 Pellet 2x 1.5 ml aliquots of cell-PBS solution in 1.5 ml tubes by centrifuging at maximum ^{20m} speed (13,000 RPM) for 20 minutes, then discard supernatant

 **3 mL PBS**

 **16,000 x g, Room temperature, 00:20:00**

We have found that this extended spin at the beginning of the protocol results in much better yield of bacterial gDNA, especially in samples with low bacterial abundance

Metapolyzyme & Proteinase K Lysis

- 3 Resuspend cell pellets (which might be invisible) and combine in 160 µl buffer P1 (50 mM Tris, 10 mM EDTA, pH 8.0)

 **160 µL buffer P1**

- 4 Add 20 µl metapolyzyme (3.3 mg/ml, 5 mg resuspended in 1500 µl PBS) and mix by flicking the tube

 **20 µL metapolyzyme (3.3 mg/ml)**

- 5 Incubate on a thermomixer for 60 minutes at 37°C with 900 RPM shaking

1h

 **900 rpm, 37°C, 01:00:00**

- 6 Add 20 µl MagAttract proteinase K and mix by flicking the tube

 **20 µL proteinase K**

- 7 Incubate on a thermomixer for 30 minutes at 56°C with 900 RPM shaking

30m

 **900 rpm, 56°C, 00:30:00**

MagAttract DNA isolation and washing

- 8 Add 150 µl MagAttract buffer AL and mix by pulse vortexing

 **150 µL buffer AL**

Our standard "pulse vortex" is 10 short (<1 second) pulses per tube

- 9 Add 15 µl MagAttract Suspension G and 280 µl MagAttract buffer MB and mix by pulse vortexing

 **15 µL Suspension G**

 **280 µL Buffer MB**

Make sure the magnetic beads (Suspension G) are really well mixed before adding them! The whole suspension should be black, not separated into a bead layer and a clear layer. We usually resuspended by vortexing for 10 or more seconds.

- 10 Spin down briefly, then pellet beads on magnet and remove supernatant

- 11 10. Add 700 µl MagAttract buffer MW1 and incubate on a thermomixer for 1 minute at room^{1m} temperature with 1,400 RPM shaking

 **700 µL buffer MW1**

 **1400 rpm, Room temperature , 00:01:00**

- 12 11. Repeat steps 9 and 10 1m

- 13 12. Spin down briefly, then pellet beads on magnet and remove supernatant

- 14 13. Add 700 µl MagAttract buffer PE and incubate on a thermomixer for 1 minute at room^{1m} temperature with 1,400 RPM shaking

 **700 µL buffer PE**

 **1400 rpm, Room temperature , 00:01:00**

- 15 14. Repeat steps 9 and 10 1m

- 16 15. Spin down briefly, then pellet beads on magnet and remove supernatant

- 17 16. Rinse the pelleted beads on the magnetic rack with 700 µl distilled water by pipetting down the opposite wall of the tube

 **700 µL distilled water**

18 17. Repeat step 16

19 18. Add 50 µl nuclease-free water off the magnet, to resuspend the bead pellet

🧴 50 µL nuclease-free water

20 19. Incubate on a thermomixer for 3 minutes at room temperature with 1,400 RPM shaking^{3m}

🌀 1400 rpm, Room temperature , 00:03:00

21 20. Spin down briefly, then pellet beads on magnetic rack and **keep supernatant** in a low-DNA binding 1.5 mL tube (e.g. [Eppendorf](#) or [Sarstedt](#))

Qubit Pre-clean-up quantification

22 Quantify DNA using Qubit dsDNA HS kit. If DNA concentration is an appropriate concentration for your experiment (for us, this means at least 0.2 ng/µl), continue to clean-up steps.

🔗 [Qubit® dsDNA HS Assay Kit Thermo Fisher](#)

Scientific Catalog #Q32854

🧴 1 µL DNA

🧴 199 µL Qubit dsDNA HS working solution

ProNex DNA clean-up

23 Add 150 µl room temperature ProNex beads to your entire tube of DNA (49 µl)

🧴 200 µL ProNex beads

Like the magnetic beads in Suspension G, make sure the ProNex beads are really well mixed (10+ seconds of vortexing) immediately before you use them.

- 24 Mix well by slowly pipetting up and down 10 times
- 25 Incubate at room temperature for 10 minutes (no shaking needed) 10m
- 🕒 00:10:00
- 🌡 Room temperature
- 26 Spin down briefly, then pellet beads on magnet and remove supernatant
- 27 Rinse the pelleted beads on the magnetic rack by pipetting 200 µl ProNex Wash Buffer down ^{1m}
the opposite wall of the tube, then incubate at room temperature for 60 seconds (no shaking), then remove Wash Buffer
- 📄 200 µL Wash Buffer
- 🌡 Room temperature
- 🕒 00:01:00
- 28 Repeat step 26
- 29 **Air-dry** (lid open) the sample on the magnetic rack for 5 minutes (longer is OK, no more than ^{5m}
60 minutes)
- 🌡 Room temperature
- 🕒 00:05:00
- 30 Add 20 µl nuclease-free water off the magnet. Resuspend the pellet by **flicking the tube**, ^{5m}
then incubate at room temperature for 5 minutes (no shaking needed)
- 📄 20 µL nuclease-free water
- 🌡 Room temperature
- 🕒 00:05:00

- 31 Spin down briefly, then pellet the beads on magnet and **keep supernatant** in a low DNA-binding tube

Qubit post-clean-up quantification

- 32 Quantify DNA using Qubit dsDNA HS kit. If DNA concentration is an appropriate concentration for your experiment (for us, this means at least 0.2 ng/μl), continue to library preparation.

 **Qubit® dsDNA HS Assay Kit Thermo Fisher**

Scientific Catalog #Q32854

 **1 μL DNA**

 **199 μL Qubit dsDNA HS working solution**