



Apr 29, 2022

🌐 Thawing of CT-2A cell line

Guido Krähenbühl¹¹MDC Berlin

1



protocol .

AG Gerhardt



Guido Krähenbühl

This protocol describes the thawing of mouse CT-2A glioma cells

Guido Krähenbühl 2022. Thawing of CT-2A cell line. **protocols.io**
<https://protocols.io/view/thawing-of-ct-2a-cell-line-b8gurtww>



protocol ,

Apr 29, 2022

Apr 29, 2022

61684

GMO safety level 1

 **DMEM Thermo Fisher****Scientific Catalog #41966** **FBS Invitrogen - Thermo Fisher**Make sure to warm up media to  **37 °C**

- 1 Remove cryovial from liquid nitrogen storage and place on dry-ice for transfer to cell culture lab.

Cryovial IDs	<input type="text"/>
Additional Identifiers	<input type="text"/>

- 2 Remove cryovial from dry-ice and hold in a **37 °C** waterbath until thawed



Avoid submerging the lid or rim of the cryovial in the waterbath as this may result in contamination

Dry the cryovial thoroughly, spray with 70% ethanol and transfer to cell culture hood

- 3 Transfer contents of cryovial to a [Falcon 15 mL Polystyrene Conical Tube Fisher Scientific Catalog #352095](#)

containing **7 mL DMEM, 10% FBS, no antibiotics** .

- 3.1 Rinse the cryovial with **1 mL** of the cell suspension and add it back to the same

[Falcon 15 mL Polystyrene Conical Tube Fisher Scientific Catalog #352095](#)

4



5m

Centrifuge at **300 rcf, Room temperature, 00:05:00**

- 5 Discard supernatant, resubstitute with fresh **7 mL DMEM, 10% FBS, no antibiotics**

- 6 Transfer cell suspension to a

 Petri dish, 10cm, polystyrene **Fisher**

Scientific Catalog #FB0875712


and agitate

gently back and forth and side to side to evenly distribute cells

7 Place the petri dish in a  **37 °C Incubator** (5%CO₂)

8 After 24 hours, inspect adherence and confluency.



If a major amount of cells did not attach, remove the media and replace with
 **7 mL DMEM, 10% FBS, no antibiotics** and check for surviving cells after 24 hours