

Jul 11, 2024

## Tile Scan Imaging and Cell Counting

DOI

**[dx.doi.org/10.17504/protocols.io.bp2l6231kgqe/v1](https://dx.doi.org/10.17504/protocols.io.bp2l6231kgqe/v1)**

Shiyi Wang<sup>1</sup>

<sup>1</sup>Duke University

ASAP Collaborative Rese...



Shiyi Wang

Duke University

---

OPEN  ACCESS



DOI: **[dx.doi.org/10.17504/protocols.io.bp2l6231kgqe/v1](https://dx.doi.org/10.17504/protocols.io.bp2l6231kgqe/v1)**

**Protocol Citation:** Shiyi Wang 2024. Tile Scan Imaging and Cell Counting. **protocols.io**

**<https://dx.doi.org/10.17504/protocols.io.bp2l6231kgqe/v1>**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 11, 2024

**Last Modified:** July 11, 2024

**Protocol Integer ID:** 103200

**Keywords:** ASAPCRN

**Funders Acknowledgement:**

**Aligning Science Across**

**Parkinson's (ASAP) initiative**

**Grant ID:** ASAP-020607



## Disclaimer

### DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## Abstract

Tile Scan Imaging and Cell Counting



## 1 **\*\*Sample Preparation\*\***

- 1.1 - Prepare 30 µm thick coronal sections containing the anterior cingulate cortex (ACC) and primary motor cortex (MOp) from P21 WT and LRRK2 G2019Ski/ki; Aldh1l1-eGFP mice.

## 2 **\*\*Confocal Microscopy\*\***

- 2.1 - Acquire tile scan images using a confocal Leica SP8 STED microscope equipped with a galvo scanner and a 20x objective.
- 2.2 - Ensure consistent settings for image acquisition across all samples.

## 3 **\*\*Double Positive Cell Identification\*\***

- 3.1 - Identify cells labeled with ALDH1L1-EGFP (astrocyte marker) and SOX9 (astrocyte marker) in the ACC and MOp regions.
- 3.2 - Use FIJI software with the cell counter tool for manual cell counting.

## 4 **\*\*Data Collection\*\***

- 4.1 - Count ALDH1L1-EGFP and SOX9 double positive cells manually in 2-4 sections per brain.
- 4.2 - Analyze sections from 3 sex- and age-matched mice per genotype.