



Jul 13, 2022

# Immunohistochemistry

Alexandra Nelson<sup>1</sup><sup>1</sup>University of California San Francisco1 *Works for me* Share[dx.doi.org/10.17504/protocols.io.14egn7nezv5d/v1](https://dx.doi.org/10.17504/protocols.io.14egn7nezv5d/v1)

kelsey.barcomb

## ABSTRACT

This protocol describes immunohistochemical staining of fixed brain sections.

## DOI

[dx.doi.org/10.17504/protocols.io.14egn7nezv5d/v1](https://dx.doi.org/10.17504/protocols.io.14egn7nezv5d/v1)

## PROTOCOL CITATION

Alexandra Nelson 2022. Immunohistochemistry. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.14egn7nezv5d/v1>



## KEYWORDS

Immunohistochemistry, IHC, Immunostaining, Brain, Mouse, ASAPCRN

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

May 23, 2022

## LAST MODIFIED

Jul 13, 2022

## PROTOCOL INTEGER ID

63075

## 1 Sectioning

1.1 Using freezing microtome section brain at 30-35  $\mu\text{m}$



1.2 Place sections in PBS in a 24-well plate until ready for immunohistochemistry; if not performing immunohistochemistry that day, cover 24 well plate and place at 4°C overnight

2 Blocking.

Incubate sections in blocking buffer at  **Room temperature** for 1-2 hours.

Blocking buffer is: 3% NDS (normal donkey serum)/ 0.1% triton (also called NDST) in PBS.


3 Primary Antibody.

Prepare primary antibody in 3% NDS at desired concentration (make sure this is the final concentration in the well if you are adding more than one primary) and let sections sit on shaker  **Overnight** at  **4 °C**

4 Wash.

Wash sections 3-5 times in PBS for 10-15 minutes each at  **Room temperature**

5 Secondary Antibody.

Prepare secondary antibody in 3% NDS at desired concentration and let sections sit on shaker for 2-4 hours at  **4 °C**

1:500 is default concentration for secondary antibodies.

6 Wash.

Wash sections 3-5 times in PBS for 10-15 minutes each at  **Room temperature**

7 Mounting.

Mount sections onto slides with Vectashield hardset mounting medium (with or without DAPI).

8 Imaging.

- 8.1 For low-magnification display of viral expression patterns or tyrosine hydroxylase staining, we typically use 10X stitched epifluorescence images.
- 8.2 For evaluation of individual neurons and overlap of fluorescent markers, we typically use 40X epifluorescence or confocal Z-stacks.