



Oct 07, 2020

# Tissue delipidation

yang1832<sup>1</sup>, Julia Laskin<sup>1</sup><sup>1</sup>Purdue University

1

Works for me

[dx.doi.org/10.17504/protocols.io.bmw6k7he](https://dx.doi.org/10.17504/protocols.io.bmw6k7he)

NanoDESI MSI Julia Laskin group Purdue



yang1832

## ABSTRACT

This step is to wash lipids off the tissue to prevent their interferences for imaging of proteins.

## DOI

[dx.doi.org/10.17504/protocols.io.bmw6k7he](https://dx.doi.org/10.17504/protocols.io.bmw6k7he)

## PROTOCOL CITATION

yang1832, Julia Laskin 2020. Tissue delipidation. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bmw6k7he>

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Sep 30, 2020

## LAST MODIFIED

Oct 07, 2020

## PROTOCOL INTEGER ID

42686

## ABSTRACT

This step is to wash lipids off the tissue to prevent their interferences for imaging of proteins.

- 1 Immerse the tissue section in 70%, 90%, 100% ethanol, 20s each.
- 2 Immerse in 99.8% chloroform for 25s.
- 3 Dry in desiccator for 10min.